**Engineered for you**

## **Vacuum Aspiration Systems**

**AA – AC – AZ – AE – AK**

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[info@hettichbenelux.com](mailto:info@hettichbenelux.com)

# Content

<b>1. Basic Data</b> .....	<b>4</b>
Safety .....	4
Warranty .....	5
Disclaimer of warranty.....	5
Components.....	5
<b>2. First Steps</b> .....	<b>6</b>
Assembly of systems AA and AC .....	6
Assembly of system AZ.....	7
Setup .....	8
Recommended operating times.....	8
Disposal of liquids .....	9
<b>3. Automatic systems AA</b> .....	<b>9</b>
Pump operation .....	9
Handle .....	10
Connection of additional units.....	11
<b>4. Basic systems AC</b> .....	<b>13</b>
Pump operation .....	13
<b>5. Compact systems AZ</b> .....	<b>14</b>
Pump operation .....	14
<b>6. Aspiration systems AE &amp; AK</b> .....	<b>15</b>
Pump operation .....	15
<b>7. Maintenance</b> .....	<b>15</b>
Cleaning.....	15
Service.....	16
<b>8. Article description accessories</b> .....	<b>18</b>
<b>9. Explanations</b> .....	<b>19</b>

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## Basic Data

The HLC by Hettich Benelux Aspiration Systems are made for the controlled aspiration of liquids from typical lab containers like cuvetts, tubes, microplates, petri dishes or roller bottles. They are also suited for the disposal of liquid biological or chemical waste and filtration of liquids by means of filter systems. If there are any uncertainties regarding the liquids, which should be aspirated, please contact Hettich Benelux. The HLC by Hettich Benelux devices are meant for indoor usage only. Please use only HLC by Hettich Benelux or accessories recommended by HLC by Hettich Benelux.

## Safety

The HLC by Hettich Benelux devices possess an EC declaration, which Hettich Benelux will provide on request. The devices comply with the standards and directives mentioned in the EN DIN 61010-1 norm. IQ, OQ and PQ documents are available on request.

### **Please take note of the following safety measures:**

- » Do only connect the system to a Schuko socket of 230 V, 50 Hz.
- » Attach the safety filter before the vacuum pump. No liquid may get into the pump.
- » Empty the collection bottle in time (max. filling level 75%).
- » The ventilation of the device must not be blocked. To ensure the proper ventilation at all times please consider a minimum distance of 10 centimetres on each side of the device.
- » Damaged mains cable should only be replaced with equal ones.

### **Operational conditions of the system**

- » The device does not own an admission to be placed within an explosiv atmosphere and to extract contents within explosiv atmospheres.
- » The device is not designed for the aspiration of
  - unstable materials
  - materials, which might react explosive under impact (mechanical duress) and / or under elevated temperature without air supply

- self-igniting materials
  - materials, which might explode without air supply
  - explosive material
- 
- » The devices are not authorized for usage in underground mining.
  - » The devices are not designed for the aspiration of dusts.

## Warranty

All functions of the systems have been tested thoroughly. After that, the system and the accessories were packed carefully in perfect condition. If, however, any damages or defects are detected during installation or setup, please contact your local dealer or the Hettich Benelux service department. The Hettich Benelux contact information can be found at [www.hettichbenelux.com](http://www.hettichbenelux.com) (see also our Terms and Conditions at [www.hettichbenelux.com](http://www.hettichbenelux.com)). If used under normal laboratory conditions according to the Operating Manual, Hettich Benelux grants a warranty of 2 years, starting with the day of shipment.


## Disclaimer of warranty

In case of unspecified use of the device, the manufacturer assumes no responsibility.

## Components

The following components are part of every HLC by Hettich Benelux aspirations system:

- » 1 x Autoclavable collection bottle (2 or 4 litres) incl. screw cap
- » 2 x Self-locking fast lock clutches
- » 1 x Ventilation insert
- » 2 x Hose clips
- » 1 x Safety filter
- » 3 m Silicone hose (∅ inside 6mm, ∅ outside 9mm)
- » Y-manifold for connecting several outlets (standard only for AA 24)

## First Steps


### Assembly of systems AA and AC

Place the bottle in the safety bottle holder **(1)**. Cut off a piece of the silicone hose as required for connecting bottle and pump. Cut this hose in two and place the safety filter between bottle and pump, relatively close to the bottle **(2)**.

Insert a free hose end into the hose olive at the back of the pump **(3)**. Put a hose clip on the other end of the hose which is inserted into the self-locking clutch on the cap **(4)**.

Cut off another piece of the silicone hose, the length depending on the desired distance between bottle and workplace. Put the second hose clip on a free hose end which is inserted into the second self-locking clutch on the cap **(5)**.

Put the other hose end in the olive of the handle **(6)**.


## Assembly of system AZ

Cut off a piece of the silicone hose as required for connecting bottle and safety filter. Put the hose clip on a free hose end which is inserted into one of the self-locking clutches on the cap **(1)**. Cut this hose in two and place the safety filter between bottle and pump, relatively close to the bottle **(2)**.

Insert the hose end, which leads away from the safety filter, into the lateral nipple of the separate pump **(3)**.

Cut off another piece of the silicone hose, the length depending on the desired distance between bottle and workplace. Put a hose clip on the end of the hose, which is inserted into the second self-locking clutch on the cap **(4)**. Insert another free hose end into the hose olive and then connect it to your aspiration instrument (e.g. Pasteur pipette) **(5)**.

Connect the mains cable with plug **(6)** and the cable switch **(7)** to a socket with 230 V, 50 Hz.

## Setup

For the systems AA and AC: Plug the mains cable in the IEC socket on the back of the system and connect it with the mains power 230 V, 50 Hz via a Schuko socket. For the system AZ: Connect the mains cable with the cable switch with the mains power 230 V, 50 Hz via a Schuko socket.

## Recommended operating times

### AA 02 / 04 (Pump AP 30)

The pump turns off automatically after reaching its maximum vacuum, therefore the system is suited for continuous operation.

### AC 02 / 04 (Pump AP 15)

If the system gets turned off regularly or a foot switch is used, the system is suited for continuous operation. The maximum approved uninterrupted operating time is 60 minutes, after that, a break for cooling down is required.

### AZ 02 / 04 (Pump AP 04)


The system is suitable for continuous operation, there is no danger of failure or overheating.


## Disposal of liquids

For the quick and easy separation and removal of the collection bottle, it is provided with fast-lock clutches. For unlocking the clutches, push on the lateral metal plate and pull the hose clip out of the clutch on the cap. Now, the bottle is hermetically sealed since the integrated clutches lock automatically. The bottle can be transported safely.

If you want to autoclave the bottle in sealed state, insert the provided ventilation insert into one of the fast-lock clutches on the cap. This opens the bottle and the pressure equalization required for autoclaving is ensured.

If liquid gets into the hose leading to the pump during aspiration, a 5 cm hose can be put on the olive of a fast-

lock clutch inside the cap. This increases the distance between ingress of liquid and vacuum.

Take care that the collection bottle does not exceed app. 75% of its filling volume. If the vacuum gets noticeably weaker, check the permeability of the filter and replace it if required.

## Automatic systems AA

### Pump operation


Switch on the pump with the mains switch "Start / Stop"; the green mains control light turns on, the pump runs and generates a vacuum in the bottle. The pump runs until a vacuum is generated in the bottle or as long as the valve of the handle is opened. It is normal that the pump switches off shortly if you work with very thin aspiration instruments since the vacuum in the

bottle is sufficient for them. Even if the handle is closed, the pump runs temporarily in irregular intervals since it continuously provides the optimal vacuum. During long operation stops you should switch off the pump in order to avoid unnecessary running noises and prevent unnecessary consumption of electricity. Take note that the collection bottle does not exceed 75% of its filling volume if possible. If the vacuum gets noticeably weaker, check the permeability of the filter and replace it if required.

## Handle

The handle AH 01 can be used by right-handed as well as left-handed persons comfortable with one hand. The thumb is e.g. placed on the even surface, the valve is rotated by the index finger self-retaining between 0° and 90°.

If the mark on the rotary knob points to the direction of the flow, the valve is completely open. The valve can be rotated clockwise up to app. 90° until it is completely closed.

If the valve is completely closed, the pump switches off automatically as soon as the desired vacuum is generated in the bottle. In this way, required vacuum is available continuously even without the pump running. For multi-user systems, an unused handle must always be closed since otherwise there is a continuous flow there and now vacuum can be generated. If you use a Pasteur pipette or a suction rake, please cut off a short piece of the enclosed silicone hose and connect your aspirator and the pump with it. High quality steel aspirators can be screwed on the handle. Pipette tips can be connected to the handle with the adapter AR 40, which can be bought separately.

## Cleaning


After the completion of the aspiration process, the handle should be rinsed to avoid a blocking of the valve. During rinsing, slowly close and open the valve. If the valve is blocked, please contact our service department or remove the little lower knurled screw and try to loosen the valve by means of dissolving liquid.

<b>Technical Data</b>		<b>AA 02 / 04, AA 24</b>
<b>Vacuum pump</b>		
Aspiration capacity		17 l / min
Final vacuum		100 mbar abs.
Material of pump head		aluminium
<b>Housing</b>		
Dimensions (W x L x H)		170 x 220 x 155 mm
Height with bottle 02 l		440 mm
Height with bottle 04 l		510 mm
Weight without accessories		5.5 kg
Material housing		aluminium
Power input		100 W
Electr. supply		230 V / 50 Hz
Fuse		2 AT
<b>Environmental conditions</b>		
ambient air temperature		4 – 40°C
in service (non-condensing)		
relative humidity		max. 80%


## Connection of additional units

**This description refers to the connection of additional units to the systems AA 02, AA 04 and AA 24**

Up to 4 additional units can be connected to the aspiration systems. This can be done with 4 bottles with one handle each.


Or with 2 bottles with two handles each - or with a part of each option.

In both cases it is required that the Y-manifold AY 01 of HLC by HETTICH BENELUX is used for connecting the bottles since this ensures a straight-line connection to the fast-lock caps of the bottles (see description “Assembly AA24 for 2 Workplaces”).

### Assembly of the system AA 24 for 2 workplaces

Extend the three connection of the Y-manifold with one short piece of hose each. Then insert the Y-manifold into a fast-lock clutch on the cap **(1)**. Connect each handle with the Y-manifold **(2)**.

If more than one handle is used – for the multi-user systems AA 02 and AA 04 as well as for the two-user system AA 24 – the unused handles must always be closed since otherwise no vacuum can be generated in the bottle.

## Basic systems AC

### Pump operation

Switch on the pump with the mains switch "Start / Stop"; the green mains control light turns on, the pump runs and generates a vacuum in the bottle. During long operation stops you should switch off the pump in order to avoid unnecessary running noises and prevent unnecessary consumption of electricity. Take note that the collection bottle does not exceed 75% of its filling volume if possible. If the vacuum gets noticeably weaker, check the permeability of the filter and replace it if required.

Technical Data	AC 02 / 04
<b>Vacuum pump</b>	
Aspiration capacity	25 l / min
Final vacuum	250 mbar abs.
Material of the pump head	aluminium
<b>Housing</b>	
Dimensions (W x L x H)	170 x 220 x 155 mm
Height with bottle 02 l	440 mm

Height with bottle 04 l	510 mm
Weight without accessories	4.5 kg
Material housing	aluminium
Power input	100 W
Electr. supply	230 V / 50 Hz (115 V / 60 Hz available)
Fuse	2 AT
<b>Environmental conditions</b>	
ambient air temperature in service (non-condensing)	4 – 40°C
relative humidity	max. 80%

## Compact systems AZ

### Pump operation

Connect the hose to the lateral nipple of the pump and turn it on by using the cable switch. The pump runs and generates a vacuum in the bottle. It is reasonable to use the bottle holder and the foot switch AF 01 in addition, these two can be bought as an option. During long operation stops you should switch off the pump in order to avoid unnecessary running noises and prevent unnecessary consumption of electricity. Take note that the collection bottle does not exceed 75% of its filling volume if possible. If the vacuum gets noticeably weaker, check the permeability of the filter and replace it if required.

Technical Data	AZ 02 / 04
<b>Vacuum pump</b>	
Aspiration capacity	8 l / min
Final vacuum	700 mbar abs.
Material of the pump head	aluminium
<b>Housing</b>	
Dimensions of the pump (W x L x H)	75 x 170 x 60 mm
Height bottle 02 l	290 mm
Height bottle 04 l	360 mm
Weight without accessories	0.8 kg

Power input	5 W
Electr. supply	230 V / 50 Hz (115 V / 60 Hz available)
<b>Environmental conditions</b>	
ambient air temperature in service (non-condensing)	4 – 40°C
relative humidity	max. 80%

## Aspiration systems AE & AK

### Pump operation

The AE and AK systems are suited for labs where a vacuum port is present. Connect the hose of the bottle with the vacuum port in your lab. If the vacuum gets noticeably weaker, check the permeability of the filter and replace it if required.

Dimensions	AE & AK 02 / 04
Height bottle 02 l	290 mm
Height bottle 04 l	360 mm

## Maintenance

### Cleaning

Regularly clean the housing of the vacuum safety aspiration systems.

#### Precautions for avoiding electric shock

Electronic devices can cause electric shocks in case of an operating error. Never try to repair electric parts. Never open the housing.

- Switch off the system and disconnect it from the power supply before starting with cleaning or disinfection works.
- Never let get liquids inside the housing (ventilations

slit).

- Do not perform spray disinfection.
- Do only connect the system with the power supply if it is completely dry.

The repair service may only be performed by staff authorized and trained by the manufacturer. A modification of the system is not permitted.

### **Caution when handling aggressive chemicals**

Do not use aggressive chemicals like e.g. strong and weak bases, strong acids, formaldehyde, acetone, halogenated hydrocarbons or phenol for cleaning the system and its accessories.

- In case of contamination with aggressive chemicals, clean the system with a neutral detergent immediately.
- Use neither corrosive detergents nor aggressive solvents or abrasive polishing agents.

### **Cleaning**

1. Disconnect the system from the power supply before you start cleaning.
2. Clean all outer parts of the system with a mild detergent and a lint-free cloth.
3. Wipe off the detergent with Aqua dest.
4. Dry all cleaned parts.

### **Disinfection**

1. Disconnect the system from the power supply before you start disinfecting.
2. Let the system cool down.
3. Clean the system as described above.
4. Select a disinfection method compliant to the applicable local legal regulations and directives.
5. Wipe off all outer parts of the system with the disinfectant and a lint-free cloth.

## Service

If a technical problem arises, please contact your local dealer or the HETTICH BENELUX service department. The HETTICH BENELUX contact information can be found at [www.Hettich Benelux.com](http://www.Hettich Benelux.com). If required, the system will be returned for repair – please take note of our service guidelines, which can be found at [www.Hettich Benelux.com](http://www.Hettich Benelux.com) as well.

### **Decontamination before shipment**

If you send the system to the authorised technical service for repair or to your distributor for disposal, please decontaminate all parts you want to send. Document the decontamination in a Decontamination Certificate inclusive serial number.

## Article description accessories

### Bottles / bottle accessories

Art. No.	Description
800014600	Spare bottle 2 l with fast-lock clutches
800014700	Spare bottle 4 l with fast-lock clutches
800010100	Bottle holder for 2 l bottle
800010200	Bottle holder for 4 l bottle
100021028	Hose clip
800015900	Y-manifold, for connecting several collection units
800015800	Safety filter (shipping unit 3 pcs.)
800014400	Cap for 2 l bottle with fast-lock clutches
800014500	Cap for 4 l bottle with fast-lock clutches


### Handle / adapter

Art. No.	Description
800013200	Handle with flow regulation
800013300	Stainless steel aspiration adapter, length 60 mm
800013400	Stainless steel aspiration adapter, length 150 mm
800013500	Stainless steel aspiration adapter, length 200 mm
800013600	Aspiration adapter for yellow pipette tips
800013700	Adapter with aspiration rake, 8-channel, autoclavable
800015600	Aspiration rake, 8-channel, as spare part

### Pumps & foot-operated switches

Art. No.	Description
800015300	Membrane vacuum pump
800015000	Membrane vacuum pump
800015100	Membrane vacuum pump
800014000	Foot-operated switch for switching the pump on and off

## Explanations

	Caution! Risk of electric shock!
	Caution!
	Important information
	Earth conductor
	Fuse

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