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

# Directions for Use AccuOne | EnergyOne

Art. No. 5900-1050, -1070, -1100, -2050, -2070, -2100

Read carefully before using.

## Identification of dangers in the Directions for Use:

**Warning!**

This symbol identifies actions that may pose dangers to life and limb of persons if they are not performed properly.

**Caution!**

This symbol identifies actions that could result in possible dangers for objects.

Identification on the device:

This sign is used to indicate that depending on the liquid being pumped, personal protective equipment must be worn. The operator must determine requirements for specific cases and ensure the requirements are implemented (safety data sheet, hazard analysis).

Comply with the warnings on the device. Always ensure that the warnings and mandatory signs are in proper condition. Replacement signs can be obtained from the manufacturer.

## 1. Safety

The AccuOne/EnergyOne decanting and emptying system is designed, constructed and tested for operating safety according to the latest state of the art. It left the factory in flawless safety condition. Nevertheless, the device may pose hazards for persons and objects if it is not operated properly. It is therefore essential to read the operating manual from cover to cover and to comply with the safety instructions. In case of non-designated use, the manufacturer shall assume no liable and warranty claims may not be accepted.

### 1.1 Designated use

- ▶ The AccuOne decanting and emptying system must only be used for decanting and emptying liquids in barrels and canisters.
- ▶ The pump is not designed for use in a potentially explosive atmosphere or with readily flammable liquids. Any such use is prohibited!
- ▶ The device must not be used to pump organic solvents.
- ▶ The device must only come in contact with media against which the materials are resistant. (For materials see the technical data;

for resistance see the list entitled "Chemical resistance" under [www.buerkle.de](http://www.buerkle.de)

- ▶ The pump is a non-hand-held, portable motorised device. It must only be used under direct personal supervision.
- ▶ The position of the pump in accordance with designated use is vertical.
- ▶ Continuous operation is limited to 20 minutes.
- ▶ The drive must not be immersed in the pumping liquid.
- ▶ The pump must only be used in buildings or in an area where it is protected from the elements.

### 1.2 Working and operational safety

**The pump is not designed for use in a potentially explosive atmosphere or with readily flammable liquids. Any such use is prohibited!**

- ▶ The operator is responsible for ensuring that only persons who have read and understood the operating instructions are permitted to use the device.
- ▶ Check the device for damage every time

before commissioning. Any damage must be repaired before continuing with operation.

- ▶ If the pump station is disconnected, the drive must not be turned on (danger of finger injuries).
- ▶ The device must only be switched on to pump liquids. In particular the pump must not be in operation if the pump tube has been removed from the container.
- ▶ When the pump tube is removed from the filling medium, make certain that a suitable container is available for the residual liquid that will flow out.
- ▶ If the discharge hose is connected, fluid may continue to run out even if the pump is not in operation. To prevent this, the nozzle must be closed and the end of the hose must be positioned above the container that is being filled.
- ▶ The discharge tube must always be pointed into the opening of the collecting container while the pump is in operation.
- ▶ Conversions or modifications to the device are prohibited. Use only the manufacturer's original spare parts and accessories.

### 1.3 Hazardous materials

If hazardous materials are used during operation or maintenance or while the device is being cleaned observe the current safety data sheets.

Devices that are contaminated with hazardous materials must not be returned to the device manufacturer for repairs.



## 2. Commissioning/operation


**Danger of injury due to hazardous materials! Always wait to switch on the device until it is correctly mounted in the barrel or canister as described in the mounting instructions. Ensure that all safety precautions have been taken.**

The pump must be switched off during mounting. Follow these steps:


1. Connect the discharge tube (Fig. 1/E) or discharge hose (Fig. 1/F) to the discharge threads (Fig. 2 and 3). Make certain the seal (Fig. 1/H) is inserted and tighten the union nut finger tight.




2. Move the adapter reduction (Fig. 4/A) and suitable and/or required thread adapters (Fig. 4/B) or barrel plugs (Fig. 5) over the pumping tube of the pump station. Place the pump station in the container (Fig. 6) and screw the thread adapter together (Fig. 7).

3. The barrels and canisters must be stable.
4. Power is supplied to the drive either via the rechargeable battery (Fig. 1/A) or through the power supply unit (Fig. 1/C). Standard plug-in adapters for the relevant country must be used for the power supply unit and charger. Connect the drive with the charged rechargeable battery or power supply unit (Fig. 8) on the pump station.


5. Place a suitable, stable container under the discharge tube or hose.
6. Move the centre On/Off key  to switch the state of the device to ready for operation (Fig. 9/B).
7. Press the  key (Fig. 9/A) to increase the speed of the motor. Press and hold the key until the required volume flow is reached. Use the  key (Fig. 9/C) to reduce the volume flow.
8. Press the centre On/Off key  again to switch off the drive.

**! Caution!**

**The drive must always be switched off before it is disconnected from the pump station!**

To remove the drive from the pump station, simultaneously press the two keys on the discharge section of the pump station and pull out the drive (Fig. 10).


**! Caution!**

**The rechargeable battery or power supply unit must not be disconnected from the drive while the pump is in operation (risk of damaging the electronics).**

**The pump must not run dry. Switch off the pump immediately when the container is empty (risk of damaging the pump).**

**Never turn the pump tube and the connection head against each other (Fig. 11)!**

## 2.1 The rechargeable battery

The rechargeable battery consists of lithium ion cells and can be charged at any time. The charger included with delivery prevents the rechargeable battery from overcharging. To prevent deep discharge, the rechargeable battery must be charged at least once every 12 months.

### To charge the rechargeable battery:

1. Press the keys on the side to remove the rechargeable battery from the drive (Fig. 12).

2. Connect the charging cable to the charging socket of the rechargeable battery (Fig. 13).

3. Connect the charger to the socket. The green lamp flashes during the charging process. The rechargeable battery is completely charged when the lamp is continuously lit. The charging time is about 1 hour at a mains voltage of 230V.

**! Caution!**

**Only the charger included with delivery may be connected to the plug connection (charging socket) on the rechargeable battery. An unsuitable charger could damage and overcharge the rechargeable battery (danger of fire).**

**The rechargeable battery must not be exposed to shocks and impacts. If it is, it could be damaged and high heat could develop as a result (danger of fire).**

## 3. Cleaning, maintenance and upkeep

The device must only be maintained and cleaned as specified in these directions for use. The pump is maintenance-free.

Please use suitable cleaning fluids for cleaning depending on the medium that is being pumped. Make certain the materials of the device are not detrimentally affected.

To clean the device, rinse and flush out the pump station. If necessary the inside of the pump station can be cleaned by pumping cleaning fluid through it. Observe the safety instructions noted previously.

The motor and rechargeable battery must not be immersed in liquid.

The coupling in the connection head must not be exposed to any moisture and therefore must not be cleaned wet.

## 4. Troubleshooting

Troubleshooting and fault elimination must only be performed when the power supply unit or rechargeable battery is disconnected.

### The device does not start

- ▶ Device is not switched on: Press the On/Off key; LED lamp indicates ready for operation
- ▶ Rechargeable battery fully discharged: Charge rechargeable battery
- ▶ Power supply unit not connected: Connect to socket
- ▶ Drive faulty: Send the device in to the manufacturer
- ▶ Drive blocked: Send the device in to the manufacturer

### Drive running, but not pumping any liquid.

- ▶ Viscosity of the liquid too high (see Technical data)
- ▶ Pump station clogged: Rinse or clean pump station

### Fluid escaping between drive and pump station

- ▶ Switch-off the pump immediately. The shaft seal is probably faulty. Replace the pump station or send it in to the manufacturer.
- ▶ If repairs are necessary, the device can be sent in to the manufacturer or repairs can be performed by an electrician.

## 7. Disposal

Devices that are contaminated with hazardous materials must be cleaned or disposed of as hazardous materials!

The materials of the cleaned device are recyclable. They can be presented to appropriate collecting stations.

## 6. Accessories, additional and spare parts AccuOne/EnergyOne

	Art. No
Ⓐ Rechargeable battery without charger	5900-3010
Ⓑ Charger	5900-3015
Ⓒ Power supply unit incl. connector	5900-3020
Ⓓ Drive without rechargeable battery, power supply unit and pump station	5900-3030
Ⓔ Discharge tube	5900-3040
Ⓕ Discharge hose with nozzle	5900-3045
Ⓖ Pump station with immersion pipe length 500 mm	5900-3050
Ⓖ Pump station with immersion pipe length 700 mm	5900-3070
Ⓖ Pump station with immersion pipe length 1000 mm	5900-3100

### Barrel connections

R2", steel barrel thread, internal	5600-3130
Container barrel screw connection for 59.5-mm opening Ø	5600-3140
Barrel plug PVC for 40 – 70-mm opening Ø	5600-3150
Mauser 2", coarse thread, internal	5600-3170
Tri-Sure, coarse thread, internal	5600-3180

## 7. Technical data

### Pump station

Parts in contact with medium	PP, PTFE, stainless steel (Hastelloy C)
Discharge tube	PP
Discharge hose	PVC with PP MS nozzle
Immersion pipe length	50, 70, 100 cm
Immersion tube Ø	25 mm
Connecting threads	¾"
Delivery volume in reference to water (density 1 g/cm <sup>3</sup> )	12 l/min
Max. viscosity of the pumping medium	200 mPas
Max. density of the pumping medium	1.5 g/cm <sup>3</sup>

### Drive

Rechargeable battery	14.4 V 1.5 Ah lithium ion
Capacity of rechargeable battery when pumping water	Up to 600 l
Charger	14.4V; 1400 mA; In 230V; Out 16.8V; Protection class II
Power supply unit	30 W; In 100 V/230V; 50 – 60 Hz; Out 15V; Protection class II
Drive	12 V max. 2.5A IP 54
Airborne noise emission	≤ 55 dB (A)
Operating ambient temperature	+5 ... +40 °C average temperature/24 h 35 °C
Relative humidity	Max. 90 % at 20 °C, 50 % at 40 °C
Storage and transport temperature	-25 °C to +55 °C
Altitude	≤ 2000 m above sea level
Weight	max. 1.2 kg

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