winkler®

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Laborheizhauben

OPERATING INSTRUCTIONS

Laboratory heating mantles

WHG1.../G2...
WHG2R.../G4R...
WHG2H...
WHLG2...
WHLG/ER...
WHLMG3/ER...
WHLSG2/ER...
WHU2/ER...
WHRE...

winkler®

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Laborheizhauben

**WHG1.../G2...
WHG2R.../G4R...
WHG2H...
WHLG2...
WHLG/ER...
WHLMG3/ER...
WHLG2/ER...
WHU2/ER...
WHRE...**

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IMPORTANT SAFETY INSTRUCTIONS REGARDING USE AND OPERATION

Warning:
Dangerous electric voltage!

Important safety instructions regarding use and operation

Please read the operating instructions very carefully before using the device. Observe the information on the type plate and the warnings on the device. For subsequent use and operation of the device, the operating instructions should be kept in a safe place where they must always be visible and available to the operator.

PILZ® heating mantles are to be operated in compliance with the relevant legislation on occupational health & safety, the country-specific requirements and the stipulations of the Employers' Liability Insurance Associations (e.g. TRGS 526 – Technical Rules for Hazardous Substances in Laboratories). Please include the present instructions into your Quality Management Manual as part of the work instructions. They are to be treated as operating instructions according to TRGS 526, section 4.2.

Warning:
Dangerous area!

Warnings on the product must not be removed!

The heating mantle is an electrical appliance.

Operation and maintenance of PILZ® heating mantles must only be carried out by a qualified electrician or a person with electrotechnical training.

Electric shocks may result in serious injuries.

The PILZ® heating mantle is provided with a residual-current device (RCD) in the power plug in order to prevent electric shock hazards. This device will switch off the heating mantle in the event of fault occurrence (fault current).

The correct functioning of the residual-current device needs to be checked before putting the heating mantle into operation and such checks must be repeated at regular intervals (see page 6, functional test of the residual-current device).

Do not touch!

No moisture protection

The PILZ® heating mantle is not protected against moisture penetration. Moisture penetrating into the heating mantle – such as water, liquid chemical substances or other electroconductive media used during the operating of the PILZ heating mantle – can result in dangerous electric shocks. In the event of spillage, flask breakage or overflowing liquid the heating mantle must no longer be touched.

The residual-current device (RCD) integrated into the connecting mains lead or power plug will switch off the PILZ® heating mantle in response to a fault current flow of > 30 mA, so that a dangerous touch voltage at accessible parts cannot occur. The PILZ heating mantle must then be disconnected from the mains by pulling the plug out of the grounded socket (Schuko).

Warning:
Flammable

The PILZ® heating mantle must afterwards not be put back into operation!

The defective PILZ® heating mantle must be marked accordingly and stored in a safe place. It can only be released for reuse by the electrician in charge or a person with electrotechnical training.

The PILZ® heating mantle may be dried in a drying cupboard. However, before being put back into operation it has to be completely dry and must first be checked by an electrician.

Other cases when the PILZ® heating mantle must be taken out of operation immediately:

Visible heat conductor in internal section > mechanical damage

Brittle or broken glass fabric > overheating of glass fabric

In the event of this kind of damage or any other similar damage, the PILZ® heating mantle has to be taken out of operation immediately and reuse is no longer permitted. It must then be disposed of so that any further operation is excluded.

Warning:
Hot surface!

Precautions against overheating

Take care to avoid overheating of the PILZ heating mantle. Overheating can be avoided by using suitable temperature controllers during operation or – where the device is equipped with several heating zones – by switching off individual heating circuits. When the liquid level goes down, any heating circuits that are not required need to be switched off or else the total heating power must be reduced.

Warning:
Dangerous electric voltage!

PILZ® heating mantles must not be operated with empty or encrusted/scaly flasks or with no flask at all.

Do not touch the PILZ® heating mantle when it is switched on or ready for operation. The heating surface must not be touched with metal or other electroconductive objects.

Do not use any metal vessels, vessels with metal parts or other electroconductive vessels in the area of the heating surfaces. Do not stick any needles or similar into the PILZ® heating mantle. This could damage the heating coil or other electrical conductors and/or the corresponding protective systems. Furthermore, it could result in electric shocks and serious injuries.

The PILZ® heating mantle is fitted with a heating element along its internal section. This heating element must not be touched when the PILZ® heating mantle is switched on or when it is still hot. The upper edge of the PILZ® heating mantles and the fastening eyelets can reach high temperatures. There is a risk of injury due to burns.

Warning:
Explosive substances

This type of PILZ® heating mantle is not explosion protected! It must not be used for heating processes where there is a hazard due to explosive media or explosive gas-air mixtures.

You should also take into account that flammable liquids can cause a fire in the event of flask breakage, overfilling or spillage.

GENERAL OPERATING INSTRUCTIONS

PILZ® heating mantles for laboratory and industrial applications will function to your satisfaction provided you observe the following information. Our products are manufactured to a high quality standard; they are systematically monitored and checked during all stages of production and only leave our works after a final safety test. As a rule, PILZ® heating mantles must be employed in accordance with the valid regulations.

PILZ® heating mantles are glass-silk insulated heating systems designed for use in industry and laboratories. The heating mantles are not protected against moisture because of their high temperature range. The flexible heating units can be adapted to the shapes and sizes of the corresponding vessels and thus ensure uniform heat transfer as well as careful treatment of the laboratory vessels.

All PILZ® heating mantles are tested and fully heated up before leaving our works. Any resulting discoloration of the glass fabric is technically induced and does not represent a deterioration in quality.

Intended use of PILZ® heating mantles

PILZ® heating mantles are designed for heating non-combustible, liquid or solid substances or mixtures in electrically non-conductive glass vessels (flasks).

Note: In line with the relevant regulations, PILZ® heating mantles are not suitable for heating electroconductive vessels.

Only use glass vessels (flasks) that match the corresponding heating mantles in size and shape. PILZ® heating mantles may only be operated by a qualified electrician or a person with electrotechnical training. The accepted rules and standards (VDE, SEV, ÖVE, DIN standards, and particularly the accident prevention regulations) must be complied with.

Stable positioning of the PILZ® heating mantle

The operator of a PILZ® heating mantle must take appropriate measures to ensure that the heating mantle is safely positioned (e.g. fitting into a support, tripod). PILZ® heating mantles that have been built into corresponding housings by the manufacturer do not require any other positioning system. In the case of PILZ® heating mantles of type series WHLG.. / WHLMG.. / WHLSG.. / WHU.. (metal housings) the screwed-on rod clamp may only be loaded in such a way that there is no risk of tilting or mechanical damage (e.g. use an additional stand with a broad base to guard against tilting).

Caution: PILZ® heating mantles must always be in an upright, level position during operation. Any objects on top of the glass vessel endanger the stability and need to be specially secured.

PILZ® heating mantles are suitable for carrying the deadweight of the glass vessels. With further setups, the additional weight must be taken up by suitable supports. Heating mantles of series WHG... are provided with special hooks for mounting on a tripod with support ring. Heating mantles of series WHLG... , WHLMG... , WHLSG... , WHU... are built into self-standing metal housings.

Transport

PILZ® heating mantles have to be suitably packed and transported with care. The connecting lead should never be pulled or used for carrying the heating mantle. Sharp objects that could cause damage should be kept away from the device. To avoid damage (fire risk) the PILZ® heating mantle must have cooled down sufficiently before being packed for transport.

Operation of PILZ® heating mantles

PILZ® heating mantles may only be operated in suitable and well-ventilated rooms. During operation all combustible and easily flammable objects must be no closer than 1 m from the PILZ® heating mantle.

Grounding-type plug (Schuko) as a circuit breaker (residual-current device)

In the case of PILZ® heating mantles with a built-in power plug, a residual-current device in the lead serves as a suitable circuit breaker.

GENERAL OPERATING INSTRUCTIONS

Temperatures

The temperature stated of the type plate of the PILZ® heating mantle represents the max. permissible surface temperature of the heating appliance. The operator must take appropriate measures to ensure that this temperature is not exceeded at any point.

Caution: In the case of operation without temperature control there is a risk that the max. operating temperature will be exceeded.

This can be caused by:

Heat concentration because of empty or partially filled vessels, varying filling levels and changes in ambient conditions (temperatures). PILZ® heating mantles are damaged or even destroyed by excessive temperatures (fire risk). Furthermore, it should be realized that the stipulated temperature range will also be exceeded if the heating systems are operated with unsuitable or incorrectly set temperature controllers and temperature sensors. For safety reasons, PILZ® heating mantles should always be operated together with a temperature controller.

Important: Both the vessel and the medium to be heated must be suitable for the applicable max. operating temperature. The specifications of the vessel and medium manufacturers must be adhered to.

Where sensitive or dangerous media are involved, the operator must take suitable measures in line with the technical requirements to prevent excessively high temperatures and damage to the medium or the PILZ® heating mantle. As a rule, the temperature sensor must be placed at the hottest point according to the particular application. Where a power selector is employed, the process must be constantly monitored by the operator using a temperature measuring instrument and adjusted if necessary.

Check if the operating temperatures are being correctly maintained during operation and if the total system functions normally.

Penetration of medium into the PILZ® heating mantle

If a medium penetrates into a PILZ® heating mantle because of spillage, glass flask breakage or overflow, this can produce a dangerous touch voltage. The residual-current device RCD (FI) > 30mA, which is integrated into the power plug, will switch off the heating mantle in response to fault current flow, so that a dangerous touch voltage at accessible parts cannot occur. Nevertheless, the PILZ® heating mantle must then be disconnected from the mains by pulling out the power plug. The heating mantle must not be put back into operation and should be stored in a suitable place without allowing reuse. The PILZ® heating mantle may only be put back into operation after a safety test by a qualified electrician (operator) or the manufacturer. The type of overloading that has occurred must be determined and the cause eliminated.

Contamination

In the event of contamination, the PILZ® heating mantle must be disconnected from the network immediately. The type of contamination load must then be determined. In cooperation with a competent specialist, the PILZ® heating mantle must immediately be stored in a suitable place and secured against reuse.

Maintenance

Like all electrical equipment, PILZ® heating mantles must be checked by the operator at regular intervals to ensure that they are in a safe condition (VDE test according to BGV A3 and BetrSicherV).

Maintenance and service are to be carried out in compliance with the standards stated under "Safety Instructions" and the relevant stipulations of the Employers' Liability Insurance Associations, as well as the regulations applicable to the particular application. Any damaged appliances must immediately be taken out of operation and may not be reused.

Repairs

PILZ® heating mantles are made of special materials which have to meet high requirements in terms of temperature as well as electrical and chemical stability. This also applies to the connecting lead. Repairs may only be performed by the manufacturer or an authorized company.

Electric heating systems are to be installed and operated in such a way that unattended operation or accidental switching on does not represent a danger. This makes it necessary to put suitable safety measures into place, e.g. operation using a controller and/or cut-out.

Where other rules are being applied, it is the duty of the operator to provide appropriate protective measures.

Switching options of PILZ® heating mantles

Depending on the type of design, PILZ® heating mantles are fitted with an ON-OFF switch, a series/parallel switch or a heating zone switch, step switch or power selector (see following table). Please check in advance what type of PILZ® heating mantle you are going to use and ensure that the relevant switch is off, so that the device does not start heating accidentally. After plugging in and switching on the residual-current device (RCD), the heating mantle can be turned on at the appropriate switch. The PILZ® heating mantle features different heating zones, depending on design.

GENERAL OPERATING INSTRUCTIONS

Heating mantles with metal housings

This heating mantle series is equipped with an ON-OFF switch.

WHLG1	1 heating zone	2 power steps	min/max switch 50% / 100%	
WHLG2	2 heating zones	2 power steps	min/max switch 50% / 100%	
WHLG1/ER	1 heating zone	2 power steps	min/max switch 50% / 100%	power selector 0-100% electronic
WHLG2/ER	2 heating zones	2 power steps	min/max switch 50% / 100%	power selector 0-100% electronic
WHLMG3/ER	3 heating zones	3 steps	zone 1; zone 2; zone 3	power selector 0-100% electronic
WHLG2/ER	2 heating zones	2 power steps	min/max switch 50% / 100%	power selector 0-100% electronic
WHU/ER	2 heating zones	2 power steps	min/max switch 50% / 100%	power selector 0-100% electronic

Heating mantles with flexible outer jacket

Zone switch and series parallel switch with ON-OFF position

WHG1	1 heating zone	2 power steps	zone switch 50% / 100%
WHG2	2 heating zones	3 power steps	series/parallel switch 0 / 50 / 25 / 100%
WHG2R	2 heating zones	3 power steps	series/parallel switch 0 / 50 / 25 / 100%
WHG4R	4 heating zones	4 steps	zone 1, zone 2, zone 3, zone 4
WHG2H	2 heating zones	3 power steps	series/parallel switch 0 / 50 / 25 / 100%

Multi-place heating mantles

The devices are equipped with an ON-OFF switch.

WHRE4	1 heating zone per heating place	4 heating places	1 power selector per heating place	(OFF) 0-100% electromechanical
WHRE6	1 heating zone per heating place	6 heating places	1 power selector per heating place	(OFF) 0-100% electromechanical

Heating mantle types

WHLG1; WHLG1/ER; WHG1	WHG2; WHG2R; WHG2H	WHLG2; WHLG2/ER; WHLSG2; WHU/ER	WHLMG3/ER; WHG4R
Number of heating zones			
1 	2 	2 	3 und 4 
Zone switch and min/max switch	Series/parallel switch	Min/max switch	Step switch
Power			
0 1/2 1/1	0 1/4 1/2 1/1	0 1/2 1/1	0 1/3 2/3 3/3 0 1/4 2/4 3/4 4/4

To check prior to putting the heating mantle into operation:

1. Do your ordering data correspond to the type plate?
2. Is the PILZ® heating mantle suitable for your application according to the technical specifications (type plate, description)?
3. Are there any measures in place to prevent excessive temperatures?
4. Is the PILZ® heating mantle safely positioned?
5. Has it been ensured that you can disconnect the PILZ® heating mantle from the mains in the event of danger?
6. Has the flask or vessel been cleaned? Incrusted or scaly flasks possess poor thermal conductivity and can cause flask breakage.

Important: Make sure to protect the PILZ® heating mantle (incl. switch) against moisture and wet conditions.

GENERAL OPERATING INSTRUCTIONS

Residual-current device (RCD) / functional test

Putting into operation: Put plug into socket and press green "RESET" button. The device is operative when the red functional indicator lights up. After unplugging or power failure, the device will switch itself off automatically.

A functional test must always be carried out before putting the device into operation. In the event of repeated failure, the connected device must be checked. Please observe that the residual-current device is no substitute for fundamental safety measures.

Functional test procedure of the residual-current device:

Press blue "TEST" button; the device will switch off.

Press green "RESET" button; the device will again be operative when the red functional indicator lights up.

Operation

First place a suitable glass vessel into the PILZ® heating mantle and check that the glass vessel is correctly seated in the heating mantle. Put the mains plug (with residual-current device) into the socket and switch on the heating mantle.

The table under "Switching options of PILZ® heating mantles" shows the different switch settings.

Important: Do not leave the PILZ heating mantle unattended during operation. Observe all the rules and stipulations of the present operating instructions.

The technical data of the PILZ heating mantles are stated on the type plate of the heating mantle.

Safety class of housing (except heating surface):	I	protective conductor
Safety class of heating surface / crocheted outer jacket:	IP 00	No moisture protection
Safety class of housing (series with metal housing):	IP 43	
Safety class of power plug:	IP 44	Residual-current device
Max. heat conductor temperatures	450 °C	(series WHG2H up to 900 °C)

Dimensions/ further technical data

> See Laboratory Catalogue

The operating instructions must be fully and correctly followed, otherwise the PILZ® heating mantle may be not protected as intended.

Please contact the manufacturer if you have any questions.

GENERAL OPERATING INSTRUCTIONS

Our safety connection system

The new generation of PILZ® heating mantles offers the operator optimal protection during day-to-day operation thanks to our combinable safety connection system. As a rule all PILZ® heating mantles (laboratory heating mantles) are equipped with a multipole plug.

Fig. 1 – WHG2 crocheted heating mantle

Fig. 2 – WHLG2 metal-cased heating mantle

Fig. 3 – WHLG2/ER metal-cased heating mantle

The heating mantles are connected to the mains supply network via a suitable connection system that plugs into a grounded socket (Schuko).

The **WHZA 1000** residual current device (RCD) is connected to the PILZ heating mantle via a multipole coupling as a standard (Fig. 4).

If the RCD is not required because the heating mantle is operated on a protected circuit (e.g. in a laboratory), the PILZ® heating mantle may also be connected directly to the supply network via the **WHZA 2000** Schuko plug (Fig. 5, optional accessory).

Fig. 4 – Connection WHZA 1000 (standard)

Fig. 5 – Connection WHZA 2000 (optional accessory)

Important:
When using the WHZA 2000 connection kit, ensure that the power supply circuit is protected by a residual current device (RCD) before putting the PILZ® heating mantle into operation.

Safety connection system and temperature control

PILZ® heating mantles can be connected to appropriate control equipment, e.g. a power selector or temperature controller.

The illustrations below show different safety connection options.

Please contact us if you have any questions; our specialists will be happy to assist you.

Fig. 6 – Power selector WRT00117

Fig. 7 – Microprocessor-based controller WRT2010X

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are in conformity with

the European Directives

EC Directive 2006/95/EC "Low Voltage Directive" (codified version)
"Electrical equipment for use within certain voltage limits"
(formerly DIR 73/23/EEC of 19.02.1973 amended by DIR 93/68/EEC)

EC Directive 89/336/EEC "Electromagnetic compatibility"
amended by DIR RL 91/263/EEC, 92/31/EEC and 93/68/EEC.

The following national standards are adhered to with respect to the products described:

DIN EN 60519-1/-2 Safety in electroheat installations
(VDE 0721-1/-2) Part 1: General requirements / Part 2: Particular requirements...

DIN EN 61010-1/-2 Safety requirements for electrical equipment for measurement, control and laboratory use
(VDE 0411-1/-2) Part 1: General requirements / Part 2: Particular requirements...

This declaration does not constitute a warranty of properties within the meaning of the product liability law.

The safety instructions of the product information must be observed!

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**Winkler GmbH