◀ DUPONT ▶™

Tyvek®

For greater good™

500

Cat. III PROTECTION LEVEL

ACCESSORIES

2 | ◀ DUPONT ▶ | 8

1 | Tyvek® 500 |

3 | Hood model PH30L0 Boot cover model POBO
Jacket model PP33 Boot cover with antislip model POBA
Trousers model PT31L0 Shoe cover model POSO
Labcoat with press studs and pockets model PL30 Shoe cover with antislip model POSA
Labcoat with press studs model PL30NP Sleeve model PS32LA
Labcoat with zipper and pockets model PL309 Apron model PA30L0
Labcoat with zipper model PL309NP

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© 2021 DuPont. All rights reserved. DuPont™, the DuPont Oval Logo, and all trademarks and service marks denoted with ™, SM or ® are owned by affiliates of DuPont de Nemours, Inc. unless otherwise noted.

Internet: dpp.dupont.com

DuPont de Nemours (Luxembourg) s.à r.l.
L-2984 Luxembourg

Cert. Ref: Tyvek® 500 Accessories
January 2019/26/V2

DuPont Ref: IFUAT_012

Overview of the accessories

Tyvek® 500 Hood model PH30L0

- Hood with flange
- Elasticated face and neck

Seams: Stitched internal seams

Fabric: **Tyvek® 500**

Article measurements in cm

Size	Hood width	Hood length side
Single	60	57

Tyvek® 500 Jacket model PP33

- Hooded jacket
- Zipper closure

Seams: Stitched internal seams

Fabric: **Tyvek® 500**

Body measurements in cm (see table)

Sizes	M - 2XL

Tyvek® 500 Trousers model PT31L0

- Elasticated waist
- No elastic at ankles
- Without pockets

Seams: Stitched internal seams

Fabric: **Tyvek® 500**

Body measurements in cm (see table)

Sizes	M - 2XL

Tyvek® 500 Labcoat with press studs and pockets model PL30

- Labcoat with collar
- 5 press studs
- 3 pockets

Seams: Stitched internal seams

Fabric: **Tyvek® 500**

Body measurements in cm (see table)

Sizes	M - 2XL

Tyvek® 500 Labcoat with press studs model PL30NP

- Labcoat with collar
- 5 press stud closures
- Without pockets
- Elasticated cuffs (not tunneled)

Seams: Stitched internal seams

Fabric: **Tyvek® 500**

Body measurements in cm (see table)

Sizes	M - 2XL

Tyvek® 500 Labcoat with zipper and pockets model PL309

- Labcoat with collar
- Zipper closure
- 2 pockets
- Elasticated cuffs (tunneled)

Seams: Stitched internal seams

Fabric: **Tyvek® 500**

Body measurements in cm (see table)

Sizes	S - 2XL

Tyvek® 500 Labcoat with zipper model PL309NP

- Labcoat with collar
- Zipper closure
- Without pockets
- Elasticated cuffs (tunneled)

Seams: Stitched internal seams

Fabric: **Tyvek® 500**

Body measurements in cm (see table)

Sizes	S - 2XL

Tyvek® 500 Boot cover model POB0

- Knee-length overboot
- Elasticated top
- Fixation ties

Seams: Stitched internal seams

Fabric: Tyvek® 500

Article measurements in cm

Size	Boot length	Boot height
Single	35	48

Tyvek® 500 Boot cover with antislip model POBA

- Knee-length overboot
- Elasticated top
- Fixation ties
- Slip-retardant sole

Seams: Stitched internal seams

Fabric: Tyvek® 500

Article measurements in cm

Size	Boot length	Boot height
Single	35	48

Tyvek® 500 Shoe cover model POS0

- Overshoe
- Elasticated ankle

Seams: Stitched internal seams

Fabric: Tyvek® 500

Article measurements in cm

Size	Shoe length	Shoe height
Single	38	22

Tyvek® 500 Shoe cover with antislip model POSA

- Overshoe
- Elasticated ankle
- Slip-retardant sole

Seams: Stitched internal seams

Fabric: Tyvek® 500

Body measurements

	UK	EU
Shoe size	3 1/2 - 7 1/2	36-42
Shoe size	7 1/2 - 11	42-46

Tyvek® 500 Sleeve model PS32LA

- Adjustable arm opening
- Upper-arm in blue-coloured thread for identification purposes

Seams: Stitched internal seams

Fabric: Tyvek® 500

Article measurements in cm

Size	Length	Upper arm opening diameter	Wrist opening diameter
Single	50	13-24 (elasticated)	8-15 (elasticated)

Tyvek® 500 Apron model PA30L0

- Shin-length apron
- Neck and waist ties

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
ENGLISH

INSTRUCTIONS FOR USE

INSIDE LABEL MARKINGS ① Trademark. ② PPE manufacturer. ③ Model identification for accessories made of Tyvek® 500 fabric. This instruction for use provides information on these accessories. ④ CE marking - The Partial Body (PB) accessories comply with the requirements for category III personal protective equipment according to European legislation, Regulation (EU) 2016/425. Type-examination and quality assurance certificates were issued by SGS Fimko Oy, Takomitie 8, FI-00380 Helsinki, Finland, identified by the EC Notified Body number 0598. ⑤ Indicates compliance with European standards for chemical protective clothing. ⑥ Partial body protection "types" achieved by these accessories defined by the European standards for chemical protective clothing: EN 13034:2005 + A1:2009 (Type PB [6]). These accessories also fulfill the requirements of EN 14126:2003 Type PB [6-B]. ⑦ Country of origin. ⑧ Date of manufacture. ⑨ Flammable material. Keep away from fire. These accessories and/or fabrics are not flame resistant and should not be used around heat, open flame, sparks or in potentially flammable environments. ⑩ Sizing pictogram indicates body or article measurements (cm) & correlation to letter code. Check your body, arm or foot measurements and select the correct size. ⑪ Wearer should read these instructions for use. ⑫ Do not re-use. ⑬ Other certification(s) information independent of the CE marking and the European notified body (see separate section at end of the document).

TYVEK® 500 ACCESSORIES:

FABRIC PHYSICAL PROPERTIES			
Test	Test method	Result	EN Class*
Abrasion resistance	EN 530 Method 2	> 100 cycles	2/6**
Flex cracking resistance	EN ISO 7854 Method B	> 100 000 cycles	6/6**
Trapezoidal tear resistance	EN ISO 9073-4	> 10N	1/6
Tensile strength	EN ISO 13934-1	> 30N	1/6
Puncture resistance	EN 863	> 10N	2/6

* According to EN 14325:2004 ** Visual end point

FABRIC RESISTANCE TO PENETRATION BY LIQUIDS (EN ISO 6530)		
Chemical	Penetration index - EN Class*	Repellency index - EN Class*
Sulphuric acid (30%)	3/3	3/3
Sodium hydroxide (10%)	3/3	3/3

* According to EN 14325:2004

FABRIC RESISTANCE TO PENETRATION OF INFECTIVE AGENTS		
Test	Test method	EN Class*
Resistance to penetration by blood and body fluids using synthetic blood	ISO 16603	3/6
Resistance to penetration by blood-borne pathogens using bacteriophage Phi-X174	ISO 16604 Procedure C	no classification
Resistance to penetration by contaminated liquids	EN ISO 22610	1/6
Resistance to penetration by biologically contaminated aerosols	ISO/DIS 22611	1/3
Resistance to penetration by biologically contaminated dust	ISO 22612	1/3

* According to EN 14126:2003

SEAM PROPERTIES		
Test method	Test result	EN Class*
Seam strength (EN ISO 13935-2)	> 50N	2/6

* According to EN 14325:2004

For further information about the barrier performance, please contact your supplier or DuPont: dpp.dupont.com

RISKS AGAINST WHICH THE PRODUCT IS DESIGNED TO PROTECT: These accessories only provide partial body protection and are intended to protect those parts of the body exposed to risk. They are designed to help protect sensitive products and processes from contamination by people and to protect workers from some hazardous substances. They may be used separately or in combination with other Personal Protective Equipment to increase the protection level needed. Fabric used for these accessories has been tested according to EN 14126:2003 (protective clothing against infective agents) with the conclusion that the material offers a limited barrier against infective agents (see above table).

LIMITATIONS OF USE: These accessories and/or fabric are not flame resistant and should not be used around heat, open flame, sparks or in potentially flammable environments. Tyvek® melts at 135°C. It is possible that a type of exposure to bio hazards not corresponding to the tightness level of these accessories may lead to a bio-contamination of the user. Stitched seams of these accessories do not offer a barrier to infective agents. According to the European standard (EN 13034:2005 + A1:2009 (Type PB[6-B])) for chemical protective clothing Category III: these accessories offer partial body protection and have not been tested according to the whole suit low level spray test (EN ISO 17491-4, Method A). Exposure to certain very fine particles, intensive liquid sprays and splashes of hazardous substances may require full body personal protective equipment (e.g. coverall) of higher mechanical strength and barrier properties than those offered by these accessories. The user must ensure suitable reagent to garment compatibility before use. For enhanced protection in certain applications, whole suit chemical protective garments underneath with taping of cuffs, ankles, separate hood and zipper flap will be necessary. The user shall verify that tight taping is possible in case the application would require doing so. For best results in taping, the pieces of tape should overlap themselves. Care shall be taken when applying the tape, that no creases appear in the fabric or tape since those could act as channels. Ensure, that the ties on the overboots, made of these fabrics are securely tied and do not constitute a tripping hazard. Despite the slip retardant sole on the overboots, care should be taken to avoid slipping,

especially on wet surfaces. Ensure the overboots provide adequate mechanical resistance for the surface to be walked upon and that the sole is not damaged. The sole of the overboots is not liquid tight. The overboots, offering limited protection in spray form (Type PB [6]), are unsuitable for walking or standing in liquid spills and must be worn in conjunction with appropriate chemical resistant footwear. These accessories do not comply with standard EN 1149-5 (surface resistance) and are unsuitable for use in explosive zones. Please ensure that you have chosen the accessories suitable for your job. For advice, please contact your supplier or DuPont. The user shall perform a risk analysis upon which he shall base his choice of PPE. He shall be the sole judge for the correct combination of the partial body protective accessories and ancillary equipment (gloves, boots, respiratory protective equipment etc.) and for how long these accessories can be worn on a specific job with respect to their protective performance, wear comfort or heat stress. DuPont shall not accept any responsibility whatsoever for improper use of these accessories.

PREPARING FOR USE: In the unlikely event of defects, do not wear the accessory.

STORAGE AND TRANSPORT: These accessories may be stored between 15°C and 25°C in the dark (cardboard box) with no UV light exposure. DuPont has performed naturally and accelerated ageing tests with the conclusion that this fabric retains adequate physical strength over a period of 10 years. Product shall be transported and stored in its original packaging.

DISPOSAL: These accessories can be incinerated or buried in a controlled landfill without harming the environment. Disposal of contaminated garments is regulated by national or local laws.

DECLARATION OF CONFORMITY: Declaration(s) of conformity can be downloaded at: www.safespec.dupont.co.uk

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Additional information for other certification(s) independent of CE marking

Regulation 2016/425 on personal protective equipment as brought into UK law and amended.

For the purpose of these instructions for use, all BS EN or BS EN ISO standards are identical to the EN or EN ISO standards, including the date of publication, mentioned in the English text of these user instructions.

Manufacturer: DuPont de Nemours (Luxembourg) s.à r.l.
L-2984 Luxembourg

Importer of record: Du Pont (U.K.) Limited
Kings Court, London Road
Stevenage, Hertfordshire
United Kingdom, SG1 2NG

Approved Body address: SGS United Kingdom Limited
Rossmoor Business Park
Ellesmere Port, South Wirral
Cheshire, CH65 3EN

**UK
CA 0120**

Eurasian Conformity (EAC) - Complies with Technical Regulations of the Customs Union TR TS 019/2011.

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2 | Manufactured by DuPont de Nemours (Luxembourg) s.à r.l. L-2984 Luxembourg *DuPont registered trademark Ref.: XXX_XXX Made in XXX Произведено в XXX | 6