EBA 280**EBA 280 S**

DE	Bedienungsanleitung	13
EN	Operating Manual.....	51
FR	Mode d'emploi	90
IT	Istruzioni per l'uso	129

Andreas Hettich GmbH & Co. KG
Föhrenstraße 12, D-78532 Tuttlingen / Germany
Phone +49 (0)7461 / 705-0
Fax +49 (0)7461 / 705-1125
info@hettichlab.com, service@hettichlab.com
www.hettichlab.com

Single Registration Number:

DE-MF-000010680

© 2014 by Andreas Hettich GmbH & Co. KG

All rights reserved. No part of this publication may be reproduced without the prior written permission of the copyright owner.

Änderungen vorbehalten! , Modifications reserved! , Sous réserve de modifications ! , Con riserva di modifiche!

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We hereby declare under our responsibility without involvement of a notified body that the designated device:

Type of device	Small centrifuges
Name	EBA 280
	EBA 280 S
Basic UDI-DI	4050674010008QE
GMDN	36465
Classification	in vitro diagnostic, class A (Annex VIII, Rule 5)
according to	Regulation (EU) 2017/746 Annex IX

and its accessories, which are listed in the related technical documentation and whose conformity has been assessed together with the device, complies with the relevant provisions of the Regulation (EU) 2017/746 on in vitro diagnostic devices.

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The device also complies to the applicable provisions of the following European directives, ordinances and standards

- 2006/42/EC "Directive on machinery"
- 2014/30/EU "EMC Directive"
- 2014/35/EU "Low Voltage Directive"
- 2011/65/EC "RoHS Directive"
(without involvement of a notified body)
- (EC) 1907/2006 "Regulation on REACH"
(without involvement of a notified body)

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This declaration of conformity is valid from 01.01.2023 until 31.12.2023

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Basic UDI-DI **4050674010008QE**

GMDN **36465**

Classification **in vitro diagnostic, class A
(Annex VIII, Rule 5)**

according to **Regulation (EU) 2017/746
Annex IX**

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The centrifuge may only be used by qualified personnel in closed laboratories. The centrifuge is only intended for the use referred to above.

Intended use also includes observing all instructions in the Operating Manual and compliance with the required inspection and maintenance work.

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- 2006/42/EC "Directive on machinery"
- 2014/30/EU "EMC Directive"
- 2014/35/EU „Low Voltage Directive“
- 2011/65/EC "RoHS Directive"
(without involvement of a notified body)
- (EC) 1907/2006 „Regulation on REACH“
(without involvement of a notified body)

Applied relevant standards

DIN EN ISO 13485:2016; DIN EN ISO 14971:2020;
DIN EN 61010-1:2011-07; DIN EN 61010-2-020:2017;
DIN EN 61010-2-101:2002; DIN EN 61010-2-011:2017

Tuttlingen, 09.12.2022

Klaus-Günter Eberle
Gérant, Chief Executive Officer

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- **USA:** QSR, 21CFR 820 "CFR Title 21 - Food and Drugs: TITLE 21- FOOD AND DRUGS, CHAPTER I - FOOD AND DRUG ADMINISTRATION DEPARTMENT OF HEALTH AND HUMAN SERVICES, SUBCHAPTER H - MEDICAL DEVICES, Part 820 QUALITY SYSTEM REGULATIONS"
- **Kanada:** CMDR, SOR/98-282 "Medical Devices Regulations"

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Standards and regulations which apply to this device

The device is a high-end technical product. It is subject to extensive testing and certification procedures according to the following standards and regulations in their respectively valid version:

Electrical and mechanical safety for design and final testing:

Standard series: IEC 61010 (conform to standards of DIN EN 61010)

- IEC 61010-1 "Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 1: General requirements" (Pollution Degree 2, Excess-voltage category II)
- IEC 61010-2-010 „Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-010: Particular requirements for laboratory equipment for the heating of Materials" (only valid for centrifuges with heating)
- IEC 61010-2-011 „Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-011: Particular requirements for refrigerating equipment" (only valid for centrifuges with cooling)
- IEC 61010-2-020 "Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 2-020: Particular requirements for laboratory centrifuges"
- IEC 61010-2-101 "Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-101: Particular requirements for in vitro diagnostic (IVD) medical equipment"

Electromagnetic Compatibility:

- EN 61326-1 "Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 1: General requirements"

Risk management:

- DIN EN ISO 14971 "Application of risk management to medical devices"

Restriction of Hazardous Substances (RoHS II):

- EN 50581 "Technical documentation for assessing electric and electronic devices with regard to the restriction of hazardous substances"

European directives applied for conformity assessment procedures:

- Regulation (EU) 2017/746 on in vitro diagnostic devices.
- Directive 2011/65/EU for the restriction of use of certain hazardous substances in electric and electronic devices. Carrying out the EC conformity assessment process is the sole responsibility of the manufacturer, without the involvement of a notified body.

Applied medical device regulations outside Europe:

- **USA:** QSR, 21CFR 820 "CFR Title 21 - Food and Drugs: TITLE 21- FOOD AND DRUGS, CHAPTER I - FOOD AND DRUG ADMINISTRATION DEPARTMENT OF HEALTH AND HUMAN SERVICES, SUBCHAPTER H - MEDICAL DEVICES, Part 820 QUALITY SYSTEM REGULATIONS"
- **Canada:** CMDR, SOR/98-282 "Medical Devices Regulations"

Certified quality management system according to

- ISO 9001 "Quality management systems – Requirements"
- ISO13485 "Medical devices - Quality management systems - Requirements for regulatory purposes"

Environmental management system according to

- ISO 14001 "Environmental management systems - Requirements with guidance for use"

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Restrictions relatives aux substances dangereuses (RoHS II):

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- **USA** : QSR, 21CFR 820 "CFR Title 21 - Food and Drugs : TITLE 21- FOOD AND DRUGS, CHAPTER I - FOOD AND DRUG ADMINISTRATION DEPARTMENT OF HEALTH AND HUMAN SERVICES, SUBCHAPTER H - MEDICAL DEVICES, Part 820 QUALITY SYSTEM REGULATIONS"
- **Canada**: CMDR, SOR/98-282 "Medical Devices Regulations"

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Contents

1	How to use these operating instructions	53
2	Symbol meanings	53
3	Intended use	54
4	Remaining risks	54
5	Technical specifications	55
6	Notes on safety	56
7	Transport and storage	57
7.1	Transport	57
7.2	Storage	57
8	Scope of delivery	57
9	Removing the transport securing device	58
10	Initial operation	58
11	Opening and closing the lid	59
11.1	Opening the lid	59
11.2	Closing the lid	59
12	Inserting and removing the rotor	60
13	Loading the rotor	61
14	Operating and display elements	62
14.1	Displayed symbols	62
14.2	Control panel keys	62
14.3	Setting options	63
15	Programming	64
15.1	Inputting or changing programs	64
15.2	Calling up programs	66
15.3	Saving programs	66
16	Centrifugation	67
16.1	Centrifugation with preset time	67
16.2	Continuous run	68
16.3	Short-term centrifugation	69
17	Emergency Stop	69
18	Rotor detection	70
19	Cycle counter	71
20	Settings and queries	72
20.1	Querying system information	73
20.2	Acoustic signal	74
20.3	Optical signal after ending the centrifugation run	75
20.4	Automatic unlocking of the lid after the centrifugation run	76
20.5	Backlighting of the display	77
20.6	Displayed centrifugation parameters after switching on the device	78
20.7	Activating/deactivating the "Dual time" function	79
20.8	Querying the hours of operation and the number of centrifugation runs	80
20.9	Resetting the cycle counter to zero	81
21	Relative centrifugal force (RCF)	82

22	Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm ³	82
23	Emergency unlocking	82
24	Maintenance and servicing.....	83
24.1	Centrifuge (housing, lid and centrifuging chamber).....	83
24.1.1	Surface cleaning and care	83
24.1.2	Surface disinfection	83
24.1.3	Removal of radioactive contaminants	83
24.2	Rotors and Attachments	84
24.2.1	Cleaning and care.....	84
24.2.2	Disinfection	84
24.2.3	Removal of radioactive contaminants	84
24.2.4	Rotors and accessories with limited service lives	84
24.3	Autoclaving	85
24.4	Centrifuge containers.....	85
25	Faults	86
26	Change mains input fuses.....	88
27	Returning Devices	88
28	Disposal	89
29	Anhang / Appendix	168
29.1	Rotoren und Zubehör / Rotors and accessories.....	168

1 How to use these operating instructions

- Before using the centrifuge, read the operating instructions and observe them.
- The operating instructions are a part of the device. They must always be kept readily available.
- If the device is set up at a different location, the operating instructions must be provided with it.

2 Symbol meanings

Symbol on the device:
Attention, general hazard area.

Symbol on the device:
Observe operating instructions.
This symbol indicates that the user must observe the operating instructions provided.

Symbol in this document:
Attention, general hazard area.
This symbol refers to safety relevant warnings and indicates possibly dangerous situations.
The non-adherence to these warnings can lead to material damage and injury to personal.

Symbol on the device and in this document:
Beware of biohazard.

Symbol in this document:
This symbol refers to important circumstances.

Symbol on the device:
The symbol is on the housing, near the hinge.
Direction of rotation of the rotor.

Symbol on the device:
The symbol is located on the front of the housing, below the borehole.
The direction of rotation of the Allen key to unlock the lid in an emergency.

Symbol on the device and in this document:
Symbol for the separate collection of electric and electronic devices according to the guideline 2012/19/EU.
Applies in the countries of the European Union, as well as in Norway and Switzerland.

Symbol on the shipping carton label:
This way up.

Symbol on the shipping carton label:
The shipping packaging must be transported and handled within the indicated humidity range (10% - 80%).

Symbol on the shipping carton label:
The shipping packaging must be stored, transported and handled within the temperature range shown (-20°C - +60°C).

Symbol on the shipping carton label:
The shipping packaging must be kept away from rain and kept in a dry environment.

Symbol on the shipping carton label:
Fragile, handle with care.

Symbol on the shipping carton label:

Stack limit. Maximum number of identical packages which may be stacked on the bottom package, whereby "n" stands for the number of permissible packages. The bottom package is not included in "n".

3 Intended use

The centrifuge **EBA 280 / EBA 280 S** is an in vitro diagnostic medical device according to the In Vitro Diagnostic Medical Devices Regulation (EU) 2017/746.

The device is used for centrifuging and enriching sample material of human origin for subsequent further processing for diagnostic purposes. The user can set each of the variable physical parameters within the limits set by the device.

The centrifuge may only be used by qualified personnel in closed laboratories. The centrifuge is only intended for the use referred to above. Intended use also includes observing all instructions in the Operating Manual and compliance with the required inspection and maintenance work.

Any other use or use beyond this is considered improper. Andreas Hettich GmbH & Co. KG shall not be liable for any damage arising from this.

4 Remaining risks

The device is built according to the state-of-the-art and the recognized safety regulations. If used and handled improperly, there could be life-threatening danger to the user or third parties, or the device could be impaired or there could be other property damage. The device is only to be used for its intended purpose and only when it is in safe working condition.

Malfunctions which could affect safety must be corrected immediately.

5 Technical specifications

Manufacturer	Andreas Hettich GmbH & Co. KG D-78532 Tuttlingen					
Model	EBA 280		EBA 280 S			
Basic-UDI-DI	4050674010008QE					
Type	1101	1101-01	1102	1102-01		
Mains voltage ($\pm 10\%$)	200 - 240 V 1~	100 - 127 V 1~	200 - 240 V 1~	100 - 127 V 1~		
Mains frequency	50 - 60 Hz	50 - 60 Hz	50 - 60 Hz	50 - 60 Hz		
Connected load	185 VA	185 VA	330 VA	330 VA		
Current consumption	0.85 A	1.75 A	1.6 A	3.0 A		
Max. capacity	6 x 50 ml					
Allowed density	1.2 kg/dm ³					
Speed (RPM)	6000					
Force (RCF)	4146		5071			
Kinetic energy	2700 Nm					
Obligatory inspection (DGUV Regel 100 - 500)	no					
Ambient conditions (EN / IEC 61010-1)	Indoors only Up to 2000 m above sea level 2°C to 40°C Maximum relative humidity 80% for temperatures up to 31°C, linearly decreasing to 50% relative humidity at 40°C. II 2					
Device protection class	I					
Not suitable for use in explosion-endangered areas.						
EMC	– Emitted interference, Interference immunity	EN / IEC 61326-1, Class B	FCC Class B	EN / IEC 61326-1, klasse B	FCC Class B	
Noise level (dependent on rotor)		≤ 51 dB(A)		≤ 56 dB(A)		
Dimensions	– Width 326 mm – Depth 400 mm – Height 242 mm					
Weight	approx. 11 kg					

6 Notes on safety

 **No claim of warranty will be considered by the manufacturer unless ALL instructions in this manual have been followed.**

 **Reports of serious incidents involving the device**

Report any serious incidents involving the device to the manufacturer and, if necessary, to the competent authority.

-  • The centrifuge should be installed on a good, stable base.
 • Before using the centrifuge, make sure you check that the rotor is engaged on the motor shaft.
 • When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.
 • Rotors, suspensions and accessories that possess traces of corrosion or mechanical damage or if their term of use has expired may not be used any longer.
 • The centrifuge may no longer be put into operation when the centrifuging chamber has safety-related damages.
 • For centrifuges without temperature control, when the room temperature is increased and/or if the device is frequently used, the centrifuging chamber could be heated up. Therefore, it can't be ruled out that the sample material might be changed due to the temperature.

- Before the initial operation of your centrifuge you should read and pay attention to the operating instructions. Only personnel that has read and understood the operating instructions are allowed to operate the device.
- Along with the operating instructions and the legal regulations on accident prevention, you should also follow the recognised professional regulations for working in a safe and professional manner. These operating instructions should be read in conjunction with any other instructions concerning accident prevention and environmental protection based on the national regulations of the country where the device is to be used.
- This centrifuge is a state-of-the-art piece of equipment which is extremely safe to operate. However, it can lead to danger for users or others if used by untrained staff, in an inappropriate way or for a purpose other than that it was designed for.
- The centrifuge must not be moved or knocked during operation.
- In case of fault or emergency release, never touch the rotor before it has stopped turning.
- To avoid damage due to condensate, when changing from a cold to a warm room the centrifuge must either heat up for at least 3 hours in the warm room before being connected to the mains, or run hot for 30 minutes in the cold room.
- Only the rotors and accessories approved by the manufacturer for this device may be used (see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories"). Before centrifuge vessels are used which are not listed in the chapter "Appendix, Rotors and accessories", the user must make sure they can be used by asking the manufacturer.
- The centrifuge rotor may only be loaded in accordance with the chapter "Loading the rotor".
- When centrifuging with maxim revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm³.
- The centrifuge may only be operated when the balance is within the bounds of acceptability.
- The centrifuge may not be operated in explosion-endangered areas.
- The centrifuge must not be used with:
 - inflammable or explosive materials
 - materials that react with one another producing a lot of energy.
- If users have to centrifuge hazardous materials or compounds contaminated with toxic, radioactive or pathogenic micro-organisms, they must take appropriate measures.

For hazardous substances centrifuge containers with special screw caps must strictly be used. In addition to the screw cap centrifuge containers, for materials in hazard category 3 and 4 a biosafety system must be used (see the World Health Organisation's "Laboratory Biosafety Manual").

Without the use of a biosafety system the centrifuge is not microbiologically sealed in the sense of the EN / IEC 610101-2-020 standard.

No biosafety systems are available for this centrifuge.

- The centrifuge must not be operated with highly corrosive substances which could impair the mechanical integrity of rotors, hangers and accessories.
- Repairs must only be carried out by personnel authorised to do so by the manufacturer.
- Only original spare parts and original accessories licensed by the Andreas Hettich GmbH & Co. KG company are allowed to be utilised.
- The following safety regulations apply:
EN / IEC 61010-1 EN / IEC 61010-2-020 as well as their national deviations.
- The safe operation and reliability of the centrifuge can only be guaranteed if:
 - the centrifuge is operated in accordance with the operating instructions,
 - the electrical installation on the site where the centrifuge is installed conforms to the demands of EN / IEC stipulations.
- Meeting the country-specific requirements concerning occupational safety with regard to the use of laboratory centrifuges at the workplaces provided for this purpose by the user is the responsibility of the user.

7 Transport and storage

7.1 Transport

Before transporting the device, the transport securing device must be installed.

When the device and accessories are transported, the following ambient conditions must be complied with:

- Ambient temperature: -20°C to +60°C
- Relative humidity: 10% to 80%, non-condensing

7.2 Storage

The device and the accessories may only be stored in closed and dry rooms.

When the device and accessories are stored, the following ambient conditions must be complied with:

- Ambient temperature: -20°C to +60°C
- Relative humidity: 10% to 80%, non-condensing

8 Scope of delivery

The following accessories are delivered with the centrifuge:


- 1 connecting cable
- 2 fuses
- 1 hex head wrench
- 1 operating manual
- 1 information sheet, transport securing device
- 1 information sheet, emergency unlocking

The rotor(s) and associated accessories are included in the delivery in the quantity ordered.

9 Removing the transport securing device

! It is imperative that the transport securing device be removed.
 Keep the transport securing device in a safe place since it must be installed again before transporting the device.
 The device may only be transported with the transport securing device installed.

To protect the device from damage during transport, the motor is fixed in place.
 This transport securing device must be removed when the device is put into operation.

- Remove the two screws (a) and spacer sleeves (b).

! The transport securing device is installed in the opposite order.

10 Initial operation

- Remove the transport securing device from the housing floor; see the "Removing the transport securing device" chapter.
- Position the centrifuge in a stable and level manner in a suitable place. During set-up, the required safety margin of 300 mm around the centrifuge is to be kept according to EN / IEC 61010-2-020.**

When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.

- Ventilation openings may not be blocked.
 A distance of 300 mm must be maintained from the ventilation slots and openings of the centrifuge.
- Check whether the mains voltage tallies with the statement on the type plate.
- Connect the centrifuge with the power cord to a standard mains socket. For connection ratings refer to Chapter "Technical specifications".
- Switch on the mains switch (switch position "I"). The LEDs will briefly light up.
 The following displays appear one after the other:
 1. The centrifuge model
 2. The type number and program version
 3. The rotor code (R) and the maximum speed of the rotor (maxRPM) of the last rotor detected by the rotor detector
 4. The centrifugation parameters of the last used program or program 1

If the lid is closed, the message "Open the lid" is displayed.

In this case, open the lid to display the centrifugation parameters.

11 Opening and closing the lid

11.1 Opening the lid

The lid can only be opened if the centrifuge is switched on and the rotor is stationary.
If this shouldn't be possible, see the "Emergency unlock" chapter.

When the cycle counter is activated, after a centrifugation run, while opening the lid, the remaining number of running cycles (centrifugation runs) is briefly displayed.

Example:

PROG | >RCF< | RPM | t/min:s
RemCycles = 16703

- Press the key.
 - The lid is unlocked by the motor.
- : Lid unlocked.
The left LED in the key will go out.

Example:

PROG | >RCF< | RPM | t/min:s
1 6 4500 5:00

11.2 Closing the lid

Do not reach with your fingers between the lid and housing.
Do not slam the lid closed.

- Put the lid on and lightly press the front edge of the lid down.
 - The lid is locked by the motor.
- : Lid locked.


The left LED in the key is illuminated.

Example:

PROG | >RCF< | RPM | t/min:s
1 5 4500 5:00

12 Inserting and removing the rotor

Rotor with unlocking knob

Inserting the rotor:


Dirt particles between the motor shaft and rotor prevent the rotor from having a perfect seat and cause it to run unsteadily.
If necessary, clean the motor shaft (B) and the rotor hub before inserting the rotor.

- Lift up the rotor by the unlocking knob (A) and place horizontally on the motor shaft (B). The rotor will automatically engage on the motor shaft.
- Hold the rotor firmly at the left and right and pull up lightly to check whether it has properly engaged on the motor shaft.

Removing the rotor:

- Pull up the unlocking knob (A) on the rotor, hold onto it and lift the rotor from the motor shaft.

Rotor without unlocking knob

Inserting the rotor:

Dirt particles between the motor shaft and rotor prevent the rotor from having a perfect seat and cause it to run unsteadily.
If necessary, clean the motor shaft (B) and the rotor hub before inserting the rotor.

- Lift up the rotor by the handle (C), place it horizontally on the motor shaft (B) and press down as far as it can go.

Removing the rotor:


- Hold the rotor firmly by the handle (C) and lift it off of the motor shaft.

13 Loading the rotor


Standard centrifuge containers of glass will not stand RCF values exceeding 4000 (DIN 58970, pg. 2).

- Hold the rotor firmly at the left and right and pull up lightly to check whether it has properly engaged on the motor shaft.
- With swing-out rotors all rotor positions must be lined with **identical** hangers.
- The rotors and hangers may only be loaded symmetrically. The centrifuge containers have to be distributed evenly on all rotor positions. For authorised combinations see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".
In the case of angle rotors all possible rotor positions must be loaded, see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

Rotor is evenly loaded


Not permitted!

Rotor is not evenly loaded


- In containers with rubber inserts, the same number of rubber inserts must always be among the centrifuge containers.
- The centrifuge containers may only be filled outside of the centrifuge.
- The maximum filling quantity for the centrifuge containers specified by the manufacturer must not be exceeded.

In the case of angle rotors, the centrifuging vessels may only be filled so far that no fluid can be expelled from them while the centrifuge is running.

- When loading the angle rotors, no liquid may get into the angle rotors or centrifuging chamber.
- When loading the hangers and swing-out rotors and when swivelling out the hangers while the centrifuge is running, no liquid may enter the hangers or the centrifuging chamber.
- In order to maintain the weight differences within the centrifuge container as marginal as possible, a consistent fill level in the containers is to be heeded.
- The weight of the permissible filling quantity is specified on each rotor. This weight may not be exceeded.

14 Operating and display elements

14.1 Displayed symbols

- Lid unlocked.
- Lid locked.
- Rotation display. The rotation display lights up, rotating counter clockwise, as long as the rotor is turning.

14.2 Control panel keys

- Calling up programs
 - Save program.
-
-
- Switching between the RPM display (RPM) and RCF display (>RCF<). RCF values are displayed in parentheses ><. The LED above the key will light up if the RCF display is activated.
- RPM Speed
RCF : Relative centrifugal acceleration

- Selection key for activating the individual parameters. Every time the key is pressed, the next parameter is activated.
- Call up the "MACHINE MENU" (keep the key pressed for 8 seconds).
- In the "Machine Menu", select the menus "→ Info", "→ Settings" and "→ Time & Cycles".
- Scroll forward in the menus.

- Change the value of a previously selected parameter.
If the key is kept pressed, the value changes with increasing speed.

- Start the centrifugation run.
The LED in the key will be illuminated during the centrifugation run as long as the rotor is turning.
- Short-term centrifugation.
Centrifugation is run as long as the key is kept pressed. The LED in the key will be illuminated during the centrifugation run as long as the rotor is turning.
- Save inputs and changes.
- Call up the menus "→ Info", "→ Settings" and "→ Time & Cycles".
- Finish the centrifugation run.
The rotor runs down with a pre-selected brake stage. The right LED in the button is illuminated until the rotor is still. After the rotor stops, the left LED in the key will light up.
Pressing the key twice triggers the Emergency Stop.
- Unlock the lid.
The left LED in the key will go out.

14.3 Setting options

PROG RCL Program location of the program to be called

The program locations 1 to 9 and # can be called.

t/min Runtime. Adjustable from 1 - 99 min, in steps of 1 minute.

t/sec Runtime. Adjustable from 1 - 59 s, in steps of 1 second.

Continuous run "---". Set the parameters **t/min** and **t/sec** to zero.

Begins at Counting the runtime.

START The runtime begins counting immediately after the centrifugation run begins.

Begins at The runtime only begins counting after reaching the set speed.

 This setting is only possible when the "Dual time" function is activated. See the "Activating/deactivating the 'Dual time' function" chapter. This function is activated ex works.

RPM Speed

A number value from 200 rpm to the maximum speed of the rotor can be set.

Settable in increments of 10.

For the maximum rotor speed, see the chapter "Appendix, Rotors and accessories".

RAD/mm Centrifuging radius

Adjustable from 10 mm to 330 mm, in steps of 1 mm.

For the centrifuging radius, see the chapter "Appendix, Rotors and accessories".

 It is only possible to input the centrifuging radius if the RCF display (>RCF<) is activated.

RCF Relative centrifugal acceleration

A number value can be set which results in a speed between 200 rpm and the maximum rotor speed.

Settable in increments of 1.

 It is only possible to input the relative centrifugal acceleration (RCF) if the RCF display (>RCF<) is activated.

The relative centrifugal acceleration (RCF) depends on the centrifuging radius (RAD). The relative centrifugal acceleration RCF and the set centrifuging radius (RAD) are displayed.

Starting steps 1 - 9. Step 9 = shortest starting time, Step 1 = longest starting time.

Brake steps 0 - 9. Step 9 = shortest run-down time, Step 1 = long run-down time, Step 0 = longest run-down time (brakeless run-down).

PROG STO Program location where the program is stored.

9 different programs can be stored (program locations 1 to 9).

 No programs should be saved at program location #.

The program location # serves as intermediate memory for changed centrifugation parameters.

The centrifugation parameters at program location # are overwritten every time a centrifugation run is carried out with changed centrifugation parameters which were not saved to a program location.

15 Programming

 No programs should be saved at program location #.

The program location # serves as intermediate memory for changed centrifugation parameters.

If centrifugation parameters are changed, but are not saved to a program location afterwards, the display shows a dash "-" instead of the program location number. After starting the centrifugation run, these centrifugation parameters are then automatically saved at program location #.

The centrifugation parameters at program location # are overwritten every time a centrifugation run is carried out with changed centrifugation parameters which were not saved to a program location.

15.1 Inputting or changing programs

 In order to set the continuous run, the parameters **t/min** and **t/sec** must be set to zero. In the time display (t/min:s), "----" appears.

Example:

The relative centrifugal force (RCF) depends on the centrifuging radius (RAD). A correct calculation is only possible if the right centrifuging radius is set. During the input of the RCF, the set centrifuging radius is displayed.

Example:

It is only possible to display and input the relative centrifugal force (RCF) and the centrifuging radius (RAD) if the RCF display (>RCF<) is activated.

If no key is pressed for 8 seconds after selection or during parameter input, the previous values are shown on the display. The parameters must then be input again.

 Entering parameters can be cancelled at any time by pressing the **STOP OPEN** key. In this case, the settings are not stored.

- Press the key.
PROG RCL : Program location of the displayed centrifugation parameters.

- With the keys, set the program location whose centrifugation parameters are to be displayed.

- Press the key.
The centrifugation parameters are displayed.

Example:
RPM display (RPM)

Example:
RCF display (>RCF<)

Continued on next page

- Press the key to activate the RPM display (RPM) or the RCF display (>RCF<) as needed. The LED above the key is illuminated if the RCF display (>RCF<) is activated.

- Keep pressing the key until the desired centrifugation parameter is displayed.
Example: t/min : Runtime, minutes.

- Set the desired value with the keys.

- Input all the desired centrifugation parameters in succession.

- Press the key until the following is displayed:
PROG STO : Program location where the centrifugation parameters are stored.

- Set the desired program location with the keys.

Example:
RPM display (RPM)

PROG	>RCF<	RPM	t/min:s
1	a	4500	5:00

Example:
RCF display (>RCF<)

PROG	>RCF<	RPM	t/min:s
1	a >	2853<	5:00

PROG	>RCF<	RPM	t/min:s
	t/min	=05:00	

PROG	>RCF<	RPM	t/min:s
	t/min	=05:00	

PROG	>RCF<	RPM	t/min:s
	t/min	=04:00	

PROG	>RCF<	RPM	t/min:s
	t/min	=04:00	

PROG	>RCF<	RPM	t/min:s
	t/sec	=04:30	

PROG	>RCF<	RPM	t/min:s
	t/sec	=04:30	

PROG	>RCF<	RPM	t/min:s
	PROG STO	= 1	

PROG	>RCF<	RPM	t/min:s
	PROG STO	= 1	

PROG	>RCF<	RPM	t/min:s
	PROG STO	= 3	

PROG	>RCF<	RPM	t/min:s
	PROG STO	= 3	

PROG	>RCF<	RPM	t/min:s
	*** OK ***		

PROG	>RCF<	RPM	t/min:s
	*** OK ***		

PROG	>RCF<	RPM	t/min:s
3	a	4500	4:30

PROG	>RCF<	RPM	t/min:s
3	a >	2853<	4:30

15.2 Calling up programs

- Press the key.
PROG RCL : Program location of the displayed centrifugation parameters.

- With the keys, set the program location whose centrifugation parameters are to be displayed.

- Press the key.
The centrifugation parameters are displayed.

Example:
RPM display (RPM)

PROG RCL = 2

Example:
RCF display (>RCF<)

PROG RCL = 2

PROG RCL = 1

PROG RCL = 1

*** OK ***

*** OK ***

1 5 4500 5:00

1 > 2853< 5:00

15.3 Saving programs

- Press the key twice.
PROG STO : Program location where the centrifugation parameters are stored.

- With the keys, set the program location where the centrifugation parameters are to be stored.

- Press the key.
The centrifugation parameters are stored at the program location.

Example:
RPM display (RPM)

PROG STO = 3

Example:
RCF display (>RCF<)

PROG STO = 3

PROG STO = 1

PROG STO = 1

*** OK ***

*** OK ***

1 5 4500 5:00

1 > 2853< 5:00

16 Centrifugation

In accordance with EN / IEC 61010-2-020, no persons, hazardous substances or objects may be within a safety range 300 mm around the centrifuge during a centrifugation run.

After each centrifugation run, while the lid is unlocking, the remaining number of running cycles (centrifugation runs) will be briefly displayed.

Example:

PROG	>RCF<	RPM	t/min:s
RemCycles= 16703			

If the permissible weight difference within the rotor load is exceeded, the centrifugation run is aborted during the start-up and the following display appears:

PROG	>RCF<	RPM	t/min:s
IMBALANCE			

STOP

OPEN

A centrifugation run can be aborted at any time by pressing the **STOP** key.

During a centrifugation run, centrifugation parameters can be selected and changed. The changed values are stored at program location #.

With the **RCF** key, you can toggle between the RPM display (RPM) and RCF display (>RCF<) at any time. If the RCF display (>RCF<) is worked with, the centrifuging radius must be input.

The following is displayed:

PROG	>RCF<	RPM	t/min:s
Open the lid			

Further operation of the centrifuge is then only possible after opening the lid.

Operation errors and malfunctions will be shown (see the chapter on "Malfunctions").

- Plug in the mains switch. Switch position I.
- Load the rotor and close the centrifuge lid.

16.1 Centrifugation with preset time

- Press the key to activate the RPM display (RPM) or the RCF display (>RCF<) as needed. The LED above the key is illuminated if the RCF display (>RCF<) is activated.

- Input the desired centrifugation parameters (see the "Programming" chapter).

- Press the key.
The centrifugation run is started.
During the centrifugation run, the rotor speed or the resulting RCF value and remaining time are displayed.

Continued on next page

Example:
RPM display (RPM)

PROG	>RCF<	RPM	t/min:s
1	■	4500	5:00

Example:
RCF display (>RCF<)

PROG	>RCF<	RPM	t/min:s
1	■ >	2853	5:00

PROG	>RCF<	RPM	t/min:s
1	■	4000	4:30

PROG	>RCF<	RPM	t/min:s
1	■ >	2254	4:30

PROG	>RCF<	RPM	t/min:s
1	■	4000	3:15

PROG	>RCF<	RPM	t/min:s
1	■ >	2254	3:15

- After the time elapses or if the centrifugation run is aborted by pressing the  key, the run-out occurs with the set brake stage. The brake stage is displayed. Example **~_7**

Example:
RPM display (RPM)

Example:
RCF display (>RCF<)

16.2 Continuous run

- Press the key to activate the RPM display (RPM) or the RCF display (>RCF<) as needed. The LED above the key is illuminated if the RCF display (>RCF<) is activated.

↓

- Input the desired centrifugation parameters. Set the parameters **t/min** and **t/sec** to zero (see the "Programming" chapter).

↓

- Press the key.
The centrifugation run is started.
During the centrifugation run, the rotor speed or the resulting RCF value and elapsed time are displayed.

Example:
RPM display (RPM)

Example:
RCF display (>RCF<)

- Press the key.
The centrifugation run is ended.
The run-out occurs with the set brake stage. The brake stage is displayed.
Example **~_7**

16.3 Short-term centrifugation

- Press the key to activate the RPM display (RPM) or the RCF display (>RCF<) as needed. The LED above the key is illuminated if the RCF display (>RCF<) is activated.

- Input the desired centrifugation parameters (see the "Programming" chapter).

- Press the key and keep it pressed. During the centrifugation run, the rotor speed or the resulting RCF value and elapsed time are displayed.

- Release the key again. The centrifugation run is ended. The run-out occurs with the set brake stage. The brake stage is displayed. Example --_7 .

Example:
RPM display (RPM)

PROG	>RCF<	RPM	t/min:s
1	■	4500	5:00

Example:
RCF display (>RCF<)

PROG	>RCF<	RPM	t/min:s
1	■ >	2853<	5:00

PROG	>RCF<	RPM	t/min:s
1	■	4000	5:00

PROG	>RCF<	RPM	t/min:s
1	■ >	2254<	5:00

PROG	>RCF<	RPM	t/min:s
1	■	4000	0:24

PROG	>RCF<	RPM	t/min:s
1	■ >	2254<	0:24

PROG	>RCF<	RPM	t/min:s
1	■	3980	--_7

PROG	>RCF<	RPM	t/min:s
1	■ >	2231<	--_7

17 Emergency Stop

- Press the key twice. During the Emergency Stop, the run-out occurs at brake stage 9 (shortest run-out time). Brake stage --_9 is displayed.

Example:
RPM display (RPM)

PROG	>RCF<	RPM	t/min:s
1	■	4270	--_9

Example:
RCF display (>RCF<)

PROG	>RCF<	RPM	t/min:s
1	■ >	2568<	--_9

18 Rotor detection

The centrifuge is equipped with automatic rotor detection.

Rotor detection is carried out after starting each centrifugation run.

If the rotor code of a rotor is read in for the first time by the rotor detection, the centrifugation run then stops automatically after rotor detection. The rotor code (R) and the maximum speed (maxRPM) of the rotor permitted for this centrifuge are displayed.

Example:

After the rotor comes to a standstill (the left LED in the key is illuminated), press any key. Afterwards, the maximum permissible number of run cycles is displayed.

Example:

Now, the maximum permissible number of run cycles of the rotor/hangers must be set (see the "Cycle counter" chapter).

When using a rotor whose rotor code has already been read in by the rotor detection, after starting the centrifugation run, the maximum speed (ROTOR MAX) of the rotor permissible for this centrifuge is briefly displayed.

Example:

If the set speed is greater than the maximum speed of the rotor permissible for this centrifuge, the centrifugation run is aborted after rotor detection and the following displays alternately appear:

Example:

In this case, the speed on the display must be checked and adapted to the maximum rotor speed.

19 Cycle counter

 The period of use of certain rotors and hangers is limited to a maximum number of run cycles (centrifugation runs). The rotors and hangers in question are labeled with the maximum permissible number of run cycles.

The centrifuge is equipped with a cycle counter, which counts the running cycles (centrifugation runs) of the various rotor codes (see also the "Rotor detection" chapter).

After each centrifugation run, while the lid is unlocking, the remaining number of running cycles (centrifugation runs) will be briefly displayed.

Example:

If the maximum permissible number of running cycles has been exceeded, the following is displayed every time the centrifugation run is started, and the centrifugation run must be restarted.

 If this display appears:

the rotor or hanger in question must be immediately exchanged for a new rotor/hanger for safety reasons.

After the rotor/hanger has been exchanged, the cycle counter of the respective rotor code must be reset to zero again (see "Resetting the cycle counter to zero" chapter).

If the rotor code of a rotor is read in for the first time by the rotor detection, the centrifugation run then stops automatically after rotor detection. The rotor code (R) and the maximum speed (maxRPM) of the rotor permitted for this centrifuge are displayed.

Afterwards, the maximum permissible number of run cycles of the rotor/hanger must be set.

 In the case of swing-out rotors, the maximum permissible number of run cycles of the hangers must always be set.

The maximum permissible number of run cycles of the rotor/hanger must be set as follows.

Example:

Display after rotor detection.

R : Rotor code of the rotor.

maxRPM : The maximum speed of the rotor permitted for this centrifuge.

Press any key.

- Wait until the left LED in the  key is illuminated.
- Press any key.
The maximum permissible number of run cycles will be displayed.

Continued on next page

Example:

- Check whether the displayed maximum number of run cycles matches the specifications on the rotor/hangers.

If not, the maximum permissible number of run cycles specified on the rotor/hangers must be set.

 For rotors and hangers whose period of use is not limited to a maximum number of run cycles (centrifugation runs), zero must be set.


- Set the correct value with the keys.

- Press the key.
The maximum permissible number of run cycles is saved.

20 Settings and queries

**Query:**

- System information (section 20.1)

Setting:

- Acoustical signal (section 20.2)
- Automatic unlocking of the lid after the centrifugation run (section 20.4)
- Automatic switch-off of the background lighting (section 20.5)
- Optical signal after ending the centrifugation run (section 20.3)
- Displayed centrifugation parameters after switching on the device (section 20.6)
- Activating/deactivating the "Dual time" function (section 20.7)

Query:

- Hours of operation (section 20.8)
- Number of carried-out centrifugation runs (section 20.8)

Setting:


- Cycle counter (section 20.9)

20.1 Querying system information

The following system information can be queried:

- Centrifuge model
- The maximum speeds of the various rotor codes,
- Centrifuge program version
- Frequency converter type
- Frequency converter program version

With the rotor at a standstill, the query can proceed as follows:

20.2 Acoustic signal

The acoustic signal sounds:

- after a malfunction occurs in a 2-second interval.
- after ending the centrifugation run and the rotor is at a standstill in a 30-second interval.

The acoustic signal is ended by pressing any key.

If the rotor is at a standstill, the acoustic signal can be set as follows:

The procedure can be aborted at any time by pressing the key. In this case, the settings are not stored.

- Press and hold the button for eight seconds.

* MACHINE MENU *

- Press the key until the following is displayed:

-> Settings

- Press the key.
End beep : Acoustic signal after ending the centrifugation run.
off : Signal deactivated.
on : Signal activated.

End beep =off

End beep =on

- With the keys, set **off** or **on**.

End beep =off

End beep =on

- Press the key.
Error beep : Acoustic signal after a malfunction occurs.
off : Signal deactivated.
on : Signal activated.

Error beep =off

Error beep =on

- With the keys, set **off** or **on**.

Error beep =off

Error beep =on

- Press the key.
Beep volume : Volume of the acoustic signal.
min : soft
mid : medium
max : loud

Beep volume =min

Beep volume =max

- With the keys, set **min**, **mid**, or **max**.

Beep volume =min

Beep volume =max

Beep volume =mid

Continued on next page

- Press the key to save the setting.

- Press the key once to exit the "→ Settings" menu or press twice to exit the "* MACHINE MENU *".


Example:

20.3 Optical signal after ending the centrifugation run

The backlighting of the display flashes after the centrifugation run to visually signalize that the centrifugation run has finished.

The optical signal can be switched on or off when the rotor is at a standstill:

- Press and hold the button for eight seconds.

- Press the key until the following is displayed:

- Press the key.

- Press the key until one of the following is displayed:
End blinking : Backlighting of display flashes after the centrifugation run.
off : Backlighting does not flash.
on : Backlighting flashes.

- With the keys, set **off** or **on**.

- Press the key to save the setting.

- Press the key once to exit the "→ Settings" menu or press twice to exit the "* MACHINE MENU *".

Example:

20.4 Automatic unlocking of the lid after the centrifugation run

It can be set whether the lid should automatically unlock or not after the centrifugation run.

With the rotor at a standstill, this can be set as follows:

 The procedure can be aborted at any time by pressing the  key. In this case, the settings are not stored.

- Press and hold the button for eight seconds.

* MACHINE MENU *

- Press the key until the following is displayed:

-> Settings

- Press the key.

End beep =off

End beep =on

- Press the key until one of the following is displayed:
Lid AutoOpen : Automatic unlocking of the lid after the centrifugation run.
off : Lid does not unlock automatically.
on : Lid unlocks automatically.

Lid AutoOpen=off

Lid AutoOpen=on

- With the keys, set **off** or **on**.

Lid AutoOpen=off

Lid AutoOpen=on

- Press the key to save the setting.

Store setting...

-> Settings

- Press the key once to exit the "-> Settings" menu or press twice to exit the "* MACHINE MENU *".

* MACHINE MENU *

Example:

1 6 4500 5:00

20.5 Backlighting of the display

To save energy, it can be set that, after a centrifugation run, the backlighting of the display switches off after 2 minutes.

With the rotor at a standstill, this can be set as follows:

 The procedure can be aborted at any time by pressing the  key. In this case, the settings are not stored.

- Press and hold the button for eight seconds.

- Press the key until the following is displayed:

- Press the key.

- Press the key until one of the following is displayed:
Power save : Automatic shutdown of the backlighting.
off : Automatic shutdown deactivated.
on : Automatic shutdown activated.

- With the keys, set **off** or **on**.

- Press the key to save the setting.

- Press the key once to exit the "**-> Settings**" menu or press twice to exit the "*** MACHINE MENU ***".

20.6 Displayed centrifugation parameters after switching on the device


It can be set whether the centrifugation parameters of the last used program or of program 1 should be displayed after switching on the device.

With the rotor at a standstill, this can be set as follows:

 The procedure can be aborted at any time by pressing the  key. In this case, the settings are not stored.

- Press and hold the button for eight seconds.

- Press the key until the following is displayed:

- Press the key.

- Press the key until one of the following is displayed:
Start Pr : Displayed program after switching on.
First : Program 1.
Last : Last used program.

- Set **First** or **Last** with the keys.

- Press the key to save the setting.

- Press the key once to exit the "-> Settings" menu or press twice to exit the "* MACHINE MENU *".


Example:

20.7 Activating/deactivating the "Dual time" function

In order to be able to set the runtime counter (parameter **Begins at START** and parameter **Begins at SPEED**), the "Dual time" function must be activated. This function is activated ex works.

With the rotor at a standstill, this can be set as follows:

- Press and hold the button for eight seconds.

- Press the key until the following is displayed:

- Press the key.

- Press the key until one of the following is displayed:
Dual time : Counting the runtime.
on : Function activated.
off : Function deactivated.

- Set with the **off** or **on** keys.

- Press the key to save the setting.

- Press the key once to exit the "-> Settings" menu or press twice to exit the "* MACHINE MENU *".


20.8 Querying the hours of operation and the number of centrifugation runs

The hours of operation are divided into internal and external hours of operation.

Internal hours of operation: Total time the device was switched on.

External hours of operation: Total time of the previous centrifugation runs.

With the rotor at a standstill, the query can proceed as follows:

- Press and hold the button for eight seconds.

- Press the key until the following is displayed:

- Press the key.
TimeExt : External hours of operation.

- Press the key.
TimeInt : Internal hours of operation.

- Press the key.
Starts : Number of all centrifugation runs.

- Press the key twice to exit the "-> Time & Cycles" menu or press three times to exit the "* MACHINE MENU *".

Example:

20.9 Resetting the cycle counter to zero

After the rotor/hanger has been exchanged, the cycle counter of the respective rotor code must be reset to zero again.

The cycle counter may only be reset to zero if the rotor/hanger has been exchanged for a new rotor/hanger first.

With the rotor at a standstill, the cycle counter can be reset as follows:

The procedure can be aborted at any time by pressing the key. In this case, the settings are not stored.

- Press and hold the button for eight seconds.

- Press the key until the following is displayed:

- Press the key.

- Press the key until the following is displayed:
Cyc sum : Number of completed run cycles.

- Press the key.

- Press the key.
The number of completed run cycles is reset to zero.

- Press the key to save the setting.

- Press the key twice to exit the "-> Time & Cycles" menu or press three times to exit the "* MACHINE MENU *".

Example:

21 Relative centrifugal force (RCF)

The relative centrifugal force (RCF) is given as a multiple of the acceleration of gravity (g). It is a unit-free value and serves to compare the separation and sedimentation performance.

These values are calculated using the formula below:

$$\text{RCF} = \left(\frac{\text{RPM}}{1000} \right)^2 \times r \times 1,118 \quad \Rightarrow \quad \text{RPM} = \sqrt{\frac{\text{RCF}}{r \times 1,118}} \times 1000$$

RCF = relative centrifugal force

RPM = rotational speed (revolutions per minute)

r = centrifugal radius in mm = distance from the centre of the turning axis to the bottom of the centrifuge container. For more on the centrifugal radius see the chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

The relative centrifugal force (RCF) stands in relation to the revolutions per minute and the centrifugal radius.

22 Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm³

When centrifuging with maxim revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm³.

The speed must be reduced for materials or mixtures of materials with a higher density.

The permissible speed can be calculated using the following formula:

$$\text{Reduced speed (n}_{\text{red}}\text{)} = \sqrt{\frac{1.2}{\text{Greater density [kg/dm}^3\text{]}}} \times \text{maximum speed [RPM]}$$

e.g.: maximum speed RPM 4000, density 1.6 kg/dm³

$$\text{n}_{\text{red}} = \sqrt{\frac{1.2 \text{ kg/dm}^3}{1.6 \text{ kg/dm}^3}} \times 4000 \text{ RPM} = 3464 \text{ RPM}$$

If in doubt you should obtain clarification from the manufacturer.

23 Emergency unlocking

In the event of a power failure, the lid cannot be unlocked with the motor. Emergency unlocking must be done by hand.


To unlock in an emergency, disconnect the centrifuge from the mains.
Open the lid only when the rotor is at a standstill.

CAUTION! Damage to the lock during emergency unlocking by turning the hexagon Allen key in clockwise direction (to the right).

The Allen key may only be turned counter clockwise (to the left); see figure.

- Switch off the mains switch (switch setting "0").
- Look through the window in the lid to make sure that the rotor is at a standstill.
- Insert the Allen key horizontally in the bore (A) and turn carefully counter clockwise (to the left) until the lid opens.
- Pull the Allen key back out of the bore.

24 Maintenance and servicing

The device can be contaminated.

Pull the mains plug before cleaning.

Before any other cleaning or decontamination process other than that recommended by the manufacturer is applied, the user has to check with the manufacturer that the planned process does not damage the device.

- Centrifuges, rotors and accessories must not be cleaned in rinsing machines.
- They may only be cleaned by hand and disinfected with liquids.
- The water temperature must be between 20 – 25°C.
- Only detergents/disinfectants may be used which:
 - have a pH between 5 - 8
 - do not contain caustic alkalis, peroxides, chlorine compounds, acids and alkaline solutions
- In order to prevent appearances of corrosion through cleaning agents or disinfectants, the application guide from the manufacturer of the cleaning agent or disinfectant are absolutely to be heeded.

24.1 Centrifuge (housing, lid and centrifuging chamber)

24.1.1 Surface cleaning and care

- Clean the centrifuge housing and the centrifuging chamber regularly, using soap or a mild detergent and a damp cloth if required. For one thing, this serves purposes of hygiene, and it also prevents corrosion through adhering impurities.
- Ingredients of suitable detergents:
soap, anionic tensides, non-ionic tensides.
- After using detergents, remove the detergent residue by wiping with a damp cloth.
- The surfaces must be dried immediately after cleaning.
- In the event of condensation water formation, dry the centrifugal chamber by wiping out with an absorbent cloth.
- Lightly rub the rubber seal of the centrifuge chamber with a rubber care product after each cleaning.
- The centrifuging chamber is to be checked for damage once a year.

If damage is found which is relevant to safety, the centrifuge may no longer be put into operation. In this case, notify Customer Service.

24.1.2 Surface disinfection

- If infectious materials penetrates into the centrifugal chamber this is to be disinfected immediately.
- Ingredients of suitable disinfectants:
ethanol, n-propanol, ethyl hexanol, anionic tensides, corrosion inhibitors.
- After using disinfectants, remove the disinfectant residue by wiping with a damp cloth.
- The surfaces must be dried immediately after disinfecting.

24.1.3 Removal of radioactive contaminants

- The agent must be specifically labelled as being an agent for removing radioactive contaminants.
- Ingredients of suitable agents for removing radioactive contaminants:
anionic tensides, non-ionic tensides, polyhydrated ethanol.
- After removing the radioactive contaminants, remove the agent residue by wiping with a damp cloth.
- The surfaces must be dried directly after removing the radioactive contaminants.

24.2 Rotors and Attachments

24.2.1 Cleaning and care

- In order to avoid corrosion and changes in materials, the rotors and accessories have to be cleaned regularly with soap or with a mild cleaning agent and a moist cloth. Cleaning is recommended at least once a week. Contaminants must be removed immediately.
- Ingredients of suitable detergents:
soap, anionic tensides, non-ionic tensides.
- After using detergents, remove detergent residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after cleaning.
- Angle rotors, container and hanger made of aluminium are to be lightly greased after drying using acid-free grease, e.g. vaseline.
- The rotors and accessories have to be checked weekly for wear and corrosion.

Rotors and attachments may no longer be utilised upon indication of wear and tear or corrosion.

24.2.2 Disinfection

- If infectious material should get on the rotors or accessories, they must be appropriately disinfected.
- Ingredients of suitable disinfectants:
ethanol, n-propanol, ethyl hexanol, anionic tensides, corrosion inhibitors.
- After using disinfectants, remove disinfectant residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after disinfection.

24.2.3 Removal of radioactive contaminants

- The agent must be specifically labelled as being an agent for the removal of radioactive contaminants.
- Ingredients of suitable agents for removing radioactive contaminants:
anionic tensides, non-ionic tensides, polyhydrated ethanol.
- After removing the radioactive contaminants, remove agent residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after removing the radioactive contaminants.

24.2.4 Rotors and accessories with limited service lives

The use of certain rotors, hangers and accessory parts is limited by time.

These are marked with the maximum permitted number of operating cycles or with an expiration date and the maximum permitted number of operating cycles or just with the expiration date; e.g.:

- usable until end of month/year: 10/2024 ("usable until end of month/year: 10/2024")
- "max. cycles: 40000" (max. run cycles: 40000)

For safety reasons, rotors, hangers and accessory parts may no longer be used if either the indicated maximum number of operating cycles or the indicated expiration date has been reached.

The centrifuge is equipped with a cycle counter, which counts the running cycles (centrifugation runs) of the various rotor codes. For a description, see the "Cycle counter" chapter.

24.3 Autoclaving

The following accessory can be autoclaved at 121°C / 250°F (20 min):

- Swing-out rotors without hangers
- Reductions

No statement can be made about the degree of sterility.

Autoclaving accelerates the ageing process of plastics. In addition, autoclaving may discolour plastics.

24.4 Centrifuge containers

- With leakiness or after the breakage of centrifuging containers broken container parts, glass splinters and leaked centrifugation material are to be completely removed.
- The rubber inserts as well as the plastic sleeves of the rotors are to be replaced after a glass breakage.

Remaining glass splinters cause further glass breakage!

- If this concerns infectious material, a disinfection process is to be executed immediately.

25 Faults

If the fault cannot be eliminated with the help of the fault table, please inform Customer Service.

Please specify the type of centrifuge and the serial number. Both numbers can be found on the name plate of the centrifuge.

 Perform a MAINS RESET:	<ul style="list-style-type: none"> – Switch off the mains switch (switch position "0"). – Wait at least 10 seconds and then switch on the mains switch again (switch position "I").
--	---


Message / fault		Cause	Remedy
No display		No voltage Mains input fuses defective.	<ul style="list-style-type: none"> – Check distribution voltage. – Check mains power input fuse, refer to Chapter "Change mains input fuse". – Mains switch ON.
IMBALANCE		The rotor is unevenly loaded.	<ul style="list-style-type: none"> – Open the lid after the rotor is at a standstill. – Check the loading of the rotor, see chapter "Loading the rotor". – Repeat the centrifugation run.
RPM > ROTOR MAX		Speed in the selected program greater than the maximum speed of the rotor.	<ul style="list-style-type: none"> – Check the set speed. Reduce the set speed
MAINS INTERRUPT		Power failure during the centrifugation run. (The centrifugation run was not finished.)	<ul style="list-style-type: none"> – Open the lid after the rotor is at a standstill. – Push  button. – Repeat the centrifugation run if necessary.
Rxxx WRONG ROTOR (xxx = Anzeige des Rotorcodes)		The used rotor is not approved for this device.	<ul style="list-style-type: none"> – Use a rotor which is approved for this device.
KEYBOARD ERROR		Error / defect electronics	<ul style="list-style-type: none"> – Perform a MAINS RESET when the rotor has been stationary.
TACHO ERROR	1, 2, 96.1	Failure of speed impulses during operation. No rotor inserted. Error / defect electronics	<ul style="list-style-type: none"> – After the displayed waiting time expires (150 seconds), carry out a MAINS RESET. – Check whether the rotor is inserted.
LID ERROR	4.1 – 4.127	Error in lid locking or lid closure.	<ul style="list-style-type: none"> – Perform a MAINS RESET when the rotor has been stationary.
OVER SPEED	5	Rotation too fast	
ROTOR ERROR	10.1 – 10.6	Incorrect rotor coding	
VERSION ERROR	12	Incorrect centrifuge model recognized. Error / defect electronics	
UNDER SPEED	13	Rotation too slow	

Message / fault		Cause	Remedy
CTRL ERROR	22 – 25.4	Error / defect electronics	<ul style="list-style-type: none"> – Perform a MAINS RESET when the rotor has been stationary.
CRC ERROR	27, 27.1		
COM ERROR	31 – 36		
FC ERROR	60, 61.1 – 61.21, 61.64 – 61.142		
TACHO ERR	61.22	Speed measurement error	<ul style="list-style-type: none"> – The device must not be switched off as long as the "wait" message is being displayed. After the message "wait ..." is no longer displayed, carry out a MAINS RESET.
FC ERROR	61.23		
FC ERROR	61.153	Error / defect electronics	<ul style="list-style-type: none"> – Perform a MAINS RESET. – Check the loading of the rotor, see chapter "Loading the rotor". – Repeat the centrifugation run.
VERS. ERR	61.154	Invalid machine version	<ul style="list-style-type: none"> – Perform a MAINS RESET.

26 Change mains input fuses

Switch off the mains switch and separate the device from the mains!

The fuse holder (A) with the mains input fuses is located next to the mains switch.

- Remove the connecting cable from the machine plug socket.
- Press the snap-fit (B) against the fuse holder (A) and remove.
- Exchange defective mains input fuses.

Only use fuses with the rating defined for the type. See the following table.

- Reinsert the fuse holder until the snap-fit clicks shut.
- Reconnect the device to the mains supply.

Model	Type	Fuse	Order no.
EBA 280	1101	T 3,15 AH/250V	E997
EBA 280 S	1102	T 3,15 AH/250V	E997
EBA 280	1101-01	T 6,3 AH/250V	2266
EBA 280 S	1102-01	T 6,3 AH/250V	2266

27 Returning Devices

Before returning the device, a transport securing device has to be installed.

If the device or its accessories are returned to Andreas Hettich GmbH & Co. KG, in order to provide protection for people, the environment and materials, it has to be decontaminated and cleaned before being shipped.

We reserve the right to refuse contaminated devices or accessories.

Costs incurred for cleaning and disinfection are to be charged to the customer.

We ask for your understanding in this matter.

28 Disposal

The device can be disposed of via the manufacturer.

A Return Material Authorisation (RMA) form must always be requested for a return.

If necessary, contact the Technical Service Department of the manufacturer:

Andreas Hettich GmbH & Co. KG
Föhrenstrasse 12
78532 Tuttlingen, Germany
Phone: +49 7461 705 1400
Email: service@hettichlab.com

Disposal costs may be incurred.

WARNING

Risk of pollution and contamination for people and the environment.

When disposing of the centrifuge, people and the environment may be polluted or contaminated by incorrect or improper disposal.

- Removal and disposal may be carried out only by a trained and authorized service personnel.

The device is intended for the commercial sector ("Business to Business" - B2B).

According to Directive 2012/19/EU, the devices may no longer be disposed of with household waste.

The appliances are assigned to the following groups according to the Stiftung Elektro-Altgeräte Register (EAR (German foundation under civil law)):

- Group 5 (small appliances)

The crossed-out bin symbol indicates that the device must not be disposed of with household waste.

Regulations governing disposal of such devices may differ in individual countries. If necessary, contact the supplier.

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