◀ D U P O N T ▶™

Tychem®

For greater good™

Tychem Cat. III PROTECTION LEVEL

ACCESSORIES


2 | ◀ D U P O N T ▶ | 8

1 | Tychem®

3 | Tychem 2000 C Gown PL50 Gown 0290 Apron PA30LO Sleeve PS32LA Boot cover POBA

Tychem 6000 F Gown PL50 Gown 0290 Apron PA30LO Sleeve PS32LA Boot cover POBA

C €0598 | 4

10 |  | 9

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13 | Other certification(s) independent of CE marking | 11

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DuPont™, the DuPont Oval Logo, and all trademarks and service marks denoted with ™, SM or ® are owned by affiliates of DuPont de Nemours, Inc. unless otherwise noted. © 2019 DuPont.

Internet: www.ipp.dupont.com

DuPont de Nemours (Luxembourg) s.à r.l.
L-2984 Luxembourg

CE Ref.: Tychem Accessories
July 2019/25/V2

DuPont Ref.: IFUATCF_012

Overview of the accessories

Tychem® 2000 C Gown PL50

Tychem® 6000 F Gown PL50

- Gown style apparel for front protection with closure in back
- Elasticated cuffs
- Shin length

Seams: Stitched and overtaped

Fabrics: Tychem® 2000 C
Tychem® 6000 F

Body measurements in cm

| Size | Chest girth | Body height |
|--------|-------------|-------------|
| SM/MD | 84-100 | 162-176 |
| LG/2XL | 100-124 | 174-194 |

Tychem® 2000 C Gown 0290

Tychem® 6000 F Gown 0290

- Gown style apparel for front protection with closure in back
- Collar
- Knitted cuffs
- Double elasticated cuffs
- Shin length

Seams: Stitched and overtaped

Fabrics: Tychem® 2000 C
Tychem® 6000 F

Body measurements in cm

| Size | Chest girth | Body height |
|--------|-------------|-------------|
| SM/MD | 84-100 | 162-176 |
| LG/2XL | 100-124 | 174-194 |

Tychem® 2000 C Apron PA30L0

Tychem® 6000 F Apron PA30L0

- Neck and waist ties
- Shin length

Seams: Stitched

Fabrics: Tychem® 2000 C
Tychem® 6000 F

Body measurements in cm

| Size | Chest girth (apron width) | Body height (apron length) |
|-------------|---------------------------|----------------------------|
| Single size | 84-132 (70) | 162-200 (108) |

Tychem® 2000 C Sleeve PS32LA

Tychem® 6000 F Sleeve PS32LA

- Wide elastics (at cuffs and upper arm) for comfortable fit
- 50 cm length

Seams: Stitched and overtaped

Fabrics: Tychem® 2000 C
Tychem® 6000 F

Body measurements in cm

| Size | Length | Upper arm opening diameter | Wrist opening diameter |
|--------|--------|----------------------------|------------------------|
| Single | 50 | 13-25 (elasticated) | 10-15 (elasticated) |

Tychem® 2000 C Boot cover POBA

Tychem® 6000 F Boot cover POBA

- Knee length overboot with slip retardant sole

Seams: Stitched and overtaped
Sole not liquid-tight, only stitched

Fabrics: Tychem® 2000 C
Tychem® 6000 F

Body measurements in cm

| Size | Max. shoe size that can be worn below overboots (UK Mens' shoe size) |
|-------------|----------------------------------------------------------------------|
| Single size | 45 (11) |

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ENGLISH

INSTRUCTIONS FOR USE

INSIDE LABEL MARKINGS ① Trademark. ② PPE manufacturer. ③ Model identification for accessories made of Tychem® 2000 C and 6000 F fabric. This instruction for use provides information on these accessories. ④ CE marking - The Partial Body (PB) accessories comply with the requirements for category III personal protective equipment according to European legislation, Regulation (EU) 2016/425. Type-examination and quality assurance certificates were issued by SGS Fimko Oy, P.O. Box 30 (Särkiniementie 3), 00211 HELSINKI, Finland, identified by the EC Notified Body number 0598. ⑤ Indicates compliance with European standards for chemical protective clothing. ⑥ Partial body protection "types" achieved by these accessories defined by the European standards for chemical protective clothing: EN 14605:2005 + A1 2009 (Type PB [3]). These accessories also fulfill the requirements of EN 14126:2003 Type PB [3-B]. ⑦ Country of origin. ⑧ Date of manufacture. ⑨ Flammable material. Keep away from fire. These accessories and/or fabrics are not flame resistant and should not be used around heat, open flame, sparks or in potentially flammable environments. ⑩ Sizing pictogram indicates body or article measurements (cm) & correlation to letter code. Check your body, arm or foot measurements and select the correct size. ⑪ Wearer should read these instructions for use. ⑫ Do not re-use. ⑬ Other certification(s) information independent of the CE marking and the European notified body.

PERFORMANCE OF THESE ACCESSORIES:

| FABRIC PHYSICAL PROPERTIES | | Tychem® 2000 C | | Tychem® 6000 F | |
|-----------------------------|-------------------|----------------|-----------|----------------|-----------|
| Test | Test method | Result | EN Class* | Result | EN Class* |
| Abrasion resistance | EN 530 Method 2 | > 1500 cycles | 5/6** | > 2000 cycles | 6/6** |
| Flex cracking resistance | ISO 7854 Method B | > 5000 cycles | 3/6** | > 1000 cycles | 1/6** |
| Trapezoidal tear resistance | EN ISO 9073-4 | > 10 N | 1/6 | > 20 N | 2/6 |
| Tensile strength | EN ISO 13934-1 | > 100 N | 3/6 | > 100 N | 3/6 |
| Puncture resistance | EN 863 | > 10 N | 2/6 | > 10 N | 2/6 |

* According to EN 14325:2004 ** Pressure pot

TYCHEM® 2000 C FABRIC RESISTANCE TO PENETRATION BY LIQUIDS (EN ISO 6530)

| Chemical | Penetration index - EN Class* | Repellency index - EN Class* |
|------------------------|-------------------------------|------------------------------|
| Sulphuric acid (30%) | 3/3 | 3/3 |
| Sodium hydroxide (10%) | 3/3 | 3/3 |
| o-Xylene | 3/3 | 3/3 |
| Butan-1-ol | 3/3 | 2/3 |

* According to EN 14325:2004

TYCHEM® 6000 F FABRIC RESISTANCE TO PENETRATION BY LIQUIDS (EN ISO 6530)

| Chemical | Penetration index - EN Class* | Repellency index - EN Class* |
|------------------------|-------------------------------|------------------------------|
| Sulphuric acid (30%) | 3/3 | 3/3 |
| Sodium hydroxide (10%) | 3/3 | 3/3 |
| o-Xylene | 3/3 | 3/3 |
| Butan-1-ol | 3/3 | 3/3 |

* According to EN 14325:2004

TYCHEM® 2000 C FABRIC AND TAPED SEAMS RESISTANCE TO PERMEATION BY LIQUIDS (EN ISO 6529 METHOD A, BREAKTHROUGH TIME AT 1 µg/cm²/min)

| Chemical | Breakthrough time (min) | EN Class* |
|----------------------|-------------------------|-----------|
| Sulphuric acid (98%) | > 480 | 6/6 |

* According to EN 14325:2004

TYCHEM® 6000 F FABRIC AND TAPED SEAMS RESISTANCE TO PERMEATION BY LIQUIDS (EN ISO 6529 METHOD A, BREAKTHROUGH TIME AT 1 µg/cm²/min)

| Chemical | Breakthrough time (min) | EN Class* |
|----------|-------------------------|-----------|
| n-Hexane | > 480 | 6/6 |

* According to EN 14325:2004

FABRIC RESISTANCE TO PENETRATION OF INFECTIVE AGENTS

| Test | Test method | EN Class* |
|---------------------------------------------------------------------------------|-----------------------|-----------|
| Resistance to penetration by blood and body fluids using synthetic blood | ISO 16603 | 6/6 |
| Resistance to penetration by blood-borne pathogens using bacteriophage Phi-X174 | ISO 16604 Procedure C | 6/6 |
| Resistance to penetration by contaminated liquids | EN ISO 22610 | 6/6 |
| Resistance to penetration by biologically contaminated aerosols | ISO/DIS 22611 | 3/3 |
| Resistance to penetration by biologically contaminated dust | ISO 22612 | 3/3 |

* According to EN 14126:2003

SEAM PROPERTIES

| Test method | Tychem® 2000 C | | Tychem® 6000 F | |
|--------------------------------|----------------|-----------|----------------|-----------|
| | Result | EN Class* | Result | EN Class* |
| Seam strength (EN ISO 13935-2) | > 125 N | 4 of 6 | > 125 N | 4 of 6 |

* According to EN 14325:2004

For further information about the barrier performance, please contact your supplier or DuPont: www.ipp.dupont.com

RISKS AGAINST WHICH THE PRODUCT IS DESIGNED TO PROTECT: These accessories only provide partial body protection and are intended to protect those parts of the body exposed to risk. They are designed to protect workers from hazardous substances, or sensitive products and processes from contamination by people. They may be used separately or in combination with other Personal Protective Equipment to increase the protection level needed. They are typically used, depending on chemical toxicity and exposure conditions, for protection against certain organic (Tychem® 6000 F) & inorganic (Tychem® 2000 C) liquids. Fabrics used for these accessories have passed all tests of EN 14126:2003 (protective clothing against infective agents). Under the exposure conditions as defined in EN 14126:2003 and mentioned in the table above, the obtained results conclude that the materials offer a barrier against infective agents.

LIMITATIONS OF USE: These accessories and/or fabrics are not flame resistant and should not be used around heat, open flame, sparks or in potentially flammable environments. Tyvek® melts at 135°C, the fabric coatings melt at 98°C. It is possible that a type of exposure to bio hazards not corresponding to the tightness level of these accessories may lead to a bio-contamination of the user. According to the European standard (EN 14605:2005 + A1 2009; Type PB [3-B]) for chemical protective clothing Category III: these accessories offer partial body protection and have not been tested according to the whole suit Jet test (EN ISO 17491-3:2008). Exposure to certain very fine particles, intensive liquid sprays and splashes of hazardous substances may require full body personal protective equipment (e.g. coverall) of higher mechanical strength and barrier properties than those offered by these accessories. The user must ensure suitable reagent to garment compatibility before use. For enhanced protection in certain applications, whole suit chemical protective garments underneath with taping of cuffs, ankles and hood will be necessary. The user shall verify that tight taping is possible in case the application would require doing so. For best results in taping, the pieces of tape should overlap themselves. When wearing the gown with double cuffs the wearer must ensure that the knitted cuff attached to the inner sleeve is covered at all times by garment fabric and/or appropriate gloves. Despite the double cuff, taping is required to obtain a

tight connection between glove and sleeve. Care shall be taken when applying the tape, that no creases appear in the fabric or tape since those could act as channels. Ensure, that the ties on the overboots, made of these fabrics are securely tied and do not constitute a tripping hazard. Despite the slip retardant sole on the overboots, care should be taken to avoid slipping, especially on wet surfaces. Ensure the overboots provide adequate mechanical resistance for the surface to be walked upon and that the sole is not damaged. The sole of the overboots is not liquid tight. The overboots, offering limited protection in spray form, are unsuitable for walking or standing in liquid spills and must be worn in conjunction with appropriate chemical resistant footwear. These accessories shall not be used whilst in presence of flammable or explosive atmospheres or while handling flammable or explosive substances. Please ensure that you have chosen the accessories suitable for your job. For advice, please contact your supplier or DuPont. The user shall perform a risk analysis upon which he shall base his choice of PPE. He shall be the sole judge for the correct combination of the partial body protective accessories and ancillary equipment (gloves, boots, respiratory protective equipment etc.) and for how long these accessories can be worn on a specific job with respect to its protective performance, wear comfort or heat stress. DuPont shall not accept any responsibility whatsoever for improper use of these accessories.

PREPARING FOR USE: In the unlikely event of defects, do not wear the accessory.

STORAGE AND TRANSPORT: These accessories may be stored between 15 and 25°C in the dark (cardboard box) with no UV light exposure. DuPont has performed tests according to ASTM D-572 with the conclusion that these fabrics retain adequate physical strength over a period of 10 years. Product shall be transported and stored in its original packaging.

DISPOSAL: These accessories can be incinerated or buried in a controlled landfill without harming the environment. Disposal of contaminated garments is regulated by national or local laws.

DECLARATION OF CONFORMITY: Declaration of conformity can be downloaded at: www.safespec.dupont.co.uk

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