ROTI®ImmunoBlock

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Instructions for use

ROTI®ImmunoBlock

Blocking reagent for immunofluorescent applications

A. Introduction

ROTI®ImmunoBlock is an innovative, ready to use blocking reagent made on a polymer basis without detergent and following advantages:

- Your detection system is not contaminated by any foreign proteins. The signal intensity remains constant.
- Specific binding sites are preserved whilst nonspecific reactions are suppressed. Signal intensity increased.
- The ready to use solution saves you time and labour by eliminating the need for dissolving.

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe
P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0
Fax: +49 (0) 721/ 5606-149
info@carlroth.com • www.carlroth.com ip 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

B. Application

ROTI®Block working solution

Dilute ROTI®ImmunoBlock (10 x concentrate) before use 1:10 with H₂O (distilled). The diluted solution can be used directly.

ROTI®ImmunoBlock and samples

When blocking please cover your entire sample with ROTI®ImmunoBlock working solution.

Blocking time and temperature

Blocking should occur under gentle shaking at room temperature. A blocking period from 15 mins to 1 h has proved successful.

ROTI®ImmunoBlock and antibodies

You can also incubate your sample with primary and secondary antibody diluted in ROTI®ImmunoBlock working solution. Alternatively incubation can occur in ROTI®ImmunoBlock working solution diluted 1:10 in PBS or TBS, respectively, or PBST or TBST, respectively.¹⁴

TBST: ROTI®Stock 10 x TBST (Art. No. 1061.1)
(25 mM Tris, pH 7.4, 150 mM NaCl, 3 mM KCl, 0.05% Tween 20)

PBST: ROTI®Stock 10 x PBST (Art. No. 1059.1)
(137 mM NaCl, 2.7 mM KCl, 10 mM Na₂HPO₄, 2 mM KH₂PO₄, pH 7.4, 0.05 % Tween 20)

TBS: ROTI®Stock 10 x TBS (Art. No. 1060.1)
(25 mM Tris, pH 7.4, 150 mM NaCl, 3 mM KCl)

PBS: ROTI®Stock 10 x PBS (Art. No. 1058.1)
(137 mM NaCl, 2.7 mM KCl, 10 mM Na₂HPO₄, 2 mM KH₂PO₄, pH 7.4)

ROTI®ImmunoBlock	T144.1	250 ml
-------------------------	--------	--------