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Instructions for use

100 bp DNA-ladder, equalised

This marker was generated from different plasmides¹ by *EcoRI*-restriction. After restriction the DNA was deproteinised with phenol / chloroform, desalted, spectroscopically measured in TE-buffer and dry frozen.

¹ One mutagenesis per plasmid is legally protected. The reproduction or use of this mutagenesis requires our written consent.

Ideal for analysis of fragments ≤500 bp.

The DNA-quantity per band was set in such a way that the intensity of the individual bands is comparable (equalised).

Delivery includes:

- T833.1 - 20 µg 100 bp-DNA-ladder, equalised DNA
- 0100.1 - 1 x sample buffer ROTI®Load with glycerol (Art. No. 0100.1)

The DNA-ladder is supplied in lyophilised form.

Fragment sizes ² (in bp)	DNA-mass with 0.25 µg gel load (in ng)
1000	22
900	20
800	19
700	20
600	22
500 (3x)	58
400	18
300	20
200	22
150	15
100	14

² Taking the sequencing accuracy as a basis, we can guarantee an error rate of <1 % when indicating fragment sizes.

Application:

Depending on its application, the DNA-ladder can be absorbed directly in the supplied sterile-filtered 1x sample buffer ROTI®Load with glycerol (Art. No. 0100.1) or in TE-buffer (ROTI®Stock 100x TE, Art. No. 1052) or in sterile bidistilled water. In this case the dry frozen DNA is dissolved in the appropriate volume under occasional shaking for 10 minutes at room temperature. You will obtain a concentration of 0.05 µg/µl by dissolving the marker in 400 µl sample buffer – optimal for most applications.

Sample application:

The standard loading of MINI to MIDI gels per track is:

- After ethidium bromide staining and when visible under UV-light: 0.2 - 0.5 µg
- With signal-integrated camera systems: 0.05 - 0.3 µg

Storage:

Optimal storage temperature is -20 °C. Repeated (>10 times) thawing and freezing will damage the marker and should be avoided. We recommend storing the marker aliquoted in portions.

Please note:

The 100 bp DNA-ladder, equalised, is highly compatible with the 1 kbp DNA-ladder (Y014). Mixed 1:1 (f. e. 0.4 µg 100 bp DNA-ladder and 0.4 µg 1 kbp DNA-ladder) they result in a broad and reliable range of DNA-marker bands for all standard applications in the laboratory.


Fig.: 1.8 % agarose (NEEO, Art. No. 2267)

100 bp DNA-ladder, equalized

- T833.1 20 µg DNA + sample buffer
- T833.2 4 x 20 µg DNA + sample buffer

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe
P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0
Fax: +49 (0) 721/ 5606-149
info@carlroth.com • www.carlroth.com ip 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.