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## Instructions for use

### 100 bp DNA-ladder, equimolar

This marker was generated from different plasmids<sup>1</sup> by *EcoRI*-restriction. After restriction the DNA was deproteinised with phenol / chloroform, desalted, spectroscopically measured in TE-buffer and dry frozen.

<sup>1</sup> One mutagenesis per plasmid is legally protected. The reproduction or use of this mutagenesis requires our written consent.

#### Routine marker for small fragments.

The individual bands of the marker are present in the mixture in the same molar concentration (equimolar).

#### Delivery includes:

T834.1 - 50 µg 100 bp-DNA-ladder, equimolar DNA

0100.1 - 1 x sample buffer ROTI®Load with glycerol (Art. No. 0100.1)

The DNA-ladder is supplied in lyophilised form.

#### Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe  
P.O. Box 100121 • 76231 Karlsruhe  
Phone: +49 (0) 721/ 5606-0  
Fax: +49 (0) 721/ 5606-149  
info@carlroth.com • www.carlroth.com ip 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

Fragment sizes <sup>2</sup> (in bp)	DNA-mass with 1 µg gel load (in ng)
1000	163
900	146
800	130
700	114
600	97
500 (2x)	163
400	65
300	49
200	33
150	24
100	16

<sup>2</sup> Taking the sequencing accuracy as a basis, we can guarantee an error rate of <1 % when indicating fragment sizes.

#### Application:

Depending on its application, the DNA-ladder can be absorbed directly in the supplied sterile-filtered 1x sample buffer ROTI®Load with glycerol (Art. No. 0100.1) or in TE-buffer (ROTI®Stock 100x TE, Art. No. 1052) or in sterile bidistilled water. In this case the dry frozen DNA is dissolved in the appropriate volume under occasional shaking for 10 minutes at room temperature. You will obtain a concentration of 0.1 µg/µl by dissolving the marker in 500 µl sample buffer – optimal for most applications.

#### Sample application:

The standard loading of MINI to MIDI gels per track is:


- After ethidium bromide staining and when visible under UV-light: 0.3 - 0.8 µg
- With signal-integrated camera systems: 0.05 - 0.3 µg

#### Storage:

Optimal storage temperature is -20 °C. Repeated (>10 times) thawing and freezing will damage the marker and should be avoided. We recommend to store the marker aliquoted in portions.

#### Please note:

The 100 bp DNA-ladder, equimolar, is highly compatible with the 1 kbp DNA-ladder (Y014). Mixed 1:1 (f. e. 0.4 µg 100 bp DNA-ladder and 0.4 µg 1 kbp DNA-ladder) they result in a broad and reliable range of DNA-marker bands for all standard applications in the laboratory.

#### 100 bp DNA-ladder, equimolar

T834.1 50 µg DNA + sample buffer

T834.2 4 x 50 µg DNA + sample buffer