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Instructions for use

100 bp DNA-ladder, extended

This marker was generated from seven plasmids¹ by *EcoRI*-restriction. After restriction the DNA was deproteinised with phenol / chloroform, desalted, spectroscopically measured in TE-buffer and dry frozen.

¹ One mutagenesis per plasmid is legally protected. The reproduction or use of this mutagenesis requires our written consent.

Routine marker for small and medium sized fragments.

Fragment sizes 100 to 1000 bp are equimolar.

Delivery includes:

T835.1 - 50 µg 100 bp-DNA-ladder, extended DNA

0100.1 - 1 x sample buffer ROTI®Load with glycerol (Art. No. 0100.1)

The DNA-ladder is supplied in lyophilised form.

Fragment sizes ² (in bp)	DNA-mass with 1 µg gel load (in ng)
5000	111
4000	89
3000	67
2500	56
2000	44
1500	33
1000	98
900	88
800	78
700	68
600	58
500 (2x)	98
400	39
300	29
200	20
150	14
100	10

(40 % DNA-mass between 5000-1500 bp and 60 % DNA-mass between 1000-100 bp, equimolarly mixed in both ranges).

² Taking the sequencing accuracy as a basis, we can guarantee an error rate of <1 % when indicating fragment sizes.

Application:

Depending on its application, the DNA-ladder can be absorbed directly in the supplied sterile-filtered 1x sample buffer ROTI®Load with glycerol (Art. No. 0100.1) or in TE-buffer (ROTI®Stock 100x TE, Art. No. 1052) or in sterile bidistilled water. In this case the dry frozen DNA is dissolved in the appropriate volume under occasional shaking for 10 minutes at room temperature. You will obtain a concentration of 0.1 µg/µl by dissolving the marker in 500 µl sample buffer – optimal for most applications.


Sample application:

The standard loading of MINI to MIDI gels per track is:

- After ethidium bromide staining and when visible under UV-light: 0.5 - 0.8 µg
- With signal-integrated camera systems: 0.2 - 0.5 µg

Storage:

Optimal storage temperature is -20 °C. Repeated (>10 times) thawing and freezing will damage the marker and should be avoided. We recommend storing the marker aliquoted in portions.

100 bp DNA-ladder, extended

T835.1 50 µg DNA + sample buffer

T835.2 4 x 50 µg DNA + sample buffer

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe
P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0
Fax: +49 (0) 721/ 5606-149
info@carlroth.com • www.carlroth.com ip 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.