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PCR-Marker DNAscore

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Instructions for use

PCR-Marker DNAscore

This marker was generated by *EcoRI*-restriction of different plasmids¹. After restriction, the DNA was deproteinised, desalted, spectroscopically measured in TE-buffer and lyophilised.

¹One mutagenesis per plasmid is protected by law. Reproduction or use of this type of mutagenesis requires our written consent.

PCR-Marker DNAscore separates very quickly and is therefore ideal for analysis of PCR-fragments and for diagnostic routines. Bands 500 bp, 1000 bp and 1500 bp are intensified to enable speedy and easy orienting in the gel. The special 1000 bp band also enables the marker to be used for quantifying in gel. When loading the gel with 0.5 µg marker (5 µl marker solution, 0.1 µg/µl) the 1000 bp band contains exactly 140 ng DNA.

Delivery range:

T917.1 - 50 µg Marker DNA
0100.1 - 1 x gel loading buffer ROTI®Load DNA with glycerol

Fragment sizes² (in bp)

1500 (intensified)
1200
1100
1000 (intensified) Quantitation band
900
800
500 (intensified)


²Against the background of the sequence accuracy, we can guarantee an error rate of <1 % regarding the fragment sizes.

Application:

Depending on its use, the PCR-Marker DNAscore can be absorbed directly in the supplied sterile-filtered 1 x gel loading buffer ROTI®Load DNA with glycerol (Art. No. 0100.1) or in TE-buffer (ROTI®Stock 100 x TE, Art. No. 1052.1) or in sterile double-distilled water. In this case the lyophilized DNA is dissolved in the appropriate volume under occasional shaking for 10 minutes at room temperature. You will obtain a concentration of 0.1 µg/µl by dissolving the marker in 500 µl gel load buffer – optimal for most applications.

Storage:

Optimal storage temperature is –20 °C. Repeated (>10 times) thawing and freezing damages the marker and should be avoided. We recommend aliquotation and freezing the marker in portions.

- a) 1.0 cm: The marker shows three bands (500 bp, 1000 bp, 1500 bp).
- b) 1.5 cm: Four satellite bands flank the 1000-bp-band.
- b) 3 4 cm: All bands are resolved. Efficient size analysis of fragments with a difference of 50 bp is possible.
- b) 6-7 cm: The marker can also be used for quantifying after approx. 6 cm run: the 1000-bp-band corresponds to 140 ng DNA (when loading 0.5 µg DNA-marker).

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe
P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0
Fax: +49 (0) 721/ 5606-149
info@carlroth.com • www.carlroth.com ip 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

PCR-Marker DNAscore

T917.1 50 µg DNA + gel loading buffer