## PCR-Marker DNAscore

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# Instructions for use

## PCR-Marker DNAscore

This marker was generated by *EcoRI*-restriction of different plasmids<sup>1</sup>. After restriction, the DNA was deproteinised, desalted, spectroscopically measured in TE-buffer and lyophilised.

<sup>1</sup>One mutagenesis per plasmid is protected by law. Reproduction or use of this type of mutagenesis requires our written consent.

PCR-Marker DNAscore separates very quickly and is therefore ideal for analysis of PCR-fragments and for diagnostic routines. Bands 500 bp, 1000 bp and 1500 bp are intensified to enable speedy and easy orienting in the gel.

The special 1000 bp band also enables the marker to be used for quantifying in gel. When loading the gel with 0.5 µg marker (5 µl marker solution, 0.1 µg/µl) the 1000 bp band contains exactly 140 ng DNA.

### Delivery range:

T917.1 - 50 µg Marker DNA  
0100.1 - 1 x gel loading buffer ROTI®Load DNA with glycerol

### Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5  
76185 Karlsruhe  
Postfach 100121  
76231 Karlsruhe  
Telefon: +49 (0) 721/ 5606-0  
Telefax: +49 (0) 721/ 5606-149  
E-Mail: info@carlroth.de  
Internet: www.carlroth.de

gh 01/2020

### Fragment sizes<sup>2</sup> (in bp)

1500 (intensified)  
1200  
1100  
1000 (intensified) Quantitation band  
900  
800  
500 (intensified)


<sup>2</sup>Against the background of the sequence accuracy, we can guarantee an error rate of <1 % regarding the fragment sizes.

### Application:

Depending on its use, the PCR-Marker DNAscore can be absorbed directly in the supplied sterile-filtered 1 x gel loading buffer ROTI®Load DNA with glycerol (Art. No. 0100.1) or in TE-buffer (ROTI®Stock 100 x TE, Art. No. 1052.1) or in sterile double-distilled water. In this case the lyophilized DNA is dissolved in the appropriate volume under occasional shaking for 10 minutes at room temperature. You will obtain a concentration of 0.1 µg/µl by dissolving the marker in 500 µl gel load buffer – optimal for most applications.

### Storage:

Optimal storage temperature is –20 °C. Repeated (>10 times) thawing and freezing damages the marker and should be avoided. We recommend aliquotation and freezing the marker in portions.

- a) 1.0 cm:** The marker shows three bands (500 bp, 1000 bp, 1500 bp).  
**b) 1.5 cm:** Four satellite bands flank the 1000-bp-band.  
**b) 3.4 cm:** All bands are resolved. Efficient size analysis of fragments with a difference of 50 bp is possible.  
**b) 6-7 cm:** The marker can also be used for quantifying after approx. 6 cm run: the 1000-bp-band corresponds to 140 ng DNA (when loading 0.5 µg DNA-marker).

### PCR-Marker DNAscore

T917.1 50 µg DNA + gel loading buffer