250 bp DNA-Leiter

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Instructions for use

250 bp DNA-Ladder

This marker was generated by *EcoRI*-restriction of different plasmids. After restriction, the DNA was deproteinized, desalted, spectroscopically measured in TE-buffer, and lyophilized.

1) One mutagenesis per plasmid is protected by law. Reproduction or use of this type of mutagenesis requires our written consent.

Universal marker for size analysis of plasmids and DNA-insertions. With two distinctly intensified DNA-bands at 1500 bp and 2500 bp, the 250 bp DNA-Ladder enables speedy orienting in 1 to 1.4 % agarose gels. This enables clear and precise gel-electrophoretic size analysis of DNA-insertions in the 250 to 3000 bp range and of vector plasmids at a higher molecular range. In contrast to an equimolar distribution of bands in other DNA-markers, our 250 bp DNA-ladder has intensified lower bands (250 – 1000 bp), while the upper bands are reduced in mass and, therefore, staining intensity.

This enables not only optimum depiction with only 70 mm running distance, but also as little as 500 ng marker DNA to obtain good band signals in mini gels.

Content:

T918.1 50 µg DNA (250 bp DNA-Ladder)
0100.1 1 x ROTI®Load loading buffer (with glycerol)

Fragment sizes² (in bp)

8000
6000
5000
4000
3000
2750

2500 (intensified)

2250
2000
1750

1500 (intensified)

1250
1000
750
500
250

2) Against the background of the sequence accuracy, we can guarantee an error rate of <1 % regarding the fragment sizes.

Application:

Depending on its use, the DNA-Ladder can be dissolved directly in the supplied sterile-filtered 1x sample buffer ROTI®Load with glycerol (Art. No. 0100.1), in TE-buffer (ROTI®Stock 100 x TE, Art. No. 1052.1), or in sterile bidistilled water. The lyophilized DNA is dissolved in the appropriate volume under occasional shaking for 10 minutes at room temperature. You will obtain a concentration

of 0.1 µg/µl by dissolving the marker in 500 µl sample buffer – optimal for most applications.

Sample application/concentration:

Standard loading for mini to midi gels per lane is:

- with bands visible in UV-light after ethidium bromide staining:
0.5 – 0.7 µg
- with detection after ethidium bromide staining with signal-integrated camera systems:
0.1– 0.5 µg

Storage:

Optimal storage temperature is –20 °C. Repeated (>10 times) thawing and freezing damages the marker and should be avoided.

We recommend freezing the marker in aliquots.


Fig.: 1.8 % agarose (NEEO, Art. No. 2267)

250 bp DNA-Leiter

T918.1 50 µg DNA + sample buffer

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe
P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0
Fax: +49 (0) 721/ 5606-149
info@carlroth.com • www.carlroth.com ip 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.