
LAB

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
LAB

**N820.18 EX
TRANSLATION OF ORIGINAL
OPERATING INSTRUCTION
ENGLISH**

EN

LABOPORT®

CHEMICALLY- RESISTANT LABORATORY PUMP

Notice!

Before operating the pump and the accessories, read the operating and installation instructions and observe the safety notices!

KNF develops, produces and sells high-quality diaphragm pumps and systems for gases, vapors and fluids. They are available as built-in pumps (OEM), for process applications and for use in laboratories.

With an extensive range application-specific pumps and systems, KNF has been a technology leader for decades. Countless customers around the world trust KNF as a supplier of long-lasting and reliable products for a wide array of different applications.

Thank you for placing your trust in us and for deciding to purchase one of our products.

Quality has top priority with our products, which is why each of our products undergoes a 100% test. All processes are certified in accordance with ISO 9001, and the RoHS- and Reach guidelines are satisfied for all materials.

KNF Neuberger GmbH
Alter Weg 3
79112 Freiburg
Germany
Tel. +49 (0)7664/5909-0

E-mail: info.de@knf.com
www.knf.com

Index

1	Scope of delivery	4
2	About this document	5
3	Use	8
4	Information on the Ex-designation	16
5	Safety.....	17
6	Technical data	20
7	Design and function	23
8	Transport	26
9	Set up, installation and connection	28
10	Operation	30
11	Servicing	36
12	Troubleshooting	46
13	Spare parts and accessories	50
14	Returns	51
15	Index	52

1 Scope of delivery

- Laboport ® N820/840
- Key for hose connector (WAF 14)
- Power cable
- Silencer
- Operating instructions
- QuickStart


Fig. 1: Scope of delivery

Unpacking the pump

Inspect the pump and the included accessories for transport damage after unpacking. If the packaging is damaged, inform the responsible forwarder so that a damage report can be prepared. For further information, read Chapter 8 *Transport*.

2 About this document

2.1 Using the operating instructions

The operating instructions are part of the pump.

- In the event of uncertainties with regard to the content of the operating instructions, please contact the manufacturer (contact data: see www.knf.com). Please have the type designation and serial number of the pump ready.
- Read the operating instructions before you commission the pump.
- Give the operating instructions to the next owner.
- Keep the operating instructions within reach at all times.

Project pumps For customer-specific project pumps (pump models that begin with "PJ" or "PM"), there may be deviations from the operating and installation instructions.

- For project pumps, also observe the agreed specifications.

2.2 Exclusion of liability

The manufacturer assumes no liability for damages and malfunctions resulting from failure to observe the operating instructions.

The manufacturer assumes no liability for damages and malfunctions resulting from changes or modifications to the device and improper handling.

The manufacturer assumes no liability for damages and malfunctions resulting from impermissible spare parts and accessories.

2.3 Symbols and markings

Warning notice

A notice that warns you of danger is located here.

Possible consequences of a failure to observe the warning notice are specified here. The signal word, e.g., warning, indicates the danger level.

→ Measures for avoiding the danger and its consequences are specified here.

Danger levels

Signal word	Meaning	Consequences if not observed
DANGER	warns of immediate danger	Death or serious injury or serious damage will result.
WARNING	warns of possible danger	Death or serious injury or serious damage is possible.
CAUTION	warns of a possibly dangerous situation	Minor injuries or damage is possible.
NOTICE	Warns of possible damage	Damage is possible.

Tab. 1: Danger levels

Other notices and symbols

→ An activity to be carried out is specified here (a step).

1. The first step of an activity to be carried out is specified here.
Follow other sequentially numbered steps.

i This symbol indicates important information.

Explanations of pictograms

Pictogram	Meaning
	General warning symbol
	Warning of hot surface
	Warning of electrical voltage
	Warning of explosive materials
	Warning of poisonous substances
	Warning of hand injuries through crushing
	Observe the operating instructions
	General mandatory sign

Tab.2: Explanations of pictograms

3 Use

3.1 Proper use

The pumps are intended exclusively for transferring gases and vapors.

Responsibility of the owner

Operating parameters and conditions Only install and operate the pumps in accordance with the operating parameters and conditions described in Chapter 6 *Technical data*.

Only pumps that are fully assembled and in the condition as delivered may be operated.

Make sure that the installation location is dry and that the pump is protected against rain, splash, gushing, and drip water as well as from other contaminants.

The pump may not be used in explosive areas. The Ex-designation is valid only for the transfer chamber (media-contacting area).

Check the tightness of the connections between the pipes of the application and the pump (or the connection of the pump) at regular intervals. Leaky connections carry the risk of releasing dangerous gases and vapors from the pump system.

Requirements on the transferred medium Before transferring a medium, check whether the medium can be transferred danger-free in the specific application.

Before using a medium, check the compatibility of the media-contacting components (see 6 *Technical data*) with the medium.

Make certain that no risk of explosion arises even in extreme operating situations (temperature, pressure) and in the event of system breakdowns.

Only transfer gases that remain stable under the pressures and temperatures that arise in the pump.

Accessories Laboratory accessories or additional components that are connected to a pump must be designed for the pneumatic data of the pump (see 6 *Technical data*).

3.2 Improper use

The pumps may not be operated in potentially explosive atmospheres.

The pumps are not suitable for use below ground.

The pumps are not suitable for transferring:

- Dusts
- Liquids
- Aerosols
- Biological and microbiological substances
- Fuel
- Explosives
- Fibers
- Oxidants
- Food

The pump may not be used to simultaneously produce vacuum and operating pressure.

No operating pressure may be applied to the suction side of the pump.

The pump may not be used if reactive, explosive or otherwise dangerous mixtures could result when the gas ballast valve is open (e.g., with the medium).

3.3 Using for transferring explosive atmospheres

To transfer explosive atmospheres, only use pumps of the appropriate device category and temperature class.

These pumps have the following explosion protection designations:

Designation	Description
	Symbol for explosion-proof pumps
II	Device group (see 3.4.1 <i>Device groups</i>)
2/-G	Device category (see 3.4.2 <i>Device categories for gas</i>)
IIB + H2	Explosion groups (see 3.4.3 <i>Explosion groups</i>)
T3	Temperature class (see 3.4.4 <i>Temperature classes</i>)
Gb	Equipment protection level (See Chapter 3.4.6 <i>Equipment protection level for gas</i>)
	Special operating conditions (See Chapter Special operating conditions)
Internal atmosphere only	Special conditions (see 3.4.7 <i>Special operating conditions</i>)

Tab.3: Explosion protection designation

An ignition hazard assessment in accordance with standards DIN EN ISO 80079-36 and DIN EN ISO 80079-37 was performed for the pumps.

The explosion protection designation can also be found at the following location:

- Type plate of the pump

3.4 Explanations of the explosion protection designation

3.4.1 Device groups

- Device group I Device group I applies for devices that are used in underground plants of mines as well as their underground systems that could be endangered by methane and/or combustible dusts.
- Device group II Device group II applies for devices that are used in other areas that could be endangered by an explosive atmosphere.

3.4.2 Device categories for gas

The device category describes the frequency and the duration of the occurrence of explosive atmospheres during operation.

Device category	Description
1G	Devices of this category are designed for use in areas in which an explosive atmosphere consisting of a mixture of air and gases, vapors or mists is present constantly or for long periods of time or often.
1D	Devices of this category are designed for use in areas in which an explosive atmosphere consisting of a dust/air mixture is present constantly or for long periods of time or often.
2G	Devices of this category are designed for use in areas in which it is to be expected that an explosive atmosphere consisting of gases, vapors or mists forms occasionally.
2/2G	Devices that extract from zone 1 and are designed for use in areas in which it is to be expected that an explosive atmosphere consisting of gases, vapors or mists forms occasionally.
2/-G	Devices that extract from zone 1 but are not designed for installation in a potentially explosive atmosphere (zone).
2D	Devices of this category are designed for use in areas in which it is to be expected that an explosive atmosphere consisting of a dust/air mixture forms occasionally.
3G	Devices of this category are designed for uses in areas in which it is to be expected that an explosive atmosphere resulting from gases, vapors or mists occurs, though in all likelihood occurs only seldom and for a very short length of time.
3D	Devices of this category are designed for uses in areas in which it is to be expected that an explosive atmosphere resulting from stirred-up dust occurs, though in all likelihood occurs only seldom and for a very short length of time.

Tab.4:

3.4.3 Explosion groups

Combustible gases and vapors are classified according to explosion groups (I, IIA, IIB and IIC) and temperature classes. The following table shows the classification of the most common combustible gases and vapors.

	T1	T2	T3	T4	T5	T6
I	Methane	–	–	–	–	–
IIA	Acetone Ethane Ethyl acetate Ammonia Ethyl chloride Benzene Acetic acid Carbon monoxide Methane Methanol Methyl chloride Naphthalene Phenol Propane Toluene	i-amyl acetate n-butane n-butyl alcohol Cyclohexanone 1,2-dichloroethane Acetic anhydride	Gasoline Diesel fuel Jet fuel Heating oils n-hexane	Acetaldehyde	–	–
IIB	Town gas	Ethylene Ethyl alcohol	Hydrogen sulfide	Ethyl ether	–	–
IIC	Hydrogen	Acetylene	–	–	–	Carbon disulfide

Tab.5:

The classification of gases and vapors into groups with respect to explosion group and temperature class applies for the transferred medium.

Transferred medium The pump may only be used to transfer gases and vapors that belong to the respective explosion group and corresponding temperature class (or below) (see designation on the type plate) or that are not explosive and combustible.

Environment of the pump The pump must not be set up in potentially explosive areas. It is only suitable for transferring explosive atmospheres.

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Maximum surface temperature The maximum surface temperature is the highest temperature that is reached under the most unfavorable conditions of a surface of the pump.

Ignition temperature The maximum surface temperature of the pump must always be lower than the lowest ignition temperature of the gas or vapor/air mixture in which it is used.

Temperature class The maximum surface temperature arises from the design of the pump and is specified as temperature class.

Temperature class	Max. surface temperature [°C]	Ignition temperature [°C]
T1	450	> 450
T2	300	> 300
T3	200	> 200
T4	135	> 135
T5	100	> 100
T6	85	> 85

Tab.6:

3.4.5 Ignition protection type

Designation	Description
H	Constructional safety "c"
H	Ignition source monitoring "b"
H	Liquid immersion "k"

Tab.7:

An ignition hazard assessment in accordance with standards DIN EN ISO 80079-36 and DIN EN ISO 80079-37 was performed for the pumps. The protective goals were reached by applying ignition protection type constructional safety "c".

3.4.6 Equipment protection level for gas

The equipment protection level describes the frequency and the duration of the occurrence of explosive atmospheres in an area.

Equipment protection level	Description*	Constructional safety
Ga	Devices with very high protection level for use in potentially explosive areas. With these devices, there is no risk of ignition during normal operation or in the event of foreseeable or infrequent faults/malfunctions.	Very high
Gb	Devices with high protection level for use in potentially explosive areas in which there is no risk of ignition during normal operation or in the event of foreseeable or infrequent faults/malfunctions.	High
Gc	Device with increased protection level for use in potentially explosive areas. There is no risk of ignition during normal operation. The devices have a number of additional protection measures which ensure that, in the event of commonly foreseeable faults in the device, no danger of ignition exists.	Increased

Tab.8: *according to DIN EN ISO 80079-36

3.4.7 Special operating conditions

Designation	Description
Internal atmosphere only	Special operating conditions

Additional conditions for the diaphragm pump:

- Do not install the pump outdoors. Commissioning may only be performed with suitable weather- and corrosion-protection paneling.
- Do not install the pump in potentially explosive atmospheres. It is only suitable for transferring explosive atmospheres.

- Install up the pump so that it cannot be damaged from the outside.
- Install the pump so that it is not exposed to any UV radiation.

4 Information on the Ex-designation

This KNF pump is designated with the following device designation in accordance with the latest explosion protection directive. The designation is only valid for the compression chamber (area that comes into contact with the medium) of the pump:


Fig. 2: Ex-designation of the pumps

Vacuum pumps of category "2/-G" are designed for transferring gases, vapors or mists for which it is to be expected that an explosive atmosphere forms occasionally.

The devices are to be arranged so that they cannot be mechanically damaged from the outside.

No changes may be made to the pumps. After replacing a wear part, the original function of the pump must be tested by achieving the specified ultimate vacuum (see operating instructions, Chapter 11 *Servicing*).

5 Safety

i Observe the safety information provided in the chapters entitled Installation and connection and Operation.

The pumps are produced in accordance with the generally recognized rules of engineering, as well as the occupational health & safety and accident prevention regulations. Nevertheless, dangers can arise during their use that lead to injuries to the user or third parties or to damage to the pump or other property.

Only use the pumps in perfect technical condition, for their intended use, safely and aware of the dangers and in observation of the operating and installation instructions.

Make certain that the components that are to be connected to the pumps are designed according to the pneumatic data of the pumps.

When connecting the pumps to the electrical mains, observe the corresponding safety rules.

Personnel Make sure that only specially trained and instructed personnel or specialist personnel work on the pumps. This applies, in particular, to assembly, connection and servicing work.

Make sure that the personnel have read and understood the operating and installation instructions, particularly the Safety chapter.

Working in a safety conscious manner Observe the regulations on accident prevention and safety during all work on the pumps and during operation.

Avoid contact with the heads and housing parts of the pump, as it heats up during operation.

Do not expose any body parts to the vacuum.

Ensure that no hazards arise from flows when gas connections are open, from the effects of noise or from hot, corrosive, dangerous, and environmentally hazardous gases.

Working with hazardous media When pumping hazardous media, observe the safety regulations for the handling of these media.

Working with combustible media and explosive atmospheres

Note that the pumps are only suitable for transporting explosive atmospheres and must not be set up in potentially explosive areas.

Make certain that the temperature of the medium is always sufficiently below the ignition temperature of the medium so as to prevent ignition or explosion. This also applies for abnormal operating situations.

It is important to note here that the temperature of the medium increases when the pump compresses the medium. You should therefore ensure that the temperature of the medium also remains far enough below the ignition temperature of the medium, even when it is compressed to the maximum permissible operating pressure of the pump. The maximum permissible operating pressure of the pump is stated in the technical data (6 *Technical data*).

Where applicable, also take into account external energy sources (such as radiation sources) that could additionally heat the medium.

In case of doubt, contact KNF Customer Service.

Environmental protection

Store and dispose of the pump, as well as all replacement parts, in accordance with applicable environmental regulations. Observe both the respective national and international regulations here. This applies in particular to parts that are contaminated with toxic substances.

Dispose of packaging material that is no longer needed in an environmentally responsible manner. The packaging materials are recyclable.

Dispose of old devices in an environmentally responsible manner. Dispose of old devices via suitable collection systems. Old devices contain valuable recyclable materials.

EU/EC directives/standards

The part of the pump that comes into contact with the media complies with Directive 2014/34/EU (ATEX).

The pumps comply with directives:

- 2011/65/EU
- 2014/30/EU (EMC)
- 2006/42/EC.

The following harmonized standards are met:

- DIN EN 61326-1
- DIN EN 61010-1
- DIN EN 1012-2
- DIN EN ISO 12100

The part of the pumps that comes into contact with the media satisfies the following harmonized standards:

- DIN EN ISO 80079-36
- DIN EN 1127-1

Per IEC 664, the pumps comply with:

- Overvoltage category II
- Degree of contamination 2

Customer
service and
repairs

The pumps are maintenance-free. However, KNF still recommends inspecting the pumps periodically for obvious changes in terms of noises and vibrations.

Only have repairs to the pumps performed by the responsible KNF Customer Service.

Housings with live components may only be opened by specialist personnel.

Use only original parts from KNF during servicing work.

6 Technical data

Technical data

Pump materials

Assembly	Material
Pump head	Modified PTFE
Diaphragm	PTFE-coated
Valve	FFPM
Connection	PTFE/FFPM
Hose connector	PVDF/FFPM
Gas ballast	PTFE/FFPM

Tab.9: Pump materials

Pneumatic data

N820.3FT.29.18G

Parameter	Value
Max. permissible operating pressure [bar rel*]	0.1
Ultimate vacuum [mbar abs.]	
At min. speed:	
Gas ballast closed	≤ 6
Gas ballast open	≤17
At max. speed:	
Gas ballast closed	≤ 8
Gas ballast open	≤15
Flow rate at atm. pressure [l/min]**	
At min. speed:	10 ± 10%
At max. speed:	20 ± 10%

Tab.10: Pneumatic data for N820.3FT.29.18G

*Bar rel related to 1013 hPa

**Liters in standard state (1013 hPa, 20°C)

N840.3FT.29.18G

Parameter	Value
Max. permissible operating pressure [bar rel*]	0.1
Ultimate vacuum [mbar abs.]	
At min. speed:	
Gas ballast closed	≤ 6
Gas ballast open	≤ 17
At max. speed:	
Gas ballast closed	≤ 8
Gas ballast open	≤ 15
Delivery rate at atm. pressure [L/min]**	
At min. speed:	18 ± 10%
At max. speed:	34 ± 10%

Tab. 11: Pneumatic data for N840.3FT.29.18G

*Bar rel related to 1013 hPa

**Liters in standard state (1013 hPa, 20°C)

Pneumatic connections

Parameter	Value
Hose connection [mm] (stepped hose nipple)	ID 8 / 9.5

Tab. 12: Pneumatic connections

Electrical data

Parameter	Value N820	Value N840
Voltage [V]	100 – 240	
Frequency [Hz]	50/60	
Power consumption [W]	60	100
Max. current consumption [A]	0.66 – 0.35	1.0 – 0.6
Max. permissible mains voltage fluctuations	± 10%	

Tab. 13: Electrical data

Weight

Pump type	Weight [kg]
N820.3FT.29.18G	8.8
N840.3FT.29.18G	11.1

Tab. 14: Weight

Other parameters

Parameter	Value
Permissible ambient temperature [°C]	+ 5 to + 40
Permissible media temperature [°C]	+ 5 to + 40
Highest permissible relative air humidity of the environment	80% for temperatures to 31°C, decreasing linearly to 50% at 40°C (non-condensing).
Maximum installation altitude [m above sea level]	2000
Protection class of pump	IP30
Dimensions L x H x W [mm] N820.3FT.29.18G	259 x 220 x 163
Dimensions L x H x W [mm] N840.3FT.29.18G	289 x 240 x 177
Equipment protection	<ul style="list-style-type: none"> → Overcurrent protection → Overtemperature protection (drive) → Blocking protection (drive)

Tab. 15: Other parameters

7 Design and function

Design

- 1 Pneumatic pump inlet
- 2 Handle
- 3 Rotary/push knob for:
 - Switching the pump on and off
 - Setting the pump speed
- 4 Pneumatic pump outlet
- 5 Power switch
- 6 Status display
- 7 Interface*


Fig.3: Diaphragm pump (shown is pump N820.3FT.29.18G)

*Interface:

NOTICE

Injury to property resulting from use of non-authorized accessories

- Only ever use accessories that are authorized by KNF (see the chapter entitled *13.2 Accessories*).

Function of a diaphragm pump

- 1 Outlet valve
- 2 Inlet valve
- 3 Transfer chamber
- 4 Diaphragm
- 5 Eccentric
- 6 Connecting rod
- 7 Pump drive


Fig.4: Function of a diaphragm pump

Diaphragm pumps transfer, compress (depending on the version) and evacuate gases and vapors.

The elastic diaphragm (4) is moved up and down by the eccentric (5) and the connecting rod (6). In the downwards stroke, it aspirates the gas to be transferred via the inlet valve (2). In the upwards stroke, the diaphragm presses the medium out of the pump head via the outlet valve (1). The transfer chamber (3) is separated from the pump drive (7) by the diaphragm.

7.1 Gas ballast


Fig.5: Operating button for gas ballast

DANGER

Personal injury through poisoning or explosion and damage to the pump

- When the gas ballast valve is open, make sure that no reactive or explosive mixtures can form.
- Close the gas ballast valve if necessary.
- If inert gas is necessary, contact KNF Service.

- i** If vaporous media are transferred, the formation of condensate in the pump heads can be minimized by opening the gas ballast valve.
- i** The final vacuum that can be achieved is worse when the gas ballast valve is open (see the chapter entitled 6 *Technical data*).

8 Transport

General

CAUTION

Personal injury and/or property damage due to incorrect or improper transport of the pump

In the event of incorrect or improper transport, the pump can fall down, be damaged or injure persons.

- Use suitable auxiliary means if necessary (carrying strap, lifting gear, etc.).
- Where appropriate, wear suitable personal protective equipment (e.g., safety shoes, safety gloves).

CAUTION

Risk of injury from sharp edges on the packaging

There is a risk of injury from cutting on the sharp edges when grabbing corners or when opening the packaging.

- Where appropriate, wear suitable personal protective equipment (e.g., safety shoes, safety gloves).

- Transport the pump in the original packaging to the installation location.
- Store the original packaging of the pump (e.g., for later storage).
- Inspect the pump for transport damage after receiving it.
- Document any transport damage in writing.

→ Remove any transport safeguards on the pump prior to commissioning.

Parameter

Parameter	Value
Storage temperature [°C]	+ 5 to + 40
Transport temperature [°C]	- 10 to + 60
Permissible humidity (non-condensing) [%]	30 to 85

Tab.16: Transport parameters

NOTICE

Prior to commissioning, make sure that the pump has reached the ambient temperature (6 *Technical data*).

9 Set up, installation and connection

Only connect the pumps in accordance with the operating parameters and conditions described in Chapter 6 *Technical data*.

Observe the safety notices (see Chapter 5 *Safety*).

- Before connecting, store the pump at the installation location to allow it to reach the room temperature (no condensation may form).

Cooling air supply

WARNING

Danger of burning on hot surfaces

Hot surfaces could occur if the pump overheats.

- When installing the pump, make sure that sufficient cooling air in-feed and discharge is ensured.

Installation location

- Make sure that the installation location is dry and that the pump is protected against rain, splash, gushing, and drip water as well as from other contaminants.
- Select a secure location (flat surface) for the pump.
- Protect the pump from dust.
- Protect the pump from vibration, impact and external damage.

9.1 Connecting the pump

1. Remove the protective caps from the pneumatic connections of the pump.

Vacuum operation with silencer

1. Remove the hose connector on the pressure side with the supplied tool (see 1 *Scope of delivery*).
2. Screw the silencer or your system component to the pressure side hand tight.
3. Connect the suction line.

NOTICE

When using a silencer, KNF recommends operating the pump in a safe environment, e.g., fume cupboard, to avoid dangers posed by harmful gases.

EN

Vacuum operation without silencer

1. Connect the suction line and the pressure line.

Connected components

→ Only connect components to the pump that are designed for the pneumatic data of the pump (see Chapter 6 *Technical data*).

Pump discharge

→ When using as a vacuum pump: Safely drain the pump discharge at the pneumatic outlet of the pump.

2. Lay the suction line and the pressure line at a downward angle to prevent condensate from running into the pump.

3. Plug the plug of the power cable into a properly installed, grounded socket.

10 Operation

10.1 Pump

10.1.1 Preparing for commissioning

Before switching on the pump, ensure the following points:

	Necessary operating requirements
Pump	→ All hoses correctly connected
Pump	<ul style="list-style-type: none"> → Data of the voltage supply system consistent with the details on the type plate of the pump. → Pump outlet not closed or restricted. → When operating with gas ballast: When venting the pump through the air inlet, no explosive or poisonous mixtures can occur.

Tab. 17: Operating requirements for commissioning

10.1.2 Commissioning

WARNING

Risk of burns from hot pump parts and/or hot medium

Some pump parts may be hot during or after operation of the pump.

- Allow the pump to cool after operation.
- Take protective measures to protect against touching hot parts.

WARNING

Injury to eyes

Coming too close to the inlet/outlet of the pump may result in injury to the eyes due to the present vacuum/operating pressure.

- Do not look into the pump inlet/outlet during operation.

-
- Only operate the pump in accordance with the operating parameters and operating conditions described in Chapter 6 *Technical data*.
 - Ensure the proper use of the pump (See Chapter 3.1 *Proper use*).
 - Eliminate the possibility of improper use of the pump (see Chapter 3.2 *Improper use*).
 - Observe the safety instructions (see Chapter 5 *Safety*).

Risk of pump head bursting due to excessive pressure increase

- Do not exceed the maximum permissible operating pressure (see 6 *Technical data*).
- Monitor the pressure during operation.
- If the pressure exceeds the maximum permissible operating pressure of the pump: immediately switch off the pump and remedy the fault (see Chapter Troubleshooting).
- Only throttle or regulate the air or gas quantity on the suction line to prevent the maximum permissible operating pressure from being exceeded.
- If the air quantity or gas quantity on the pressure line is throttled or regulated, make sure that the maximum permissible operating pressure at the pump is not exceeded.
- Ensure that the pump outlet is not closed or restricted.

**DANGER**

Risk of dangerous gas mixtures during pump operation

Depending on the medium being transferred, breakage of the media-contacting components can result in a dangerous mixture if the medium mixes with the air in the compressor housing or the surroundings.

- Before using a medium, check the compatibility of the media-contacting components (see 6 *Technical data*) with the medium.

i Excessive pressure, with all of the associated hazards, can be prevented by means of a bypass line with a pressure relief valve between the pressure side and suction side of the pump. Further information is available from KNF Customer Service (contact data: see www.knf.com).

Pump standstill → Establish normal atmospheric pressure in the lines while the pump is at a standstill (relieve pump pneumatically).

10.2 Information on switching the pump on and off

Switching on the pump


Fig. 6: Switching the pump on and off with power switch


Fig. 7: Switching the pump on and off with rotary/push knob


Fig. 8: Setting the flow rate

i The pump must not be started up against operating pressure during switch-on. This also applies during operation after a brief power failure. If a pump starts up against pressure, the pump may block, thereby activating blocking protection (drive) and switching off the pump.

- Ensure that no pressure is present in the lines when switching on.
- Switch on the pump with the power switch (see Fig. 6).
- Switch on the pump by pushing the rotary/push knob (see Fig. 7).

Setting the flow rate

The pump speed can be varied with the rotary/push knob. The flow rate can be set in this way (see Fig. 8).

- The speed setting of the rotary/push knob is retained when the pump is switched off.

Switching off/decommissioning the pump

- When transferring aggressive media, flush the pump before switching off to extend the service life of the diaphragm (see Chapter 11 Servicing).
- Switch off the pump with the rotary/push knob (see *Information on switching the pump on and off*).
- Switch off the pump with the power switch (see Fig. 3/5).
- Establish normal atmospheric pressure in the lines (relieve pump pneumatically).
- Pull out power plug of the pump.

Status display (see Fig. 3/6)

- Illuminates green if the pump was switched on via the power switch.
- Illuminates blue if the pump was also switched on via the rotary/push knob.
Pump is running:
The higher the speed is set, the brighter the status display lights up.
- Illuminates red if there is a fault:

Signal duration	Error type
100% ON (continuous light)	Drive blocked
50% ON; 50% OFF	Temperature too high (drive)
90% ON, 10% OFF	Other fault

Tab.18: Fault signal via status display

For further information, see Chapter 12 *Troubleshooting*.

11 Servicing

11.1 Servicing schedule

Risk of explosion from the formation of potentially explosive atmosphere

Leaky connections can result in dangerous potentially explosive atmospheres.

- Ensure that diaphragms and valve plates/seals are installed undamaged, cleanly and correctly.
- Check the pneumatic connections of the pump for leaks.
- Work with care during service work.
- Replace defective parts immediately.

Risk of injury when not using original parts

Failure to use original parts will result in a loss of pump functionality and safety.

The validity of the CE conformity is rendered void if genuine parts are not used.

- Use only genuine spare parts from KNF when performing maintenance/repair work.
-

Component	Servicing interval
Pump	<ul style="list-style-type: none">→ Inspect the pump periodically for external damage or leakage.→ Periodically check for noticeable changes to noises and vibrations.
Diaphragm and valve plates/seals	<ul style="list-style-type: none">→ At the latest, replace when the pump flow rate decreases.

Tab. 19: Servicing schedule

11.2 Cleaning

NOTICE

During cleaning work, ensure that no fluids enter the interior of the housing.

11.2.1 Flushing the pump

Risk of explosion by flushing the pump with air

- When using the pump with explosive media, only permit specialist to flush the pump with inert gas.

Personal injury through poisoning or explosion and damage to the pump

- When flushing the pump with inert gas, ensure that the gas ballast valve is closed and that no reactive or explosive mixtures form.

- Before switching off, flush the pump with air at atmospheric conditions (ambient pressure) for about 5 minutes (if necessary for safety reasons: with an inert gas).

11.2.2 Cleaning the pump

- Only clean the pump with a damp cloth and non-flammable cleaning agents.
- If compressed air is present, blow out the parts.

11.3 Replace diaphragm, valve plates/seals and O-rings

- Requirements
- Pump switched off and power cable pulled out of socket
 - Pump cleaned and free of hazardous materials
 - Hoses removed from pneumatic pump inlet and outlet

Material and tool	Quantity	Material
	1	TORX® screwdriver T20
	1	TORX® T25 screwdriver (only for N840.3FT.29.18G)
	1	Spare part set (see Chapter Spare parts)
	1	Felt-tip pen

Tab.20:

Information on the procedure

- Always replace diaphragms, valve plates/seals, and O-rings together to maintain the performance of the pump.
- Replace the diaphragms and valve plates/seals of the individual pump heads one after the other.

WARNING

Health hazard due to dangerous substances in the pump

Depending on the medium being transferred, caustic burns or poisoning is possible.

- Wear protective equipment if necessary, e.g., protective gloves, goggles.
- Clean the pump with suitable measures.


Preparatory steps

Fig.9: Dismount cover

1. Move the handle (Fig. 9/1) into a vertical position.
2. Slacken the screws (Fig. 9/2) on the cover (Fig. 9/3).
3. Remove the cover (Fig. 9/3).
4. Slacken the cap screws (Fig. 9/4).

Removing pump head

- 1 Screw
- 2 Pressure plate
- 3 Head plate
- 4 Valve plates/seals
- 5 Intermediate plate
- 6 Diaphragm
- 7 Shim rings


Fig.10: Pump head (pump N820 shown)

i The following item numbers refer to Fig. 10 unless specified otherwise.

1. Use a felt-tipped marker to draw a continuous line on the pressure plate (2), head plate (3), and intermediate plate (5). This prevents the parts from being incorrectly mounted later on.
2. Slacken the external screws (Fig. 9/4) of the pump heads.
3. Carefully remove the pump heads.

Replacing the diaphragm

- i** With the pump, the diaphragms (6) are changed successively to ensure that the shim rings (7) are used in the same quantity as previously.
1. Press down one diaphragm (6) so that the other diaphragm is in the upper change point.
 2. Carefully turn the upper diaphragm (6) counterclockwise by hand and remove it.
- i** Make sure that the shim rings located between the diaphragm and connecting rod do not fall into the pump housing.
Remove any shim rings that may be adhering to the diaphragms and place on the corresponding connecting rod thread.
Prerequisite for ensuring the pneumatic performance of the pump is that the same number of shim rings be mounted as before.
3. Screw in the new diaphragm (6) by hand and tighten it by hand.
 4. Perform steps 1 to 4 for the second pump head.

Changing O-rings

- 8 Gas ballast valve
- 9 Connection tube
- 10 Hose connections


Fig.11: Changing O-rings

1. Use the tool supplied to unscrew the hose connections (Fig. 11/10) and the gas ballast valve (Fig. 11/8) from the head plate (3).
2. Pull the connection tube (Fig. 11/9) out of the head plate (3).
3. Replace the O-rings.

Replacing valve plates/seals

i With the pump, the valve plates/seals (4) are replaced successively.

1. Only slacken the screw (1) of one pump head.
2. Remove the pressure plate (2) and head plate (3) from the intermediate plate (5).
3. Remove the old valve plates/seals (4).
4. Carefully clean the intermediate plate (5) (if there are deposits on it).
5. Insert the new valve plates/seals (4) into the corresponding seats on the intermediate plate (5).
6. Perform steps 1 to 5 for the second pump head.
7. Dispose of the old diaphragm (6) and valve plates/seals (4) properly.

8. Insert the connection tube (Fig. 11/9) into the head plate (3).
9. Screw the gas ballast (Fig. 11/8) into the corresponding head plate (3), so that the surface faces upward.

Fitting the pump head

1. Press on the entire perimeter of the diaphragm (6). Only press on the membrane (6) on one pump head.
2. Place the intermediate plate (5) with the valve plates/seals (4) onto the mounting fixture, using the marker pen marking for alignment.
3. Place the head plate (3) onto the intermediate plate (5), aligning it with the aligning pin.
4. Place the pressure plate (2) onto the head plate, using the marker pen marking for alignment.
5. Tighten the screws, working to a cross pattern (tightening torque: N820: 4 Nm; N840: 5 Nm).
6. Tighten the screw (1) in the center of the printing plates (tightening torque: 1 Nm).
7. Perform steps 1 to 7 for the second pump head.
8. Use the tool supplied to screw in the hose connections (Fig. 11/10) hand-tight.

Final steps

Risk of explosions from leaks

- Before recommissioning the pump, check the pump heads and pneumatic connections for leaks. Leaks may lead to a risk of explosion.

1. Function test
 - Connect the suction and discharge pressure sides on the pump.
 - Electrically connect the pump.
 - Test the pump for proper function (among other things, ultimate vacuum).
 - Again disconnect the pump electrically and pneumatically.
2. Pump integrated in application
 - Connect the suction and discharge pressure sides on the pump.
 - Electrically connect the pump.
 - Test the functionality of the pump.

12 Troubleshooting

Danger to life from electric shock

- All work on the pump may only be performed by an authorized specialist.
- Disconnect the pump power supply before working on the pump.
- Check and ensure that no voltage is present.

→ Check the pump (see following tables).

Pump not delivering	
Cause	Troubleshooting
No voltage in the electrical mains.	→ Check the circuit breaker for the room and switch it on if necessary.
Overtemperature protection of the pump has tripped	<ul style="list-style-type: none"> → Disconnect the pump from the electrical mains. → Allow the pump to cool. → Determine the cause of the overheating and rectify.
Connections or lines are blocked.	<ul style="list-style-type: none"> → Check the connections and lines. → Remove the blockage.
External valve is closed or filter is clogged.	→ Check external valves and filters.
Condensate has collected in the pump head.	<ul style="list-style-type: none"> → Separate the source of the condensate from the pump. → Flush the pump with air at atmospheric pressure for a few minutes (if necessary for safety reasons: with an inert gas).
Diaphragms or valve plates/seals are worn.	→ Replace the diaphragms and the valve plates/seals (see Chapter Replacing diaphragm and valve plates).

Tab.21: Troubleshooting: Pump not delivering

Flow rate, pressure or vacuum too low	
The pump does not reach the output stated in the technical data or data sheet.	
Cause	Fault remedy
Condensate has collected in the pump head.	<ul style="list-style-type: none"> → Separate the source of the condensate from the pump. → Flush the pump with air at atmospheric pressure for a few minutes (if necessary for safety reasons: with an inert gas). → If present, open the gas ballast and flush the pump head.
There is overpressure on the pressure side and at the same time vacuum or pressure above atmospheric pressure on the suction side.	<ul style="list-style-type: none"> → Change the pneumatic conditions.
Pneumatic lines or connection parts have insufficient cross-sections or are throttled.	<ul style="list-style-type: none"> → Disconnect the pump from the system to determine the output values. → Eliminate any throttling (e.g. valve) if necessary. → Use lines or connection parts with a larger cross section if necessary.
Leaks occur at connections, lines or pump head.	<ul style="list-style-type: none"> → Ensure the correct seating of the hoses on the hose connectors. → Ensure that the connections are correctly mounted. → Replace the leaky hoses. → Eliminate the leaks.
Connections or lines are completely or partially clogged.	<ul style="list-style-type: none"> → Check the connections and lines. → Remove any parts and particles that are causing blockages.
Head parts are soiled.	<ul style="list-style-type: none"> → Clean the head components.
Diaphragms or valve plates/seals are worn.	<ul style="list-style-type: none"> → Replace the diaphragms and the valve plates/seals (see Chapter Replacing diaphragm and valve plates).

Flow rate, pressure or vacuum too low	
The pump does not reach the output stated in the technical data or data sheet.	
Cause	Fault remedy
Replaced diaphragm and valve plates/seals.	<ul style="list-style-type: none"> → Ensure that shim rings were fitted on the diaphragm thread. → Check the hoses for leaks. → If necessary, carefully tighten the outer screws of the pressure plate crosswise.
Rotary/push knob is not set to max. speed.	<ul style="list-style-type: none"> → Set the rotary/push knob to max. speed.
Gas ballast still open.	<ul style="list-style-type: none"> → Connect the gas ballast.

Tab.22: Troubleshooting: Flow rate, pressure or vacuum too low

Pump is switched on and not running; status display is not illuminated	
Cause	Fault remedy
Pump is not connected to the electrical mains.	<ul style="list-style-type: none"> → Connect the pump to the electrical mains.
No voltage in the electrical mains.	<ul style="list-style-type: none"> → Check the circuit breaker for the room and switch it on if necessary.

Tab.23: Troubleshooting: Pump is switched on and not running; status display is not illuminated

Pump is switched on and not running; status display flashes red (50% ON, 50% OFF)	
Cause	Fault remedy
Pump has overheated, overtemperature protection has tripped.	<ul style="list-style-type: none"> → Pull power cable of the pump out of the socket. → Allow the pump to cool. → Determine the cause of the overheating and rectify.

Tab.24: Troubleshooting: Pump is switched on but not running; status display is flashing red

Pump is switched on and not running; status display illuminates red (100% ON)	
Cause	Fault remedy
Drive of the pump has blocked.	<ul style="list-style-type: none"> → Pull power cable of the pump out of the socket. → Allow the pump to cool. → Determine the cause of the blockage and rectify.

Tab.25: Troubleshooting: Pump is switched on but not running; status display lights up red

Pump is switched on and not running; status display flashes red (90% ON, 10% OFF)	
Cause	Fault remedy
Other fault	<ul style="list-style-type: none"> → Pull power cable of the pump out of the socket. → Allow the pump to cool. → Contact KNF Customer Service.

Tab.26: Troubleshooting: Pump is switched on but not running; status display is flashing red

Fault cannot be rectified

If you are unable to identify any of the specified causes, contact KNF Customer Service or send the pump to KNF Customer Service (contact data: see www.knf.com).

Preparing for return

1. Flush the pump with air for a few minutes (if necessary for safety reasons: with inert gas) to free the pump head of dangerous or aggressive gases (see Chapter 11.2.1 *Flushing the pump*).
2. Remove the pump.
3. Clean the pump (see Chapter 11.2.2 *Cleaning the pump*).
4. Send the pump together with completed Health and Safety Clearance and Decontamination Form to KNF, stating the nature of the transferred medium.

13 Spare parts and accessories

i To order spare parts and accessories, please contact your KNF sales partner or KNF Customer Service (contact data: see www.knf.com).

13.1 Spare parts

Parts	Quantity
Diaphragms	2
Valve plates/seals	4
O-ring connection tube (Ø 10 x 1.8)	2
O-ring hose connection and gas ballast valve (Ø 8 x 1.8)	3

Tab.27: Spare parts

**see Chapter 11.3 Replace diaphragm, valve plates/seals and O-rings*

Spare parts set	Ordering number
N820.3FT.29.18G	317435
N840.3FT.29.18G	317436

Tab.28: Spare parts set

13.2 Accessories

Accessories	Ordering number
Silencer	322528
Hose connection with O-ring	317278
Key for hose connector	316279

Tab.29: Accessories

14 Returns

Prerequisite for repairing a pump by KNF is a completed Decontamination Form.

This is made available on the KNF website as a download.

- To find the form, select your country on the overview page (www.knf.com).

You can find the Decontamination Form in the download area.

If you have questions, please contact your sales partner (contact data: see www.knf.com).

15 Index

A	
Additional conditions	
Conditions	14
C	
Category 2/-G	16
Cleaning	37
Commissioning	31
Connecting rod	24
Customer service	19
D	
Decontamination Form	49, 51
Device category	11
Device designation	16
Device group	10
Diaphragm	24
Disposal	18
E	
Eccentric	24
Electrical data	21
Environment of the pump	12
Environmental protection	18
Equipment protection	22
Equipment protection level	14
Explosion groups	12
F	
Fault signal	35
Flow rate	20, 34
Setting the flow rate	34
H	
Health and Safety Clearance Form	49
I	
Inlet valve	24
Installation location	28
Interface	23
Internal atmosphere only	14
K	
KNF Customer Service	49
M	
Media	
Combustible media	18
Hazardous media	17
Medium	8, 12
Mounting	28
O	
Outlet valve	24
P	
Parameters	
Operating parameters	8, 31
Transport parameters	27
Personnel	17
Pneumatic connections	21
Pneumatic data	20
Power cable	38
Power switch	23, 34
Protection type	13
Pump drive	24
Pump head	41
Pump inlet	23
Pump materials	20
Pump outlet	23

R

Recycling	18
Repair	51
Rotary/push knob	23, 34

S

Silencer	4, 28
Status display	23, 35

T

Temperature	
Ambient temperature	22
Ignition temperature	13
Maximum surface temperature	13
Media temperature	22
Storage temperature	27
Transport temperature	27
Temperature class	13
Temperature classes	12
Tool	39
Transfer chamber	24
Transport	26
Troubleshooting	46

U

Ultimate vacuum	20
-----------------	----

V

Vacuum operation	28
------------------	----

W

Warning notice	6
Weight	22

KNF worldwide

You can find our local KNF partners at: www.knf.com

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