# HIRSCHMANN®

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## Instruction manual

Please read before use!

# ceramus® classic

## EM-dispenser PP

Flaschenaufsatz-Dispenser  
Bottle-top Dispenser

Hirschmann Laborgeräte GmbH & Co. KG  
Hauptstraße 7-15 · 74246 Eberstadt Germany  
Fon +49 7134 511-0 · Fax +49 7134 511-990  
[www.hirschmannlab.com](http://www.hirschmannlab.com) · [info@hirschmannlab.com](mailto:info@hirschmannlab.com)

X2124520004

## Inhalt

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## Contents

	Page
1. Warning and safety instructions .....	4
2. Individual components .....	8
3. Functional principle .....	9
4. Prohibited use .....	10
5. Restrictions on use .....	11
6. Preparing the unit .....	12
6.1 Securing the suction tube .....	12
6.2 Securing the ejection unit .....	14
6.3 Screwing the unit on to the bottle .....	15
7. Dispensing .....	16
7.1 Closure cap removal .....	17
7.2 Venting the unit .....	18
7.3 Volume configuration .....	19
7.4 Dispensing .....	20
7.5 After use, locking of dispensing mechanism .....	21
8. Cleaning .....	22
8.1 Service/Cleaning .....	23
8.2 Intensive cleaning .....	28
9. Sterilisation .....	30
9.1 Preparation .....	30
9.2 Removal of the piston .....	30
9.3 Loosening the suction/ejection valve .....	32
9.4 Sterilisation .....	33
9.5 Reassembly of the piston .....	34

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## Contents

	Page
9.6 Tightening the suction/ejection valve .....	35
9.7 Reassembly of the dispenser .....	35
10. Valve replacement .....	36
10.1 Suction valve replacement .....	36
10.2 Ejection valve replacement .....	38
11. Repair by manufacturer .....	39
12. Volumetric inspection .....	40
13. Trouble-shooting .....	44
14. Technical data/Materials .....	46
15. Accessories and spare parts.....	49
16. Clearance certificate form .....	52
17. Declaration of conformity.....	54

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---

### **Attention!**

Before use please read the safety instructions and observe them carefully. Before use every member of staff must be familiar with these safety guidelines and they must be kept available at all times. This instruction manual cannot describe all possible safety hazards. It is every user's responsibility to observe all health and safety precautions and to determine and observe restrictions before use of the unit.

---

- Observe all general safety practices and guidelines for laboratories, such as the use of protective clothing, goggles and gloves when handling hazardous chemicals.

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#### 1. Warning and safety instructions

- Observe the recommendations of the manufacturer of the chemicals.
- The device may only be used for dispensing liquids. Defined prohibitions and restrictions for use must be observed (if in doubt contact the manufacturer).
- Before every use always check the tightness of all connections and safe installation of the ejection unit and suction tube. A loose ejection unit might cause splashing when the piston is pushed down.
- Before use it is essential to remove the closure cap from the ejection tube. Never press the piston down when the closure cap is fitted to the tube.

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## 1. Warning and safety instructions


- Make sure, that neither you nor other persons are endangered. Never dispense towards any person, avoid splashing, use suitable vessels. When dispensing, always direct the liquid on to the wall of the vessel, never to its bottom or the content.
- Never use force!
- After use close the ejection tube with the closure cap. Attention! Liquid might drip off.
- After use lock the mechanism (zero position, see point 7.5). Never leave the dispenser with the piston raised (except for cleaning purposes).

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## 1. Warning and safety instructions

- Only use genuine manufacturer's accessories and spare parts.
- After sterilisation (see point 9.) let the dispenser cool down to room temperature.
- If problems occur e. g.
  - stiff or sticking piston
  - sticking or leaking suction or ejection valves
  - leaking parts
  - cracks in componentsstop dispensing immediately. Before further use observe the instructions for cleaning (see point 8.) and Troubleshooting (see point 13.). If necessary, return it to manufacturer for repair, see point 11.

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### 3. Functional principle

- This instrument is a bottle top dispenser with digital volume selection for the dispensing of liquids. The basic model plus the various adapters (see accessories) enable its use with commonly used reagent bottles.
- When the piston is raised the ball of the suction valve is lifted. Liquid is aspirated through the suction tube from the reservoir bottle into the cylinder. After finishing the suction phase the suction valve closes.
- Pushing down the piston opens the ejection valve and dispenses the volume which was selected on the volume scale.

To select the required volume the piston must be in the down position. Turn the adjusting ring in either direction, so that the required volume mark is set above the adjustment arrow.

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#### 4. Prohibited use

---

##### The unit should not be used with the following:

- Hydrofluoric acid (HF), it attacks the glass (exception: version HF).
  - Liquids that attack ECTFE, FEP, PTFE, borosilicate glass 3.3,  $Al_2O_3$  or Hastelloy (version HF: platinum-iridium).
  - Suspensions (e. g. of activated carbon solvent), because solid particles might block the valves.
  - Liquids with low ignition temperatures and explosive liquids (e.g. carbon disulfide).
-

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## 5. Restrictions on use

- Avoid electrostatic charging when using flammable liquids (e.g. only dispense in glass vessels).
- Concentrated saline solution, highly concentrated alkaline solutions may cause the piston to stick or move stiffly.
- The dispenser is suitable for operating temperatures between 4°C and 50°C.
- The unit should only be used in areas where it can be assured that no hazardous potentially-explosive atmospheres occur.
- Do not deposit in potentially-explosive areas.

### **Attention - important hint!**

Never use force if the piston sticks or become stiff. Stop use immediately and carry out a thorough cleaning (see point 8.2).

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## 6. Preparing the unit


---

### Attention, please note before assembly!

- For safe and trouble free operation only use the genuine manufacturer's accessories and parts.
- Never use damaged or other than genuine manufacturer's parts.

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
### 6.1 Securing the suction tube

- Hold the suction tube beside the storage bottle. The tube end with the union nut should face towards the base of the bottle.
- Cut the tube at a slant level to the lower edge of the locking ring on the storage bottle.

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
#### 6.1 Securing the suction tube

- Push the suction tube on up to the stop.
- Then tighten the union nut onto the valve mantle.

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

## 6.2 Securing the ejection unit

- Ejection unit with:  
tube guide  
union nut  
ejection tube  
closure cap

- Put the ejection unit on up to the stop.

- Tighten the union nut firmly and ensure that it is securely fitted.

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
### Attention!

- Ensure that the ejection unit is firmly attached.
- Only use genuine manufacturer's ejection unit (see accessories and spare parts list).
- Never use deformed or damaged ejection units.

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## 6.3 Screwing the unit onto the bottle

Screw the dispenser onto the reservoir bottle, either directly or with one of the enclosed adapters, depending on the bottle's thread size.

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**Warning!**

- Observe all safety instructions, usage prohibitions and restrictions on use.
  - Wear protective clothing, safety glasses and protective gloves.
- 

**7. Dispensing**

---

**Attention!**

- Follow safety instructions (see point 1.).
  - Usage prohibitions and restrictions of use must be observed (see point 4. and 5.).
  - Dispense only when the ejection unit is turned away from the user and any other person.
  - Never press down the piston when the closure cap is fitted to the tube.
-

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
## 7.1 Closure cap removal

- Position a collecting vessel under the ejection unit.
  - Remove closure cap as shown in the drawing.
- Caution:** Liquid might drip off!
- Slide the closure cap's ring to the direction of the device's body as shown in the drawing.

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## 7.2 Venting the unit

- Turning the scale ring away from the zero position will set the unit for operation.
- Position a collecting vessel under the ejection unit.
- Check whether closure cap is opened. (precautionary).
- Eject air bubbles by pulling the piston a short distance and firmly pressing it down.
- Repeat this operation a few times until no more air bubbles are visible in the observation window (see arrows).
- Start operation (small volume) until no more air bubbles are visible in the ejection tube unit (see arrows).

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
### 7.2 Venting the unit

- If larger air bubbles can not be ejected, a valve or the suction tube is not fitting tight (see point 13. trouble-shooting).

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### 7.3 Volume configuration

- By turning the adjusting ring (in either direction) set the desired volume mark directly above the arrow.

#### **Achtung!**

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#### **Attention!**

Do this only when the piston is in down position.

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
---

### Attention!

Never use force when dispensing. If the piston becomes stiff the unit must be cleaned immediately (see point 8.).

---

- Hold the collection vessel under the ejection unit.
- Ensure that closure cap is removed.
- Gently and steadily pull up the piston to the top of its path.
- Push the piston down gently and steadily.

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## 7.5 After use, locking of dispensing mechanism

- After use the piston must always be left in the down position.
- Turning the volume scale to zero position locks the mechanism.

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### **Attention!**

The volume scale can only be turned to the zero position when the piston is in the down position.

- Place the closure cap on the ejection tube's opening.

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## 8. Cleaning

In order to ensure a long life for your dispenser, it should be cleaned regularly.

It is essential to clean it in the following cases:

- piston becoming stiff (see point 8.1)
- after use with liquids of which the use is restricted (see point 5.)
- before changing reagent
- before sterilisation
- before replacement of suction or ejection valves
- before storage

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
### Attention!

- When cleaning always wear protective clothing, goggles and gloves.
- Follow safety instructions (see point 1.).
- Suction tube, ejection unit, cylinder and valves contain liquid. Ensure that they are kept away from your body or other persons.

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## 8.1 Service/Cleaning

- After use set the scale ring to position "0" and close the cap (see point 7.5).
- Place the unit and bottle into a suitable cleaning container.

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

---

- Set the dispenser to maximum volume at the arrow (adjusting mark) through turning the volume scale (left/right).

- Unscrew the dispenser from the bottle. Ensure that protective gloves are worn.

---

**Attention!**

Liquid may drip from the suction tube.

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### 8.1 Service/Cleaning


- Hold the ejection unit above the reservoir bottle.
- Pull off the closure cap and slide it backwards.
- Move the piston up and down to eject the remaining liquid.

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## 8.1 Service/Cleaning

- Hold the suction tube into a suitable  
cleaning solution and rinse the dispenser  
completely by repeated dispensing.
- 

### **Attention:**

- Liquid may splash!
  - Always dispense away from your body!
- 
- Rinse the unit with suitable solvent by  
repetitive dispensing.
  - Remove the suction tube from the solvent  
and clear residual solvent from the dispenser  
by repeated dispensing.

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### 8.1 Service/Cleaning

- Dismantle the dispenser and clean all following parts separately after dismantling the dispenser as described following:
- Unscrew the ejection unit from the dispenser.
- Remove the union nut on the suction tube.
- Remove the suction tube from the dispenser.

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## 8.2 Intensive cleaning

- Clean the unit as described in point 8.1.
- Remove the piston as follows:
  - Set the maximum volume at the arrow by turning the adjusting ring.
  - Pull the piston right to the stop.
  - Hold the step gauge and unscrew the knob/piston.
  - Pull the piston out of the cylinder.

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
## 8.2 Intensive cleaning

- Clean the cylinder and the ceramic piston with a brush and detergent. Afterwards rinse these parts with an appropriate solvent.
- Reassemble the piston:  
Keep the unit vertical and insert the piston carefully into the cylinder, with a **smooth rotating** movement. Screw the knob (with the piston) onto the step gauge.

### **Attention!**

Take care that the piston is inserted straight and no force is used.

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---

## 9. Sterilisation

The dispenser is able to withstand autoclaving up to 121°C, 2 bar for a period of time acc. to DIN te = 20 min.

### 9.1 Preparation

Before sterilisation clean the unit as described in point 8.1.

### 9.2 Removal of the piston

---


#### **Hint!**

In order to achieve a better cleaning effect, remove the piston from the cylinder.

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## 9.2 Removal of the piston

- Set the maximum volume at the arrow, by turning the adjusting ring.
- Pull the piston right to the stop.
- Hold the step gauge and unscrew the knob/piston.
- Pull the piston out of the cylinder.

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### 9.3 Loosening the suction/ejection valve


- Loosen the suction valve with assembly wrench by turning it counterclockwise a half a revolution.

- Loosen the ejection valve with the assembly wrench by turning it counterclockwise a half a revolution.

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### **Attention!**

- Let the unit cool down to room temperature before reassembling.
- Check parts for damage and if necessary, exchange them.

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
---

## 9.5 Reassembly of piston

Keep the unit vertically and insert the piston **carefully** into the glass cylinder with a **smooth rotating** movement.

---

### **Attention!**

Take care that the piston is inserted  
straight and no force is used.

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## 9.6 Tightening the suction/ejection valve

- Tighten the suction valve firmly with the assembly wrench by turning it in a clockwise direction.
- Tighten the ejection valve firmly with the assembly wrench by turning it in a clockwise direction.

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## 9.7 Reassembly of the dispenser

Reassemble the unit (see point 6.) and check correct operation (see point 7.). In case of problems see point 13. trouble-shooting.

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## 10. Valve replacement

Leaking or damaged suction/ejection valves require immediate exchange (see point 15. Accessories and spare parts).

### 10.1 Suction valve replacement

- Clean the unit (see point 8.).
- Unscrew the suction tube.
- Hold the unit upside down, otherwise the valve ball might fall out, when changing the valve.
- Unscrew the suction valve with the assembly wrench and remove the valve ball.

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#### 10.1 Suction valve replacement

- Insert the ball again and screw on the new suction valve, tightening it with the assembly wrench.
- Put in the suction tube right to the stop and tighten the union nut on the valve (see point 6.1).

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## 10.2 Ejection valve replacement

- Realize cleaning (see point 8.).
- Unscrew the ejection unit (see point 8.).
- Screw out the ejection valve with the assembly wrench and remove.
- Screw in the new ejection valve and tighten firmly with the assembly wrench.
- Push on the ejection unit up to the stop and tighten the unit nut firmly (see point 6.2). Check to ensure a secure fitting!

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## 11. Repair by manufacturer

- Before returning a defective unit, the dispenser must be cleaned carefully by the user (see point 8.).
- Units will only be repaired if they are cleaned by the user and are no longer contaminated with any substances which could pose a risk to health.
- The clearance certificate should be filled out prior to sending in the unit (see point 16.).
- Not cleaned units will be returned to the sender without having been repaired.
- Describe the malfunction exactly.
- Indicate the dispensed liquid. (Important!)
- The warranty is invalid in case unit has been maintained other than described in this manual (e. g. use of force).

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## 12. Volumetric inspection

The terms for accuracy (R%) and the coefficient of variation (CV%), indicated in point 14., can be checked as follows:

### 1. Measuring procedure (see e.g. DIN EN ISO 8655-6):

- Set the nominal volume.
  - Dispense with de-ionized water.
  - Determine the dispensed quantity on an analytical balance.
  - Convert the measured weight to volume
- Note:** Temperature and air pressure need to be considered.

### 2. Repetitive measuring

- In order to obtain reliable measures for accuracy and coefficient of variation, the procedure described under 1. should be repeated at least 5 times.

### 3. Calculation of accuracy and coefficient of variation

- Accuracy (R%) and coefficient of variation (CV%) have to be calculated in accordance to the formulas of statistical quality control.

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### 13. Trouble-shooting

Symptom	Possible cause	Solution
Piston is stiff or does not move at all	<ul style="list-style-type: none"><li>- Closure cap not open</li><li>- Crystallisation at the piston or the valves</li></ul>	<ul style="list-style-type: none"><li>- Open the closure cap (see point 7.1) unit.</li><li>- Stop dispensing and clean the dispenser thoroughly (see point 8.2). If necessary, soak it for 24 hours in a bath (e. g. with 20% reasol® cleaning solvent).</li></ul>
Suction or dispensing is impossible	<ul style="list-style-type: none"><li>- Volume scale set on "0" position</li><li>- Suction/ejection valve sticking</li></ul>	<ul style="list-style-type: none"><li>- Turn the adjusting ring to set the required volume (see point 7.2 and 7.3)</li><li>- Clean the unit (see point 8.1)</li><li>- Shake the unit or tap the valve block against the palm of the hand.</li><li>- If necessary, change the valves (see point 10.).</li></ul>
Liquid is leaking from the gap between the valve block and the valves	<ul style="list-style-type: none"><li>- Valves loose or damaged</li></ul>	<ul style="list-style-type: none"><li>- Tighten valves with the assembly wrench <b>Note:</b> After sterilisation (see point 9.), suction/ejection valves need to be tightened.</li></ul>

### 13. Trouble-shooting

Liquid is leaking between the ejection unit and the valve

- Ejection unit is loose or damaged
  - Wrong ejection unit is used
  - Mount the original ejection unit correctly (see point 6.2).
  - Use only the original ejection unit.
- 

Unit aspirates in air

- Unit is not primed
  - Suction tube is loose
  - Suction tube or union nut damaged
  - Prime the unit (see point 7.2).
  - Tighten the suction tube (see point 6.1).
  - Replace original suction tube (see point 15. Accessories spare parts).
- 

Delivered volume is incorrect

- Ejection unit loose or damaged
- Suction tube loose or damaged
- Suction/ejection valves loose or damaged
- Install a genuine ejection unit correctly (see point 6.2).
- Tighten the suction tube (see point 6.1).
- Tighten the valves with the assembly wrench.
- Clean the unit (see point 8.1).
- If necessary, change the valves (see point 10.).

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Dispenser Volumen ml	Einstell- volumen ml	R* %	CV* %
0,2 - 1	0,05	* 0,6	* 0,2
0,4 - 2	0,05	* 0,6	* 0,2
1 - 5	0,10	* 0,6	* 0,2
2 - 10	0,25	* 0,6	* 0,2
5 - 30	0,50	* 0,6	* 0,2
10 - 60	1,00	* 0,6	* 0,2

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- \* The given values for accuracy (R) and coefficient of variation (CV) refer to the maximum volume at a:
- consistent temperature (20 °C) of the unit, surrounding and liquid (de-ionized water).
  - even and jerk-free operation.

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#### 14. Technical data/Materials

The dispenser is conformity certified acc. to the prescription of the German „Eichordnung“ (Enclosure 12).

Area of temperatures for use of the unit and the liquid between 4 °C and 50 °C.

	Kolben	Zylinder	Ventil- block	Ventil- kugel	Ventil- sitz	Ventil- feder
	piston	cylinder	valve block	valve ball	valve seat	valve spring
<b>Dispenser</b>	hochreine Al <sub>2</sub> O <sub>3</sub> -Keramik (99,7%)	Borosilikat- Glas 3.3	Halar (ECTFE)	hochreine Al <sub>2</sub> O <sub>3</sub> -Keramik (99,7%)	hochreine Al <sub>2</sub> O <sub>3</sub> -Keramik (99,7%)	Hastelloy (99,7%)
	highly pure Al <sub>2</sub> O <sub>3</sub> -ceramic (99,7%)	Borosilicate glass 3.3	Halar (ECTFE)	highly pure Al <sub>2</sub> O <sub>3</sub> -ceramic (99,7%)	highly pure Al <sub>2</sub> O <sub>3</sub> -ceramic (99,7%)	Hastelloy (99,7%)
<b>Dispenser PP</b>	Borosilikat- Glas 3.3	Borosilikat- Glas 3.3	PP	hochreine Al <sub>2</sub> O <sub>3</sub> -Keramik (99,7%)	hochreine Al <sub>2</sub> O <sub>3</sub> -Keramik (99,7%)	Hastelloy (99,7%)
	Borosilicate glass 3.3	Borosilicate glass 3.3	PP	highly pure Al <sub>2</sub> O <sub>3</sub> -ceramic (99,7%)	highly pure Al <sub>2</sub> O <sub>3</sub> -ceramic (99,7%)	Hastelloy (99,7%)

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Unit serial no. ....

Module serial no. ....

Blood, bodily fluids, patholog. samples

**YES**

**NO**

**Please define in case of YES**

\_\_\_\_\_

Other hazardous biologic substances

\_\_\_\_\_

Biologically-degradable materials which can pose a risk

\_\_\_\_\_

Chemicals containing substances harmful to health

\_\_\_\_\_

Radioactive materials

\_\_\_\_\_

Other hazardous and toxic substances

\_\_\_\_\_

### We hereby declare that

- the device was carefully emptied, cleaned and decontaminated prior to shipment,
- no risks resulting from contamination exist,
- we are aware that liability claims can be made against us in the event of damage being caused by contaminated units.

### Address of sender:

Name: \_\_\_\_\_

Date: \_\_\_\_\_

Signature: \_\_\_\_\_

Data from this form is not processed electronically. The purpose of the form is to be able to prove the correct responsibility for any liability issues caused by contaminated material. This declaration will be destroyed no later than 4 weeks after completion of the process for the goods sent in. In the event of a liability case arising, the data is processed in accordance with current legislation.

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## 17. Declaration of conformity

Declaration of conformity for measuring instruments according to the new German Federal Measures Regulation and the adaptation to European Jurisdiction of December 11, 2014, Annex 5 § 11 Paragraph 2

1. Item: Bottle top dispenser ceramus® classic and EM-dispenser 1 ml - 60 ml

2. Item number: 93X 00 .., 93X 11 .., 93X 02 ..

3. Manufacturer:

Hirschmann Laborgeräte GmbH & Co. KG  
Hauptstraße 7-15 · 74246 Eberstadt Germany  
Fon +49 7134 511-0 · Fax +49 7134 511-990  
info@hirschmannlab.com · www.hirschmannlab.com

4. Hirschmann Laborgeräte GmbH & Co. KG herewith confirms that the above mentioned item complies with the Weights and Measures Regulations as well as the statutory ordinance it is based on.

5. Indication of other technical regulations or specifications:  
ISO 8655-5, ISO 8655-6

6. The declaration of conformity is based on the following regulations:  
- Annex 12 to the Weights Regulation of its version prevailing on December 31, 2014

- PTB-requirements 12 (12/2010)


7. Indication of relevant harmonized standards or normative reference:  
Not applicable

8. Conformity assessment body: Not applicable

9. Additional indication: Not applicable

10. It is the manufacturer's sole responsibility to issue this declaration of conformity.

Issue place and date  
Eberstadt, October 9, 2017

Name, capacity, signature  
  
Roland Benz, QA-Manager

**Notizen/Notes :**

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---