# ROTOR

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**Year of manufacture:** A six digit appliance number is printed on the rating plate located on the bottom plate of the machine. The first two digits indicate the year of manufacture for that machine.

## 1. Safety tips

The Rotor Bar Blender (type RBB) is built in accordance with the relevant international standards. When using electrical appliances, certain basic safety precautions must be complied in order to prevent any possible injuries or damages. Read these instructions carefully and retain them in case of need.

**If you use the Rotor Bar Blender in an improper way or do not respect these instructions, Rotor Lips Ltd. cannot be held liable.**

- Before putting the appliance into operation, check whether the voltage value displayed on the rating plate, corresponds with that of the mains voltage, and that the appliance is equipped with the proper plug.
- Check the mixer before each use (appliance, plug and cable).
- **The appliance may not be used if it is damaged.**
- If the appliance has been dropped or has fallen, take it to an authorized service center for inspection.
- Repairs should only be made at an authorized service center (see chapter “Repairs”). Serious dangers to the user could result from repairs which are carried out improperly.
- The plug should only be plugged into a socket when the appliance is switched off. Do not touch plug contacts. Do not allow the mains cable to hang down (children!)
- The mixer should only be used by **instructed personnel**.
- The mixer should only be used when all the parts have been properly installed.
- Use only appropriate, original spare parts from Rotor Lips Ltd.
- All persons (including children), who are not in the situation to use the appliance in a proper way, due to their physical, sensoric or mental capability, their inexperience or ignorance, may only use this machine under supervision or instruction by a responsible person in charge.

**To avoid any accidents and damages on the appliance, run the machine be only with an assembled container and set lid. It is forbidden, while mixing, to remove the lid lock and to reach into the container with a hard object.**

- Never leave the appliance unattended while it is plugged in.
- Keep your hands, hair and other objects away from the moving parts.
- Never attempt to pull out the plug by the cable or with wet hands.
- Do **not** allow the mains cable to hang over a sharp edge of a table, do **not** allow it to lie on a work surface, and avoid having it come into contact with hot surfaces.
- Use only proper, approved extension cords.
- Always unplug the appliance before cleaning.
- The motor unit may **never** be submerged in water or any other liquid!
- **Never clean it with a water jet!** Wipe clean with a damp cloth.

## 2. USES

The mixer type RBB should **only** be used to stir, mix, puree and emulsify drinks.

For example : Cocktails and Milkshakes.

Recommended processing volume: max. 0,4 Liter (400 W) resp. 1,7 liter (550 W)

Larger fruits and vegetables should be cut into pieces before they are put into the bowl.


For the preparation of sauces, mayonnaise, soups, custards, cake batters, etc. we recommended to use the Rotor Gastronom GT kitchen blender.

For crushing ice cubes use the option cutter "Bar" (Nr. 1121.306).

**If the Rotor Bar Blender is used for any other purposes or altered, Rotor Lips Ltd. cannot be held responsible for any damages which result.**


### Container uses

- 2 l Polycarbonate transparent, light and manageable.
- 2 l Stainless steel solid, for longest life perfect look.
- 4 l **Must not been used!** it is not according the CE rules. →




## 3. MACHINE DESCRIPTION

- 1. Motor base
- 2. Rubber foot
- 3. Bottom plate
- 4. Rating plate
- 5. Container
- 6. Lid
- 7. Lid lock
- 8. Main switch (on/off) with pulse function
- 9. Speed selection switch
- 10. Thermal protection switch
- 11. Cutter head
- 12. Threaded ring

**4. TECHNICAL DATA**

Power supply: (version RBB2)	230 V/50 Hz/1,8 A, or 120 V/60 Hz/4,2 A (tension +/- 6%, frequency +/- 2 Hz, sinusoidal, interference free)
Power supply : (version RBB3)	230 V/50 Hz/3,2 A, or 120 V/60 Hz/5,8 A (tension +/- 6%, frequency +/- 2 Hz, sinusoidal, interference free)
Power requirement: (version RBB2)	400 W continuous operation (max. 700 W)
Power requirement: (version RBB3)	550 W continuous operation (max. 1 000 W)
Motor:	- High temperature resistant collector motor with fan - max. speed ca. 23 000 rpm
Motorcontrol: (version RBB2)	- 2 speeds
Motorcontrol: (version RBB3)	- 3 speeds (automatically regulated)
Electrical equipment:	- according to standard EN-60335-1 (grounded/protection class I) - electromagnetically shielded accordance with EMC-directives 89/336 EEC (230 V) - the motor has a thermal protection switch
Jar content :	2 Liter (processing amount up to 1,7 Liter, depending on the consistency and motor version).
Construction:	- Motor base is made out of stainless steel - Cutter blades made out of stainless steel
Weight:	4 kg net, with stainless steel container 4,6 kg
Temperature range:	operating: 10...40°C storage: -20...70°C
Humidity:	operating: 0...80% non condensing storage: 0...95% non condensing
Altitude:	0..2.000 m. (0...6600 ft)

We reserve the right to make technical changes at any time.

## 5. ELECTRICAL CONNECTION

The **power supply voltage must correspond** to the values on the rating plate on the appliance (the voltage rating is given in volts [V])

The appliance should only be connected to a **single phase grounded alternating current**. Make sure that the original mains cable and the wall outlet correspond.

The use of a **faulty-energy-protection switch**, (FI) is suggested by all the works with electro-appliances.

Appliances with a 230 V designation can also be used with a 220 V mains supply.

If at all possible, avoid using an extension cord. If it is unavoidable, only proper, approved extension cords should be used.

Do not subject the mains cable to strong tension or pressure.

**To replace the cable**, see chapter "Repairs".

## 6. PUTTING INTO OPERATION

Before using the appliance for the first time it must be **thoroughly cleaned**.

Follow the instructions in the chapter "Cleaning".

Place the Rotor Bar Blender on a **firm, flat base** and plug it into a nearby, appropriate electrical outlet.

### **Note**

The windings of the Rotor commutator motors are covered with high-quality insulating lacquer. For this reason, the motor may, at first, exude an odour when it heats up. This is normal and does not damage the motor; it disappears after the motor has been used several times.

## 7. OPERATING INSTRUCTIONS

Container assembly (see the drawings at the end)

- Insert the cutter head together with the gasket (art.no. 1121.324) into the threaded ring (art.no. 1121.112).
- Screw the bowl turning counter clockwise into the threaded ring.
- Always ensure that it is inserted straight.

**The cutter blades are sharp and may cause bad injuries.**

#### Mix

- Insert the container into the retaining ring on the motor base, taking care to keep it vertical.
- Add ingredients to be mixed.
- Place the mixing bowl on the rubber ring so that the switching rod to the bowl monitor is pressed and the blender attachment is properly and securely fixed.
- Place the lid on firmly.
- Switch on motor (main switch "on") and vary the speed according to the ingredients being mixed (selection switch).
- Turn off: at the main switch ("off").
- Remove the container from the motor base.

Liquids and fine grained materials can be added through the filler opening in the lid of the container during operation (Remove the lid lock).

By switching on and while operating at high speed it is recommended that the bowl should be held by hand (lay hand on the lid).

English

**Never use the mixer without lid!**

**Hold the lid by hand!**

To avoid heavy injuries or damages of the appliance, it is forbidden, while mixing to reach into the container with hands or with a hard solid object. This is although valid by a standing still motor and assembled container, in case of unintentional switching. Do not reach into the container with a hard solid object by opening the lid lock.

The mixer should **never run without container.** Do not touch, nor with any other object, the rotating coupling piece on top of the motorbase.

**Never put a cutter assembly on the motor base without the container (extreme dangerous!)**

**Danger of scalding with hot drinks!**

#### Container removal

First remove the complete container from the motor base and empty it totally.

**Attention:** To disengage the plastic threaded ring, turn to the right. Clockwise, like closing a water tap, as the threaded ring is **left-threaded!**  
**Do not tear the handle (risk of breakage).**

## 8. MOTOR PROTECTION (OVERLOAD PROTECTION SWITCH)

When the machine is overloaded or overheated, it switches off automatically.

**When the machine stops automatically, immediately remove the complete container from the motor unit, pull out the plug and switch off the machine.**

After about 60 seconds press the black button on the bottom plate of the motor base and switch on the machine as usual.

Charge the machine less than before.

If you notice any irregularities or defects, please contact your supplier. For more information see chapter "Repairs".

## 9. CLEANING

**Pull out the plug before cleaning!  
The motor unit must never be submerged in, or sprayed with water or any other liquid! Never attempt to clean it with a water jet!**

If necessary, clean the motor unit with a damp cloth.

**Container:** The mixing container should not be cleaned in the dishwasher while still installed. They must first be dismantled.

The container should be cleaned thoroughly after each use or once a day.

**Bowl:** Can be cleaned in the dishwasher or clean by hand using hot water and a commercial cleaning product.

**Cutter assembly:** simply rinse and, if necessary, scrub off. **Never place into water or the dishwasher.**

**Threaded ring:** should be cleaned especially inside, by using a hard brush.

## 10. MAINTENANCE

The motor bearings are permanently lubricated and therefore require no maintenance.

**11. TROUBLESHOOTING**

<b>PROBLEM</b>	<b>CAUSE/SOLUTION</b>
The appliance doesn't run	<ul style="list-style-type: none"> <li>- See chapter "Safety features and Protection"; overload switches</li> <li>- Make sure the mains cable is plugged into a wall socket</li> <li>- Check the fuse in the applicable building electrical circuit</li> <li>- Brushes worn (more than 1 000 operating hours) contact your service center</li> </ul>
The motor only runs at highest speed when the appliance is switched on	- Defective electronics (service center)
The motor stops under load	<ul style="list-style-type: none"> <li>- The overload switch has been activated</li> <li>- Brushes worn (service center)</li> <li>- Defective electronics (service center)</li> </ul>

**12. REPAIRS**

Obtain the address of your nearest authorized Rotor Lips service center, from the store where you purchased your appliance.

**Repairs on electrical appliances should only be carried out by professional technicians. Serious dangers to the user could result from improper carried out repairs.**

**Attention!**

**If the appliance's mains cable should become damaged, it should only be replaced by a professional technician with an original cable. (Do not forget to branch grounding!)**

**13. DISPOSAL/RECYCLING**

Disposal of the resulting work waste materials (e.g. food packaging, peelings etc.) must be carried out according to the local regulations.

All individual components of the unit have to be treated according to the local regulations for recycling (e.g. metal etc.) or safe disposal (e.g. P.C. board etc.) This product should not be treated as household waste. It is RoHS compliant.

**14. GUARANTEE**

The terms of the guarantee are in accordance with the sales and delivery regulations of Rotor Lips Ltd.

**15.  Declaration of conformity**

We Rotor Lips Ltd.  
Glütschbachstrasse 91  
CH-3661 Uetendorf  
Switzerland

declare under our sole responsibility that the machine

**Drink Blender Bar Blender**

**Type RBB 2** (2 speeds/400 W)

**RBB 3** (3 speeds/550 W)

to which this declaration relates, is in conformity with the following standards or other normative documents:

EN 60335-1:2002/A12:2006, 55014-1:2000, 55014-2:1997/A1:2001, EN 55104:1995,  
EN 550011:1998/A2:2002, EN 60335-2-14:2003, EN 61000-3-2:1995,  
EN 61000-3-3:1995

Following the provisions Directives:

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