Pipette / Pipette / Pipette / Pipeta

# Rotilabo<sup>®</sup>- Mikropipette

Gebrauchsanweisung  
Operating Manual  
Mode d'emploi  
Instrucciones de manejo

Carl Roth GmbH + Co. KG  
Schoemperlenstraße 3-5  
76185 Karlsruhe  
Germany

Telefon: +49 721 5606 - 0  
Fax: +49 721 5606 - 149  
E-mail: [info@carloth.de](mailto:info@carloth.de)  
Internet: [www.carloth.de](http://www.carloth.de)

---

---

---

---

## Inhalt

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# 1. Safety Instructions

## Please read the following carefully!

This instrument may sometimes be used with hazardous materials, operations, and equipment. It is beyond the scope of this manual to address all of the potential safety risks associated with its use in such applications. It is the responsibility of the user of this pipette to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

1. Every user must read and understand this operating manual prior to using the instrument and observe these instructions during use.
2. Follow general instructions for hazard prevention and safety instructions; e.g., wear protective clothing, eye protection and gloves. When working with infectious or other hazardous samples, all appropriate regulations and precautions must be followed.
3. Observe the reagent manufacturers' information.
4. Only use the instrument for pipetting liquids that conform to the specifications defined in the limitations of use and operating limitations. Observe operating exclusions (see page 5). If in doubt, contact the manufacturer or supplier.
5. Always use the instrument in such a way that neither the user nor any other person is endangered. Avoid splashes. Use only suitable vessels.
6. Avoid touching the tip orifices when working with hazardous samples.
7. Never use force on the instrument!
8. Use only original spare parts. Do not attempt to make any technical alterations. Do not dismantle the instrument any further than is described in the operating manual!
9. Before use check the instrument for visual damages. If there is a sign of a potential malfunction (e.g., piston difficult to move, leakage), immediately stop pipetting. Consult the 'Troubleshooting' section of this manual (see page 17), and contact the manufacturer if needed.

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## 2. Limitations of Use

### 2.1. Purpose

Air-displacement pipette for pipetting aqueous solutions of medium density and low to medium viscosity.

### 2.2. Limitations of Use

The instrument is intended for the pipetting of liquids within the following limitations:

- Temperature of both the instrument and solution should be between 15 °C to 40 °C (59 °F to 104 °F). Consult the manufacturer for use in temperatures outside of this range.
- Vapor pressure up to 500 mbar
- Viscosity: 260 mPa s (260 cps)

### 2.3. Operating Limitations

Viscous and highly adhesive liquids may impair volumetric accuracy. Volumetric accuracy may also be impaired when pipetting liquids that differ from ambient temperature by more than  $\pm 1$  °C /  $\pm 1.8$  °F.

### 2.4. Operating Exclusions


The user has to ensure the compatibility of the instrument with the intended application.

This instrument cannot be used:

- for liquids incompatible with polypropylene
- for liquids incompatible with polycarbonate (inspection window)
- for liquids of a very high vapor pressure
- for liquids incompatible with FKM and polyetheretherketone

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\*Tip ejection key  
The serial number is behind the tip ejection key.

\*\*Volume display  
The figures in the display are read from top to bottom, the dash represents the decimal point.

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**Note:**

Optimum analysis results can only be obtained with quality tips.

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
### 1. Fitting the tip

Use the correct tip according to the volume range or the color code.  
Ensure that the tip is securely seated.

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
### 2. Volume setting

Select the desired volume by rotating the volumes setting wheel. Avoid twisting and abrupt rotating motions during this adjustment.

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
### 3. Aspirate sample

- Press pipetting button to the first stop.
- Hold the pipette vertically and immerse the tip into the liquid.

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
- Let the pipetting button slide back slowly. In order for the liquid to reach its end position, leave the tip immersed for a few seconds.
- Touch the tip against the container wall.

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
### 4. Discharge sample

- Place the pipette tip against the wall of the vessel. Hold the pipette at an angle of 30-45° relative to the container wall.
- Press the pipetting button slowly to the first stop and hold it down. For serum and liquids of high viscosity or low surface tension, observe adequate waiting time to improve accuracy.

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
c) The blow-out stroke empties the tip completely: Press the pipetting button down to the second stop.

d) While doing this, wipe the pipette tip against the wall of the container.

e) Remove the pipette tip from the container wall and let the pipetting button slide back.

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## 5. Ejecting the tip

Hold the pipette shaft over a suitable disposal container and press the tip ejection key to the stop.

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### Important!

Don't lay the instrument horizontal when the tip is filled. Liquid may enter and contaminate the instrument.

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## 5. Checking the Volume

Depending on use, we recommend inspection of the instrument every 3 to 12 months. The cycle can, however, be adjusted to individual requirements.

The gravimetric testing of the pipette volume is performed according to the following steps and is in accordance with DIN EN ISO 8655, Part 6.

### 1. Set nominal volume

Set volume to the maximum volume indicated on the instrument (see page 7 for procedure).

### 2. Condition the pipette

Condition the pipette before testing by using a pipette tip to aspirate and discharge the test liquid (distilled  $H_2O$ ) five times. After this, discard the pipette tip.

### 3. Carry out the test

#### Note:

According to DIN EN ISO 8655-2, it is recommended to replace the tip after each measurement. This rule may be waived in accordance with the DAkkS guideline DKD-R8-1.

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- a) Aspirate liquid and pipette it into the weighing vessel.
- b) Weigh the pipetted quantity with an analytical balance. (Please follow the operating manual instructions from the balance manufacturer.)
- c) Calculate the volume, taking into account the temperature of the test liquid.
- d) At least 10 pipettings and weighings in three volume ranges (100%, 50%, 10% of nominal volume) are recommended for statistical analysis.

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### 5.1 Calculation for nominal volume $V_0$

$x_i$  = Weighing results  
 n = Number of weighings  
 Z = Correction factor (e.g. 1,0029 µl/mg at 20 °C, 1013 hPa)

#### Richtigkeit\* / Accuracy\*

$$R/A\% = \frac{\bar{V} - V_0}{V_0} \cdot 100$$

#### Variationskoeffizient\* / Coefficient of Variation\*

$$VK/CV\% = \frac{100 s}{\bar{V}}$$

\*) = Calculation of accuracy (A%) and variation coefficient (CV%): A% and CV% are calculated according to the formulas for statistical control.

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## 6. Accuracy Table

Rotilabo®-Micropipette, Digital adjustable

Final test values related to the nominal capacity (maximum volume) indicated on the instrument, and the indicated volume steps obtained when instrument and distilled water are equilibrated at ambient temperature (20 °C/68 °F) and with smooth operation. According to DIN EN ISO 8655.

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## 7. Adjustment

The instrument is permanently adjusted for aqueous solutions. If the pipette operation is clearly inaccurate, or if the instrument must be adjusted for solutions of different densities and viscosities or specially-shaped pipette tips, adjustments can be made.

1. Check the volume, determine actual value ( $\bar{V}$ ) (see page 8).

2. Remove the cover: Push the hook forward, raise it slightly and then pull it back.

3. Using a paperclip or a pipette tip, remove the protective film (this protective film can be discarded).

4. Push the red adjustment slider completely back, raise the volume-setting wheel (decoupling) and release the adjustment slider.

5. Set the adjustment value: set the volume-setting wheel to the previously determined actual value ( $\bar{V}$ ). A volume check is recommended after every adjustment.

6. Push the adjustment slider completely back again, push the volume-setting wheel downwards and release the adjustment slider. Re-insert the cover.

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### Note:

The change to the factory settings is indicated by the now visible red adjustment slider.

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## 8. Autoclaving

The pipette is completely autoclavable at 121 °C (250 °F), 2 bar absolute (30 psi) with a holding time of at least 15 minutes, according to DIN EN 285.

### **Attention!**

Prior to autoclaving, the volume adjustment must be set on an available numbered volume (e.g., 11.25 or 11.26 but not between).

1. Eject the pipette tip.
2. Autoclave the complete pipette without any further disassembling.
3. Allow the pipette to completely cool and dry.

### **Note:**

The effectiveness of the autoclaving must be verified by the user. Maximum reliability is obtained with vacuum sterilization. We recommend the use of sterilization bags.

If the pipette is autoclaved frequently, the piston and the seal should be greased with silicone grease in order to preserve smooth movement. Please use only the recommended silicone grease, see accessories page 16. If necessary after sterilization, tighten the connection between the hand grip and the pipette shaft.

## 9. UV sterilization

The unit can withstand the usual output of a UV sterilization lamp. The effects of the UV may cause some color change.

## 10. Filter pipette 5 ml / 10 ml

A hydrophobic PE filter is used as a safeguard against liquid entering the pipette.

Change the filter if it becomes wet or contaminated.

- Use a flat object such as a screwdriver
- Remove the filter without damaging the tip cone.

Remove the filter before autoclaving!

The instrument can be operated without a filter.

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## 11.1. Rotilabo®-Micropipette up to 1000 µl

### 1. Servicing

Inspect the pipette tip cone for damage. Inspect the piston and seal for contamination.

Test the instrument's piston seal. To do this affix a tip, and aspirate a sample. Hold the instrument vertically, with the sample in the tip for approximately 10 s. If a drop forms at the tip orifice, see the troubleshooting guide on page 17.

### 2. Disassembly and cleaning

1. Unscrew the pipette shaft (S) from the hand grip.
2. Unscrew the upper part of the ejector (A) from the pipette shaft.
3. Pull the shaft (B, C and D) out of the lower part (E) of the ejector.
4. Unscrew the piston unit (B).

**Note:**

Piston remains connected with piston unit (B)!

5. Remove the seal with spring (C) (this is non-removable on 10 µl pipette models).
6. Clean the parts shown with a mild soap solution or isopropanol and then rinse with distilled water.
7. Allow the parts to dry (max. 120 °C/ 248 °F).
8. Grease piston and seal with a very thin layer of silicone grease. Please use only the recommended silicone grease, see accessories page 16.
9. Assemble the ambient temperature parts in reverse order from above. Piston unit and upper part of the ejector (A, B) should only be hand-tight.

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## 11.2. Rotilabo®-Mikropipette 5 ml and 10 ml

### 1. Servicing

Inspect the pipette tip cone for damage. Inspect the piston and O-Ring-seal for contamination.

Test the instrument's piston seal. To do this, affix a tip, and aspirate a sample. Hold the instrument vertically, with the sample in the tip for approximately 10 s. If a drop forms at the tip orifice, see the troubleshooting guide on page 17.

### 2. Disassembly and cleaning


1. Remove the entire shaft (S) from the handgrip by rotating at the upper end of the ejector (F) and remove the filter (K) from the bottom part of the shaft (H).
2. Separate the bottom part of the ejector (F') by unscrewing it from the upper part of the ejector (F).
3. Unscrew and dismantle the piston unit (G) with the ejector spring (I) and the bottom part of the shaft (H).
4. Remove the O-Ring-seal from the piston unit and clean it.

#### Note:

Do not disassemble piston unit (G) any further!

5. Clean piston unit (G) and lower part of pipette shaft (H) with a soap solution or isopropanol and then rinse with distilled water.
6. Allow the parts to dry (max. 120 °C / 248 °F) and to cool down.
7. Carefully lubricate the inside and outside of the O-ring and mount it on the piston.
8. Assemble the individual components in the reverse order from that described above.

(For illustration purposes only)

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### Rotilabo®-Mikropipette

Volumen	Best.-Nr.
0,5 - 10 µl	TA21.1
2 - 20 µl	TA22.1
10 - 100 µl	TA23.1
20 - 200 µl	TA24.1
100 - 1000 µl	TA25.1
0,5 - 5 ml	TA26.1
1,0 - 10 ml	TA27.1

Regalhalter für 1 Pipette  
Best.-Nr. TA28.1

## 12. Ordering Information and Accessories

### Rotilabo®-Micropipette

Volume	Cat. No.
0,5 - 10 µl	TA21.1
2 - 20 µl	TA22.1
10 - 100 µl	TA23.1
20 - 200 µl	TA24.1
100 - 1000 µl	TA25.1
0,5 - 5 ml	TA26.1
1,0 - 10 ml	TA27.1

Shelf/rack mount for 1 pipette  
Cat. No. TA28.1

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
### 13.3. Additional Accessories

**Filter für Rotilabo®-Mikropipette 5 ml**  
VE 25 Stk. Best.-Nr. TA30.1

**Filter for Rotilabo®-Micropipette 5 ml,**  
pack of 25 Cat. No. TA30.1

**Filter für Rotilabo®-Mikropipette 10 ml**  
VE 25 Stk. Best.-Nr. TA31.1

**Filter for Rotilabo®-Micropipette 10 ml,**  
pack of 25 Cat. No. TA31.1

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**Silicone grease for Rotilabo®-Micropipette**  
up to 1000 µl Cat. No. TA32.1

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## 14. Troubleshooting

Problem	Possible cause	Corrective action
Tip dripping (instrument leaks)	<ul style="list-style-type: none"> <li>- Unsuitable tip</li> <li>- Tip not seated tightly</li> </ul>	<ul style="list-style-type: none"> <li>- Only use high-quality tips</li> <li>- Press tip on firmly</li> </ul>
The instrument does not aspirate or aspirates too little; the discharged volume is too low	<ul style="list-style-type: none"> <li>- Seal contaminated</li> <li>- The seal or cone is damaged</li> <li>- The piston is contaminated or damaged</li> </ul>	<ul style="list-style-type: none"> <li>- Clean seal</li> <li>- Replace seal or shaft</li> <li>- Clean or replace piston</li> </ul>
Aspiration is too slow	<ul style="list-style-type: none"> <li>- Shaft clogged</li> <li>- The filter in the 5 ml and 10 ml models is contaminated</li> </ul>	<ul style="list-style-type: none"> <li>- Clean shaft</li> <li>- Change the filter</li> </ul>
Discharged volume is too large	<ul style="list-style-type: none"> <li>- Pipetting button pressed too far into the blow-out position before sample uptake</li> </ul>	<ul style="list-style-type: none"> <li>- Operate properly. See ‚Pipetting‘, page 7.</li> </ul>
Piston is difficult to move	<ul style="list-style-type: none"> <li>- The piston is contaminated or needs oil</li> </ul>	<ul style="list-style-type: none"> <li>- Clean and oil piston</li> </ul>

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## 15. Repairs - Calibration Service

If a problem cannot be fixed by following the troubleshooting guide, or by replacing spare parts, then the instrument must be sent in for repair.

**For safety reasons, instruments returned for checks and repairs must be clean and decontaminated!**

### 15.1. Return for repair

- a) Clean and decontaminate the instrument carefully.
- b) Complete the „Declaration on Absence of Health Hazards“ (ask your supplier or manufacturer for the form).
- c) Send the completed form along with the instrument to the manufacturer or to the dealer with an exact description of the type of malfunction and the media used.

The return transport of the instrument is at risk and cost of the sender.

### 15.2. Calibration Service

ISO 9001 and GLP-guidelines require regular examinations of your volumetric instruments. We recommend checking the volume every 3-12 months. The interval depends on the specific requirements on the instrument. For instruments frequently used or in use with aggressive media, the interval should be shorter.

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## 16. Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the instrument or the consequences of normal wear and tear especially of wearing parts such as pistons, seals, valves and the breakage of glass as well as the failure to follow the instructions of the operating manual. We are not liable for damage resulting from any actions not described in the operating manual or if non-original parts have been used.

## 17. Disposal

For the disposal of instruments and tips, please observe the relevant national disposal regulations.

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