# KCL-Informationsbroschüre

## Kat. II

Art. 612+

<b>D</b>	Informationsbroschüre	<b>616+</b>
<b>GB</b>	Information brochure	<b>617+</b>
<b>F</b>	Brochure d'information	<b>618+</b>
<b>N</b>	Informasjonsbrosjyre	<b>620+</b>
<b>DK</b>	Informationsbrochure	<b>627+</b>
<b>FIN</b>	Tiedotuslehtisen	<b>629+</b>
<b>H</b>	Információs füzet	
<b>S</b>	Informationsbroschyr	
<b>E</b>	Folleto Informativo	
<b>I</b>	Opuscolo Informativo	
<b>NL</b>	Informatie brochure	
<b>SK</b>	Informačná brožúra	
<b>HR</b>	Informacije brošura	
<b>EST</b>	Infobrošüür	
<b>CZ</b>	Informační brožura	
<b>P</b>	Brochura de informação	
<b>BG</b>	Информационна брошура	
<b>GR</b>	Ευήμερωτικό φυλλάδιο	
<b>LT</b>	Informacinę brošiūra	
<b>PL</b>	Broszura informacyjna	
<b>LV</b>	Direktivos	
<b>RO</b>	Broşuri de Informare	
<b>SLO</b>	Informacijska brošura	
<b>TR</b>	Direktif	
<b>RUS</b>	Информационная брошюра	

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**Information brochure in section 1.4 of Annex II of Directive 89/686/EEC  
Personal Protective gloves Cat. II, Article 612+, 616+, 618+, 617+ 620+, 627+, 629+**

The CE marking on this glove means that it meets the essential requirements for in the European Directive EEC 89/686 concerning Individual Protection Equipment (IPE): Harmlessness-Comfort-Dexterity-Sturdiness.

Protective Gloves conform to the EN 420 . Approved to EN 388 according to Article 10 of Directive 89/686/EEC: CTC, 4, rue Hermann Frenkel, 69367 LYON- Cedex 07-France, Identification number: 0075 TÜV Rheinland LGA Products GmbH, Tillystr. 2, 90431 Nürnberg, Identification number: 0197

Article	Name	Size	EN 388	Testing institute
612+	Camapur® Cut 612+	6, 7, 8, 9, 10, 11	424X	0075
616+	Camapur® Comfort 616+	6, 7, 8, 9, 10, 11	3131	0075
617+	Camapur® Comfort 617+	6, 7, 8, 9, 10	1130	0197
618+	Camapur® Cut 618+	6, 7, 8, 9, 10, 11	4342	0075
620+	Camapur® Cut 620+	6, 7, 8, 9, 10, 11	4342	0075
627+	Camapur® Cut 627+	6, 7, 8, 9, 10, 11, 12	4341	0075
629+	FiroMech® 629+	6, 7, 8, 9, 10, 11	2010	0075

**EN 388, mechanic risks**

- 1st digit Abrasion resistance (min. 0; max. 4)
- 2nd digit Cut resistance (min. 0; max. 5)
- 3rd digit Tear strength resistance (min. 0; max. 4)
- 4th digit Puncture resistance (min. 0; max. 4)
- X = Cannot be tested

These gloves are suitable for activities where protection from mechanical injury is required (eg installation, dirty work).

**Warning!**

These protective gloves provide no protection against extreme cold (<-5°C), heat (>50°C), electricity, Chemical.

Do not use in the vicinity of moving machine parts. Risk of being drawn into the machine.

The specified protective index may be affected by factors such as temperature, abrasion, etc.

616+, 618+, 620+, 627+, 629+: The specified level of protection is only present on the coated surfaces.

The length of the gloves does not conform to DIN EN 420.

Unusable when the gloves are cracked, porous and stiff.

Safety glove not approved for food.

**Allergy Notice:** The gloves contain natural rubber latex that may trigger allergic reactions including anaphylactic reactions. The gloves may contain traces of dithiocarbamates.

Article 616+, 617+, 618+, 620+, 627+, 629+: Protective gloves contain traces of DMF in legally permitted ppm concentrations.

**Use:** The protective gloves you wear must be of the correct size. Note that using undergloves may result in some usage

restrictions. Check the protective gloves for damage before you use them. Damaged protective gloves must not be used under any circumstances.

**Expiry date:** If stored properly, the gloves can be used 36 months.

**Cleaning:** Article 612+, 618+, 620+, 627+: Delicate wash cycle up to 30°C with greatly reduced mechanical action. Can be dried on care programme. The specified properties are retained even when the gloves have been washed three times.

Article 616+, 617+, 629+: Delicate wash cycle up to 30°C with greatly reduced mechanical action. Can be dried on care programme. The certified index values are not guaranteed after washing/drying. Do not use any chemicals or sharp implements when cleaning (wire brushes, sandpaper, etc.).

**Storage/Transport:** Flat, dry, dark, with no additional weight load in its original packaging, at a temperature of 5 °C- 25 °C. Protect from sunlight and ozone.

**Disposal:** Gloves when used properly disposed of with household.

**For further information please contact.**

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KCL GmbH