2000014574

RV 3 eco\_062017

# IKA®

## IKA® RV 3 eco


Betriebsanleitung	DE	3	Руководство по эксплуатации	RU	73
Operating instructions	EN	17	Instruções de serviço	PT	87
Mode d'emploi	FR	31	Instrucja obsługi	PL	101
Instrucciones de manejo	ES	45	115	AR	تعليمات التشغيل
Istruzioni per l'uso	IT	59	Çalışma talimatları	TR	129

IKA-Werke, Germany  
Reg. No. 004343

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## Contents

	Page		Page
Device setup	2	Commissioning	23
Declaration of conformity	17	Maintenance and cleaning	25
Explication of warning symbols	17	Accessories	26
Safety instructions	17	Error messages	28
Correct use	19	Technical data	29
Unpacking	19	Solvent table (excerpt)	30
Useful information	20	Warranty	30
Setting up	20		

## Declaration of conformity

We declare under our sole responsibility that this product corresponds to the directives 2014/35/EU, 2006/42/EC, 2014/30/EU and 2011/65/EU and conforms with the following standards or normative documents: EN 61010-1, EN 61010-2-051, EN 61326-1, EN 60529 and EN ISO 12100.

## Explication of warning symbols

 <b>DANGER</b>	Indicates an imminently hazardous situation, which, if not avoided, will result in death, serious injury.
 <b>WARNING</b>	Indicates a potentially hazardous situation, which, if not avoided, can result in death, serious injury.
 <b>CAUTION</b>	Indicates a potentially hazardous situation, which, if not avoided, can result in injury.
 <b>NOTICE</b>	Indicates practices which, if not avoided, can result in equipment damage.
 <b>DANGER</b>	Danger - note on hazards arising from a hot surface.
 <b>CAUTION</b>	Indicates crushing risk to fingers/hand.

## Safety instructions

### ***For your protection***

- **Read the operating instructions completely before starting up and follow the safety instructions.**
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the device.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- **When working under a vacuum in particular!**
- Wear your personal protective equipment in accordance with the hazard category of the medium to be processed. There may be a risk of:

- Splashing liquids,
- Body parts, hair, clothing and jewelry getting caught,
- Injury as a result of glass breakage.

 **DANGER** Inhalation of or contact with media such as poisonous liquids, gases, spray mist, vapours, dusts or biological and microbiological materials can be hazardous to user.

- Set up the device in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- Ensure that there is sufficient space above the device as

- the glass assembly may exceed the height of the device.
- Prior to each use, always check the device, accessories and especially the glass parts for damage. Do not use damaged components.
- Ensure that the glass assembly is tension-free! Danger of cracking as a result of:
  - Stress due to incorrect assembly,
  - External mechanical hazards,
  - Local temperature peaks.
- Ensure that the stand does not start to move due to vibrations respectively unbalance.
- Beware of hazards due to:
  - Flammable materials,
  - Combustible media with a low boiling temperature.

 **CAUTION**

Only process and heat media that have a flash point higher than the adjusted safe temperature limit of the heating bath that has been set. The safe temperature limit of the heating bath must always be set to at least 25 °C lower than the fire point of the media used.

- Do not operate the device in explosive atmospheres, with hazardous substances or under water.
- Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
- Tasks with the device must only be performed when operation is monitored.
- Operation with excess pressure is not permitted (for cooling water pressure see "Technical Data").
- Do not cover the ventilation slots of the device in order to ensure adequate cooling of the drive.
- There may be electrostatic discharges between the medium and the drive which could pose a direct danger.
- The device is not suitable for manual operation (except lift movement).
- Safe operation is only guaranteed with the accessories described in the "Accessories" chapter.
- Refer to the operating instructions for the heating bath.
- Refer to the operating instructions for the accessories, e.g. vacuum pump.
- Position the positive pressure outlet of the vacuum pump under a fume hood.
- Only use the device under an all side-closed exhaust, or a comparable protective device.
- Adapt the quantity and the type of distill and to the size of the distillation equipment. The condenser must work properly. Monitor the cooling water flow rate at the condenser outlet.
- The glass equipment must always be ventilated when working under normal pressure (e.g. open outlet at condenser) in order to prevent a pressure build-up.
- Please note that dangerous concentrations of gases, vapors or particulate matter can escape through the outlet at the condenser. Take appropriate action to avoid this risk, for example, downstream cold traps, gas wash bottles or an effective extraction system.
- Evacuated glass vessels must not be heated on one side; the evaporation flask must rotate during the heating phase.

- The glassware is designed for operation under a vacuum of down to 2 mbar. The equipment must be evacuated prior to heating (see chapter "Commissioning"). The equipment must only be aired again after cooling. When carrying out vacuum distillation, uncondensed vapours must be condensed out or safely dissipated. If there is a risk that the distillation residue could disintegrate in the presence of oxygen, only inert gas must be admitted for stress relief.

 **CAUTION**

Avoid peroxide formation. Organic peroxides can accumulate in distillation and exhaust residues and explode while decomposing! Keep liquids that tend to form organic peroxides away from light, in particular from UV rays and check them prior to distillation and exhaust for the presence of peroxides. Any existing peroxides must be eliminated. Many organic compounds are prone to the formation of peroxides, e.g. dekalin, diethyl ether, dioxane, tetrahydrofuran, as well as unsaturated hydrocarbons, such as tetralin, diene, cumene and aldehydes, ketones and solutions of these substances.

 **DANGER**

The heating bath, tempering medium, evaporation flask and glass assembly can become hot during operation and remain so for a long time afterwards! Let the components cool off before continuing work with the device.

 **WARNING**

Avoid delayed boiling! Never heat the evaporation flask in the heating bath without switching on the rotary drive! Sudden foaming or exhaust gases indicate that flask content is beginning to decompose. Switch off heating immediately. Use the lifting mechanism to lift the evaporation flask out of the heating bath. Evacuate the danger zone and warn those in the surrounding area!

- The device does not automatically restart following any interruption to the power supply, which can avoid over boiling that could happen after power failure. Please make sure that the heating bath will not restart after power failure.

***For protection of the equipment***

- The voltage stated on the rating plate must correspond to the mains voltage.
- Socket must be earthed (protective ground contact).
- Removable parts must be refitted to the device to prevent the infiltration of foreign objects, liquids, etc.
- Protect the device and accessories from bumps and impacts.
- The device may only be opened by experts.

## Correct use

### • Use

Together with the accessories recommended by IKA®, the device is suitable for:

- Quick and gentle distillation of liquids,
- Evaporation of solutions and suspensions,
- Crystallization, synthesis or cleaning of fine chemicals,
- Drying of powder and granulate material,
- Recycling of solvents.

Mode of operation: Tabletop device

### • Range of use

- Laboratories
- Pharmacies
- Schools
- Universities

This device is suitable for use in all areas except:

- Residential areas;
- Areas that are connected directly to a low-voltage supply network that also supplies residential areas.

The safety of the user cannot be guaranteed:

- If the device is operated with accessories that are not supplied or recommended by the manufacturer;
- If the device is operated improperly or contrary to the manufacturer's specifications;
- If the device or the printed circuit board are modified by the third parties.

## Unpacking

### • Unpacking

- Unpack the device carefully;
- Any damage should be notified immediately to the shipping agent (post office, railway network or transport company).

### • Contents of package

	RV 3 eco Drive	Heating bath	Vertical glassware RV 10.1	Vertical glassware RV 10.10 coated	Serrated washer M5 1	Cylindrical M5x40 2	Ring spanner 3	Switching adapter 4	RV 10.8001 seal 5	HEXAGON DIN911 SW5 6	Screwed Plug D 7	Screwed-Joint cap 8	Tube clip D=9 mm 9	Base 10	Vapour tube 11	Power cable 12	Operating instructions	Warranty card
RV 3 eco V	1	1	1		4	4	1	1	1	1	1	2	2	1	1	2	1	1
RV 3 eco VC	1	1		1	4	4	1	1	1	1	1	2	2	1	1	2	1	1
RV 3 eco FLEX	1	1			4	4	1	1	1	1	1	2	2	1	1	2	1	1

## Useful information

Distillation is a thermal separating process for liquid compounds based on substance-specific, pressure-dependent boiling points through evaporation and subsequent condensation.

The boiling point temperature decreases with decreasing external pressure which means that work is usually done under reduced pressure. In this way the heating bath can be maintained at a constant temperature (e.g. 60 °C). Using the vacuum, the boiling point is set with a steam temperature of approx. 40 °C. The cooling water for the condensation condenser should not be warmer than 20 °C (60-40-20 rule).

A chemical resistant membrane pump with a vacuum controller should be used to create the vacuum. The pump is protected from solvent residue by the addition of a condensate flask and/or a vacuum separator.

Working with a jet pump to create a vacuum can only be recommended to a limited extent as the solvents may contaminate the environment when using these systems. Speed, temperature, flask size and system pressure all affect the evaporator air capacity.

The optimum capacity of the flow-through condenser is approx. 60%.


This corresponds to condensation on approx. 2/3 of the cooling coil. With larger capacities there is the risk that the uncondensed solvent vapor will be extracted.

The device is designed for operation with a cooling water supply system (e.g. laboratory thermostat), but can also be run off a water supply line. Please refer to the "Technical Data" for information on cooling water pressure, temperature constancy, and flow rate.


## Setting up

### RV 3 eco drive

#### • Mount the base to the lift

Apply the cylindrical M5x25 (2) (4x) and the serrated washer M5 (1) (4x) to mount the base (3) with the lift unit (4) after you unpack the package.

#### • Fix the switching adapter

Attach the switching adapter on the rear side and secure it with the brackets.

- Insert the switching adapter into the upper bracket. (2)
- Bend down the lower bracket and push the switching adapter into it. (3)
- Connect the low-voltage connector to the connection socket and fix it with the clamp. (4) (5)
- Plug the power cord in the switching adapter and connect it to the power supply of the valid voltage.

After connecting to the power supply, the power indicator will light up. (6)

#### • Fix the holding bracket (not included in delivery)

- Insert the condensate flask (not included in delivery) and attach the hose connectors to it

- Adjust the angle of the drive

- Remove the clamping device for the angle setting of the rotary drive by rotating the angle adjustable knob anticlockwise according to the indication mark (1).

- Set the desired angle for the rotary drive.
- Then, secure the rotary drive against rotation by tightening the angle adjustable knob clockwise.


## Glassware

### NOTICE


Please read the operating instructions of the glassware for the safe handling of laboratory glassware!

- Locking device


- There is a locking knob (1) for locking or unlocking the locking device (3).
- If a red mark (2) is shown, the locking device is unlocked, otherwise it is locked.
- To lock or unlock the locking device (3), push the locking knob (1) to the end position.

- Feed the vapour tube

- Unlock the locking device by turning it 60° anticlockwise according to the indication mark (1). The red mark should be visible.
- Feed the vapour tube (2) in until it stops.
- Then, lock the locking device by turning it clockwise by 60°.
- Push the lock knob to the end position. The red mark should be covered and invisible.
- The vapour tube is not allowed to be pulled out!
- Check the correct axial locking device on the vapour tube.
- Keep the red mark invisible.

#### • Mount & remove the evaporation flask

**Turn clockwise for locking,  
Turn anticlockwise for unlocking.**


#### **Mount the evaporation flask:**

- Place the evaporation flask (4) on the vapour tube.
- Turn the plastic nut (2) with the clip (3) anticlockwise until the clip can be put on the flange.
- Now turn the plastic nut clockwise until the clamp is hard up against the flange. When doing so, hold the locking device (1) on the vapour tube.

#### **Remove the evaporation flask:**


- Hold the locking device and loosen the plastic nut by turning it anticlockwise. This releases a tightly clamped evaporation flask.
- Hold the evaporation flask at the flange and open the clip.
- Remove the evaporation flask.
- Check that the locking of vapour tube is still closed!

#### • Mount the condenser seal


- Place the RV 10.8001 condenser seal (1) in the condenser receptacle and fit the glassware to the device according to the assembly instructions.

#### • Mount & remove the condenser

#### **Mount the condenser:**

- Place the cap nut (2) into the condenser (1), and then put the spring ring to the flange of the condenser.
- Place the condenser on the rotary drive (3) and tighten the cap nut by hand.
- Fix the receiving flask and hose connectors as shown. See also operating manual of glassware.

#### **Remove the condenser:**

- Use the ring spanner provided to loosen union nuts that are tightly fitted.
- Loosen the union nut by turning anticlockwise.



• Mount the washer

Drainage washer

**NOTICE**

Pay attention to the correct position of the washer!


• Description of special condensers

(See products/accessories at [www.ika.com](http://www.ika.com))

- RV 10.3 Vertical-intensive condenser with manifold.  
This vertical-intensive condenser features a double jacket design for particularly efficient condensation. Also available with coating (RV 10.30)
- RV 10.4 Dry ice condenser  
Dry ice condenser is for distilling low-boiling solvents. Cooling by dry ice, no cooling water is required. Maximum condensation thanks to low temperatures. Also available with coating (RV 10.40)
- RV 10.5 Vertical-condenser with manifold and cut-off valve for reflux distillation.  
Also available with coating (RV 10.50)
- RV 10.6 Vertical-intensive condenser with manifold and cut-off valve for reflux distillation  
This vertical-intensive action condenser features a double jacket design for particularly efficient reflux distillation. Also available with coating (RV 10.60)

Hose system

- Connect the water hoses (**H<sub>2</sub>O**) to the condenser according to the counter-flow principle.
- Install the vacuum connections to the condenser, condensate flask, vacuum controller with valve and vacuum pump.
- Always connect the vacuum hose (**vac**) to the condenser at the highest point to minimize solvent losses during suctioning.
- Use standard laboratory vacuum hoses with an inside diameter of 8 mm and a wall thickness of 5 mm.

Heating bath

**CAUTION**

Refer to the chapter entitled "Commissioning" in the operating instructions of heating bath.

**Commissioning**

Switch on the device

The unit is ready for service when the mains plug has been plugged in.

Use the power switch to switch on or off the device. Adjust the speed by turning the speed adjusting knob.

The scale from 1 to 10 corresponds to a speed range of 20 to 300 rpm under nominal load.

#### NOTICE

Make sure that the speed adjusting knob is in the "0" position when switching on the power switch, otherwise a warning sound will be given.

### Fill the evaporation flask

**Manual filling:** You can fill the evaporation flask manually prior to creating the vacuum. The evaporation flask should not be filled more than half its volume.

**Automatic filling:** Prior to filling the evaporation flask, a vacuum controller is used to regulate the glass apparatus to the target pressure.

- Now fill the evaporation flask using the back feed line.

#### CAUTION

The maximum allowed load (evaporation flask and contents) is 3 kg.

- Due to the vacuum present, the solvent is suctioned into the evaporation flask. This enables you to keep solvent loss due to suctioning to a minimum.

### Set up the heating bath

#### CAUTION

Refer to the operating instructions of the heating bath!

- Move the drive to the bottom position and check the position of the heating bath in relation to the evaporation flask. When using larger evaporation flasks or depending on the angle of the rotary drive, you can move the heating bath to the right.
- Fill the heating bath with the tempering medium until the evaporation flask is surrounded by tempering medium to 2/3 of its volume.

#### CAUTION

Avoid creating waves.

#### CAUTION

Avoid stress on the glass due to different evaporation flask and heating bath temperatures when lowering the evaporation flask into the heating bath!

- Switch on the rotary drive and slowly increase the speed.

- Switch on the heating bath using the power switch.

### Adjust the drive height

To adjust the drive height, grasp the height adjusting handle and hold the trigger (1).

The lift is unlocked and you can move the drive up and down. Release the trigger, the lift is locked in the adjusted height.

### Adjust the lower end stop

#### CAUTION

Depending on the size of the evaporation flask, the setting angle of the rotary drive and the position of the heating bath and lift, the evaporation flask can be in contact with the bottom of the heating bath.

#### CAUTION


**Glass may break!**

Limit the lower lift position with the variable end stop. Turn the knob (1) anticlockwise to unlock the stopper, and then you may adjust the lower end stop by moving it up and down. Turn the knob clockwise, the stopper is locked in the adjusted lower end stop.

## Assemble the stand pole RV 3.1 (Accessories)

- Assemble the condenser locking device according to the diagram.
- Mount the plate (1) using the thumb screw (2).
- Put the support rod (4) on the plate (1) and attach it with the screw nut (3).
- Attach the rubber protector (5).
- Fasten the Velcro® band (6) to the support rod (4).
- Secure the vertical glassware with the Velcro® band (6).

No stand is required if the condenser has been properly fitted and secured by tightening the condenser union nut on the rotary drive. The purpose of the stand is solely to prevent the condenser from twisting.

### CAUTION

If the stand is fitted incorrectly, the strong forces at the Velcro strap can cause stresses within the glass, which in turn can cause damage to the glass condenser.

After the condenser has been fitted, the stand is attached to the rotary drive. Make sure that the condenser is installed parallel to the lift body.


## Maintenance and cleaning

The device is maintenance-free. It is only subject to the natural wear and tear of components and their statistical failure rate.

The seal on the glass condenser should be checked at regular intervals and replaced if necessary.

### Cleaning

Remove the device from the mains before cleaning.

Only use cleaning materials recommended by **IKA®**:

Dirt	Cleaning agent
Dyes	Isopropyl alcohol
Building materials	Water containing detergent/isopropyl alcohol
Cosmetics	Water containing detergent/isopropyl alcohol
Food	Water containing detergent
Fuels	Water containing detergent
Other materials	Please consult <b>IKA®</b>

Wear protective gloves during cleaning the devices. Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.

Do not allow moisture to get into the device when cleaning.

Before using another than the recommended method for cleaning or decontamination, the user must ascertain with **IKA®** that this method does not destroy the device.

### Ordering spare parts

When ordering spare parts, please give:

- Device type,
- Serial number, see rating plate,
- Position number and description of spare part, see **www.ika.com**.


### Repairs

**Please only send devices in for repair that have been cleaned and are free of materials which might present health hazards.**

For this, use the "**Decontamination Certificate**" form which you can obtain from **IKA®** or can download a version for printing from the **IKA®** website at **www.ika.com**.

If your device requires repair, return it in its original packaging. Storage packaging is not sufficient when sending the device - also use appropriate transport packaging.

## Accessories (see [www.ika.com](http://www.ika.com))

RV 10.1	NS 29/32 Vertical glassware (1)
RV 10.10	NS 29/32 Vertical glassware, coated (1)
RV 10.2	NS 29/32 Diagonal glassware (2)
RV 10.20	NS 29/32 Diagonal glassware, coated (2)

RV 10.3	Vertical-intensive condenser with manifold (1)
RV 10.30	Vertical-intensive condenser with manifold, coated (1)
RV 10.4	Dry ice condenser (2)
RV 10.40	Dry ice condenser, coated (2)
RV 10.5	Vertical-condenser with manifold and cut-off valve for reflux distillation (no picture)
RV 10.50	Vertical-condenser with manifold and cut-off valve for reflux distillation, coated (no picture)
RV 10.6	Vertical-intensive condenser with manifold and cut-off valve for reflux distillation (3)
RV 10.60	Vertical-intensive condenser with manifold and cut-off valve for reflux distillation, coated (3)

RV 10.70	NS 29/32 Vapour tube
RV 10.74	NS 29/32 Vapour tube, short

RV 10.80	NS 29/32 Evaporation flask 50 ml
RV 10.800	NS 29/32 Evaporation flask, coated 50 ml
RV 10.81	NS 29/32 Evaporation flask 100 ml
RV 10.810	NS 29/32 Evaporation flask, coated 100 ml
RV 10.82	NS 29/32 Evaporation flask 250 ml
RV 10.820	NS 29/32 Evaporation flask, coated 250 ml
RV 10.83	NS 29/32 Evaporation flask 500 ml
RV 10.830	NS 29/32 Evaporation flask, coated 500 ml
RV 10.84	NS 29/32 Evaporation flask 1000 ml
RV 10.840	NS 29/32 Evaporation flask, coated 1000 ml
RV 10.90	NS 24/32 Evaporation flask 50 ml
RV 10.91	NS 24/32 Evaporation flask 100 ml

RV 10.100	KS 35/20 Receiving flask	100 ml
RV 10.101	KS 35/20 Receiving flask	250 ml
RV 10.102	KS 35/20 Receiving flask	500 ml
RV 10.103	KS 35/20 Receiving flask	1000 ml
RV 10.200	KS 35/20 Receiving flask, coated	100 ml
RV 10.201	KS 35/20 Receiving flask, coated	250 ml
RV 10.202	KS 35/20 Receiving flask, coated	500 ml
RV 10.203	KS 35/20 Receiving flask, coated	1000 ml

RV 10.300	NS 29/32 Powder flask	500 ml
RV 10.301	NS 29/32 Powder flask	1000 ml

RV 10.400	NS 29/32 Evaporation cylinder	500 ml
-----------	-------------------------------	--------

RV 10.500	NS 29/32 Foam brake	
-----------	---------------------	--

RV 10.600	NS 29/32 Distillation spider with	6 distilling sleeves
RV 10.601	NS 29/32 Distillation spider with	12 distilling sleeves
RV 10.602	NS 29/32 Distillation spider with	20 distilling sleeves
RV 10.610	Distilling sleeve	20 ml

RV 10.606  
RV 10.607

NS 29/32 Distillation spider with 5 flasks 50 ml  
NS 29/32 Distillation spider with 5 flasks 100 ml

RV 10.8001 Seal

RV 3.1 RV 3.1 Stand pole

Bracket

Condensate flask

## Error messages

Any malfunctions during operation will be identified by the error sound.

Proceed as follows in such cases:

- Switch off device using the main switch,
- Carry out corrective measures,
- Restart device.

If the actions described fail to resolve the fault or the other kind of error occurred, then take one of the following steps:

- Contact the service department,
- Send the device for repair, including a short description of the fault.

<b>Error sound</b>	<b>Cause</b>	<b>Effect</b>	<b>Correction</b>
Continuously beeping	- The speed has not been adjusted to 0 rpm when starting the device - PCBA failure	Motor can't run	- Turn the speed adjusting knob and set the speed to 0 rpm, after which you may set the required speed value.
Alternatively beep 1 sec and stop 1 sec	- Overload - Motor blocked - Defective cable - PCBA failure	Motor stops running	- Reduce the volume in the evaporator - Contact the service department
Alternatively beep 0.5 sec and stop 1.5 sec	- Temperature inside the device is too high - PCBA failure	Motor stops running	- Switch off the device - Let the device cool down - Switch on the device - Contact the service department
Alternatively beep 1.5 sec and stop 0.5 sec	- Overvoltage - Undervoltage - PCBA failure - Desktop power supply failure	Motor stops running	- Use the original power supply to supply the device - Contact the service department

## Technical data

Operating voltage range	<b>Vac</b>	(100-240) ± 10%
Rated voltage	<b>Vac</b>	100...240
Frequency	<b>Hz</b>	50/60
Rated power without heating bath	<b>W</b>	40
Rated power (max.) of switching adaptor	<b>W</b>	90
Speed	<b>rpm</b>	20...300
Speed tolerance	Set speed: < 100 rpm Set speed: ≥ 100 rpm	<b>rpm</b> <b>%</b> ± 1 ± 1
Speed display		Plastic part + scale
Smooth start		yes
Lift		manual
Stroke	<b>mm</b>	150
Head angle adjustable		0°...60°
Cooling surface	<b>cm<sup>2</sup></b>	1500
Cooling water flow rate min.	<b>l/h</b>	30
Cooling water flow rate max.	<b>l/h</b>	100
Cooling water pressure max.	<b>bar</b>	1
Perm. On-time	<b>%</b>	100
Perm. ambient temperature	<b>°C</b>	5...40
Perm. relative humidity	<b>%</b>	80
Protection acc. to DIN EN 60529		IP 20
Protection class		I
Contamination level		2
Weight (no glassware; no heating bath)	<b>kg</b>	8
Dimensions (W x D x H)	<b>mm</b>	440 x 280 x 590
Operation at a terrestrial altitude	<b>m</b>	max. 2000

*Subject to technical changes!*

## Solvent table (excerpt)

Solvent	Formula	Pressure for boiling point 40 °C in mbar (For HB 10 approx. 60 °C)
Acetic acid	$C_2H_4O_2$	44
Acetone	$C_3H_6O$	556
Acetonitrile	$C_2H_3N$	226
N-Amyl alcohol	$C_5H_{12}O$	11
n-Pentanol	$C_5H_{10}O$	11
n-Butanol	$C_4H_{10}O$	25
tert. Butanol	$C_4H_{10}O$	130
2-Methyl-2-Propanol	$C_4H_{10}O$	130
Butylacetate	$C_6H_{12}O_2$	39
Chlorobenzene	$C_6H_5Cl$	36
Chloroform	$CHCl_3$	474
Cyclohexane	$C_6H_{12}$	235
Dichloromethane	$CH_2Cl_2$	atm. press.
Methylenechloride	$CH_2Cl_2$	atm. Press.
Diethylether	$C_4H_{10}O$	atm. press.
1,2,-Dichloroethylene (trans)	$C_2H_2Cl_2$	751
Diisopropylether	$C_6H_{14}O$	375
Dioxane	$C_4H_8O_2$	107
Dimethylformamide (DMF)	$C_3H_7NO$	11
Ethanol	$C_2H_6O$	175
Ethylacetate	$C_4H_8O_2$	240
Ethylmethylketone	$C_4H_8O$	243
Heptane	$C_7H_{16}$	120
Hexane	$C_6H_{14}$	335
Isopropyl alcohol	$C_3H_8O$	137
Isoamyl alcohol	$C_5H_{12}O$	14
3-Methyl-1-Butanol	$C_5H_{12}O$	14
Methanol	$CH_4O$	337
Pentane	$C_5H_{12}$	atm. press.
n-Propyl alcohol	$C_3H_8O$	67
Pentachloroethane	$C_2HCl_5$	13
1, 1, 2, 2, -Tetrachloroethane	$C_2H_2Cl_4$	35
1, 1, 1, -Trichloroethane	$C_2H_3Cl_3$	300
Tetrachloroethylene	$C_2Cl_4$	53
Tetrachloromethane	$CCl_4$	271
Tetrahydrofurane (THF)	$C_4H_8O$	357
Toluene	$C_7H_8$	77
Trichloroethylene	$C_2HCl_3$	183
Water	$H_2O$	72
Xylene	$C_8H_{10}$	25

## Warranty

In accordance with IKA® warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine directly to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

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**IKA® - Werke GmbH & Co.KG**