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
Instructions for Use Liqui-Sampler – Sampler for Liquids

Art. No. 5330-2060/-2100/-2200 (PTFE), 5330-1060/-1100/-1200 (Polypropylene)

Read carefully before use and keep for later reference.

Liqui-Sampler has been developed for the safe, simple and uncontaminated sampling of liquids, and is available in PTFE (Polytetrafluorethylen) or PP (polypropylene). Liqui-Sampler PTFE is entirely made of chemically inert PTFE and is therefore suitable for the uncontaminated sampling of almost any type of chemicals. Liqui-Sampler PP has been developed for industrial, non-analytical use. Both samplers are available in lengths of 60 cm, 100 cm and 200 cm (special lengths available on request). Their design enables easy assembly before use and easy dismantling for cleaning, which is particularly important for uncontaminated sampling.

Attention! When filling and emptying the sampler with aggressive/dangerous substances, it is imperative that the general safety regulations are followed, and that suitable protective gloves, protective goggles and protective clothing are worn. Do not use Liqui-Sampler® for the storage of liquids.

Liqui-Sampler can be universally used for point sampling, cross-sectional sampling and bottom sampling, see figures below.

Sampling

- A** - ensure that the sampler is absolutely clean
 - slide the white measuring ring (2) to the desired sampling depth (from the surface of the liquid).

Point Sampling

- B** - immerse the Liqui-Sampler in the liquid until the measuring ring²⁾ reaches the surface of the liquid
 - Now press the colored push-button¹⁾ – liquid enters the sampler up to the measuring ring.
 - release the push-button – the valve at the bottom of the sampler is sealed.

Cross-Sectional Sampling

- B** - Keep the push-button¹⁾ pressed down, immerse the Liqui-Sampler in the liquid until the measuring ring²⁾ reaches the surface of the liquid – the liquid enters the sampler up to the measuring ring.
 - release the push-button – the valve at the bottom of the sampler is sealed.

Bottom Sampling

- B** - immerse the Liqui-Sampler until it reaches the bottom of the container, or until the measuring ring²⁾ reaches the surface of the liquid.
 - now press the colored push-button¹⁾ liquid enters the sampler up to the measuring ring. (Please note that slurry at the bottom of the container can block the bottom valve!)
 - release the push-button – the valve at the bottom of the sampler is sealed.
- C** - withdraw the Liqui-Sampler and let surface liquid drip off. Using the blue tool, push the wiping ring²⁾ down, thus wiping off the liquid from the sampler (use protective gloves!).
 - Attention: Do not tilt the Liqui-Sampler when it is filled, otherwise the liquid can escape through the vent at the push-button.
 - hold the Liqui-Sampler in/over the open sample-holder³⁾, (the discharge funnel ensures exact targeting) and press the colored push-button¹⁾. The collected liquid sample flows into the sample-holder.
 - seal the sample-holder.
 - clean Liqui-Sampler (see Cleaning and Assembly).
 - after use, store Liqui-Sampler in its storage case or quiver (see accessory), otherwise the extremely chemical-resistant Teflon® tube could be deformed when stored on an uneven surface!

Cleaning and Assembly

Assembly is very simple: screw off the discharge funnel (3), hold the valve cone (4) and screw off the push-button (1). All the other parts can now be removed from the sampling tube (5). Assembly is performed in the reverse sequence, taking care that the plastic threads, particularly at the discharge funnel (3) are not forced (valve-seat (4) may otherwise be deformed).

When cleaning, use cleaning agents suitable for the sampled substance. Do not clean the sampler with abrasive substances. All parts can be sterilised (using steam up to 121 °C). A good cleaning-effect can be obtained by using the cleaning brush (see accessory).

Problem

- collecting tube does not fit in the opening of the container
- Fluid is dripping from the bottom valve

Solution

- remove discharge funnel (3) which has a Ø of 33 mm and take a sample without it (tube Ø = 25 mm). If necessary, remount it before discharging the liquid.
1. Dirt particles are in the bottom valve – screw off and clean.
 2. Discharge funnel has been screwed on too tightly, valve seat is deformed – unscrew discharge funnel.

Special Accessories

- | | | | |
|----|--------------------------|--------------------|---|
| 1) | Color coded push-buttons | Art.-No. 5330-0149 | Set with 5 buttons (blue, yellow, green, white, black) |
| | Cleaning brush | Art.-No. 5304-0100 | excellent for thoroughly cleaning of the sampling tube. |
| | Storage case | Art.-No. 5618-0100 | for safe and clean storage (wall-mounting) |
| 3) | Sample-holders | | we have a wide selection of sample-holders in all kinds of shapes (narrow-necked, wide-necked, round, square) and materials (Teflon, PFA, PP, aluminum), please ask for brochure. |

²⁾ The measuring ring is used both as a sampling depth mark, and as a wiping ring for the liquid.