Mettler-Toledo AG, Analytical, Sonnenbergstrasse 74, CH-8603 Schwerzenbach, Switzerland
Phone ++41 (0)44 806 77 11, Fax ++41 (0)44 806 73 50, Internet: www.mt.com

Printed on 100% chlorine-free paper, for the sake of our environment. 11/07
Gedruckt auf 100% chlorfreiem hergestelltem Papier. Unserer Umwelt zuliebe.
Imprimé sur papier 100% exempte de chlore. Par souci d'écologie.
Stampato al 100% su carta senza cloro, per la salvaguardia dell'ambiente.
Impreso en papel totalmente exento de cloro, para preservar nuestro entorno.
この説明書は、環境に配慮し無塩素漂白紙に印刷されています。
为保护我们的环境, 此纸张为100%无氯纸

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52 003 284

Read this first

The METTLER TOLEDO InLab® pH combination electrodes are intended for accurate pH measurement in the laboratory at ambient pressure. The exact specifications of the pH and temperature range are printed on each individual electrode shaft or can be found under www.mt.com/electrodes. These operating instructions cover the following electrode models:

51 343 031	InLab® Versatile Pro	51 343 162	InLab® Micro Pro
51 343 054	InLab® Routine Pro	51 343 171	InLab® Pure Pro
51 343 056	InLab® Routine Pt1000	51 343 061	InLab® Reach Pro
51 343 057	InLab® Routine Pro-L	51 343 062	InLab® Reach Pt1000
51 343 071	InLab® Science Pro	51 343 023	InLab® Basics Pt1000

Safety precautions

The electrode is to be used only for the purpose stated above. It is advisable to wear safety glasses when handling electrolytes and cleaning agents. The electrodes contain glass parts.

Getting ready

The electrode is delivered ready for use. Please check visually for broken parts. For electrodes with detachable cable connect electrode head with cable connector by twisting clockwise. For cables and electrodes with MultiPin™ connector both marks (white/gray) must be in line before twisting. Connect other end with the pH meter. The cable for the temperature sensor must also be connected.

Next remove the wetting cap and rinse with distilled water. Eliminate any air bubbles inside the glass membrane space (internal reference electrolyte) by smoothly shaking the unit in the vertical plane (as for a fever thermometer). The electrode is now ready for calibration. Please follow the detailed calibration procedure, described in the manual of your pH meter. After calibration the electrode is ready for measurement.

Practical hint: After rinsing do not rub the glass membrane, since this can increase response time or damage the electrode.

Storage Rinse or clean electrode after use and replace the wetting cap, partly filled (1/3) with reference electrolyte.

Refilling Open the SafeLock™. Refill the electrode with the recommended reference electrolyte, if necessary. The electrolyte level of a correctly filled electrode is approximately 10 mm below the filling hole.

Trouble-shooting

- check all connections;
- temperature sensor must be connected to the corresponding pH meter input;
- for tips go to www.electrodes.net

Slow response/drift

- clean diaphragm*;
- after dry storage soak for 24 hours in reference electrolyte.

Slope not adjustable

- clean diaphragm*;
- clean connector.

Zero not adjustable

- clean diaphragm*;
- change reference electrolyte.

* Contamination by ...

- Proteins – Use cleaning solution Pepsin/HCl (Order No. 51 340 068);
 - Sulfides – Use cleaning solution Thiourea (Order No. 51 340 070);
 - Oily/organic liquids – Rinse with acetone or ethanol (flammable);
 - Acid/alkaline-soluble – Rinse with 0.1 mol/L HCl or 0.1 mol/L NaOH.
- Each time the electrode is cleaned, it must then be recalibrated.

Disposal If no relative official regulations are in force, used or defective InLab® sensors and the packaging can be disposed of as ordinary household waste.

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