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## **SONOREX SUPER und SONOREX LONGLIFE**

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# Operating Instructions

## **SONOREX SUPER and SONOREX LONGLIFE**

### **ultrasonic cleaners for aqueous cleaning agents**

Cleaning units, accessories and cleaning agents have to be strictly used according to the instruction manual or the corresponding product information.

Keep cleaning units and cleaning agents away from children and from persons not instructed for correct use by means of these instructions.

In case of improper use, the manufacturer does not take any responsibility for safety and proper functioning.

Unauthorized changes/modifications will result in expiration of the CE-conformity.

If service is required, please contact your dealer or the address mentioned below.

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# Contents

1.	Strictly observe the following general instructions for safety and correct operation! .....	13
2.	Start .....	14
2.1.	Degassing .....	14
2.2.	Cleaning .....	14
2.3.	Ultrasound ON / OFF .....	16
2.4.	Units with heating .....	16
2.5.	Emptying the ultrasonic tank .....	17
3.	Additional recommendations .....	17
3.1.	Cleaning agents .....	17
3.2.	Cleaning and care .....	18
3.3.	Treatment of medical instruments .....	18
4.	Disfunctions .....	19
5.	Dosing table for <b>SONOREX</b> – compact ultrasonic cleaners .....	20
6.	Technical data .....	21 and 32

# 1. Strictly observe the following general instructions for safety and correct operation!

- Keep cleaning units and cleaning agents away from children and also away from persons not instructed for correct use by means of these instructions.
- Do not immerse parts of the body (hand, feet) nor creatures or plants into the ultrasonic tank. During ultrasonic cleaning do not dip your hands into the cleaning liquid.
- Avoid any contact of your eyes and/or your skin with the cleaning liquid.
- Do not ingest nor breathe in the cleaning liquids.
- Never run the unit without liquid!
- Do not fill hot water into the ultrasonic tank. The filling temperature should not exceed 50 °C.
- Do not put parts to be cleaned into the tank with aggressive soiling such as acids, chloride ions or similar substances!
- Do not fill the stainless steel oscillating tanks with aggressive cleaning agents such as acids or chemicals which contain or release chloride ions (some disinfectants, dish-washing detergents, household cleaners, saline solutions).
- Do not fill the stainless steel oscillating tanks with combustible liquids.
- Use aggressive cleaning agents only in inset beakers or non-perforated plastic insert tanks.
- After removing the packing, check the unit's casing and control elements for any signs of damage. In such case, do not connect the unit to the electrical power supply, but report the damage immediately and in writing to the carrier and your supplier. Save the original packing.
- Place the unit in a horizontal position on a dry, solid base.
- Plug the unit to grounded sockets only.
- Keep the unit's surface and control elements clean and dry!
- Replace defective parts only with original **SONOREX** parts!
- During preheating, it is necessary to stir the cleaning liquid at least every 15 minutes or switch on ultrasound (retardation of boiling – danger of scalding)!
- Ultrasonic cleaning produces cavitation noise. If working within distances of 2 meters from the unit for longer periods of time, be sure to wear ear plugs to avoid damage to your hearing!
- Do not place anything on the bottom of the oscillating tank. Use appropriate **SONOREX** accessories (insert basket, positioning lid with insert vessel, etc.).
- The ultrasonic energy in the oscillating tank heats up the cleaning liquid without additional heating! Check the temperature regularly when cleaning temperature-sensitive parts.
- Do not move or carry units filled with liquid.

## Principle of ultrasonic cleaning

PZT oscillators underneath the tank convert electrical energy into mechanical vibrations.

**SONOREX** units produce oscillations of 35 kHz in the cleaning liquid. Tiniest vacuum bubbles are formed which then implode (cavitation). This cavitation principle removes dirt and grease from the deepest pores - "electronic scrubbing" - at all points that come into contact with the cleaning liquid.

## 2. Start

- Make sure that the outlet of the unit – either ball cock or closing cap (RK 100 SH) - is closed!
- Fill the oscillating tank up to the filling mark.
- Connect the unit to a grounded socket.

### 2.1. Degassing

Degassing the cleaning liquid enhances the cleaning effect.


Cleaning and contact liquids contain dissolved gases (for example, oxygen). Sonicate cleaning and contact liquid, which have just been poured into the oscillating tank or have been left in it for a longer period of time, for approx. 5 to 15 minutes before starting the cleaning process. Degas oscillating tanks with a volume of more than 10 l for roughly 30 minutes. During the degassing process the cavitation noise changes. As the loud degassing noise decreases towards the end of the process, the unit appears to run more quietly. Lower noise levels indicate the end of the degassing process and an improved cleaning impact. They do not mean a reduction of ultrasonic energy.

### 2.2. Cleaning

- Be sure that the objects in the tank are completely covered with liquid.
- Normally the objects may be cleaned **directly** in the ultrasonic tank with an insert basket.
- In the following cases clean objects **indirectly** in inset beakers or non-perforated insert tanks to protect the stainless steel oscillating tank:
  - when chemically aggressive cleaning liquids (for example acids, etc.) are used.
  - when chemically aggressive substances (for example when cleaning racks from film developing machines) or
  - when abrasive substances (for example polishing paste, quartz, sand) are to be removed.

**DIRECTLY** in the oscillating tank:


- Fill the oscillating tank to 2/3 with cleaning liquid (TICKOPUR, STAMMOPUR see dosage on label). Improperly filled units can suffer damage (see point 5 – dosing table).
- Use always an insert basket or utensil holder.
- Parts must not be placed on the tank bottom (Exemption basket K6 and SH 7).

**INDIRECTLY** in inset beakers or non-perforated insert plastic tub


⇒ **SONOREX** Accessories

- Fill the oscillating tank with a tenside-containing contact liquid (for example TICKOPUR R 33 or STAMMOPUR RD 5 diluted with water) to achieve a good transmission of ultrasound.
- Fill the inset beakers or insert plastic tub with cleaning liquid (for example TICKOPUR, STAMMOPUR).
- Suspend inset beakers in the positioning cover, do not place them on the tank bottom.
- Inset beakers or insert plastic tub must be immersed in the contact liquid to a depth of at least 2 cm. See to the tank filling to 2/3.
- Several inset beakers with different cleaning liquids can be sonicated at the same time.

**When inserting the objects to be cleaned, be sure to observe the following:**

- Place the dirtier side towards the bottom.
- Do not stack - ultrasonic energy will be absorbed.
- Sensitive pieces may not be allowed to touch one another. Avoid damage by using plastic insets, plastic insert baskets or the knobbed silicone mats in insert basket.
- Be sure that there are no air bubbles in hollow spaces (for example pocket holes).
- Check the temperature when cleaning heat-sensitive objects.

## 2.3. Ultrasound ON / OFF

### Cleaning time


Cleaning time must be kept as short as possible. Depending on the degree of soiling and the cleaning liquid used (see the data on its label), the cleaning time varies between 1 and 15 minutes.

### Timer operation

- Turn the knob clockwise to the appropriate cleaning time → the green pilot lamp lights up.
- The timer switches off automatically.
- You can shorten the time of cleaning or switch the unit off by turning the knob anticlockwise.

### Continuous operation (∞ position)

- Turn the knob anticlockwise until it clicks into position → the green pilot lamp lights up.
- The unit does not switch off automatically, turn the knob back to “0” in order to switch off.

The unit may remain plugged in when it is not running.  
To disconnect it pull the mains plug.

### Important leads:

- Warmed-up liquids improve the cleaning efficiency.
- However, the effect of the ultrasonic cavitation decreases at higher temperatures<sup>1</sup>. Experience has shown that best cleaning results are achieved by using an aqueous solution with a bath temperature of 50 to 60 °C.
- Ultrasonic energy also heats the cleaning agent. Continuous irradiation and/or the covering of the oscillating tank may cause the temperature of the liquid to rise more quickly and exceed the thermostat's default setting.
- **Do not** heat disinfectants.
- When reaching a critical temperature, the ultrasonic power is reduced to protect the electronic components inside the ultrasonic device. This way, the internal temperature is prevented from further rising.


## 2.4. Units with heating

Units marked with „H” have a built-in heating.


The heating system runs independently from the ultrasonic energy.

- Tip switch operation (fixed temperature of 65 °C):  
Heating switch on – control lamp is on.  
When the bath temperature of approx. 65 °C is reached, the heating will switch off, the control lamp remains on.
- Set the temperature knob clockwise to the selected temperature → the thermostatic regulation is in operation:
  - the yellow and white pilot lamps light up.
  - the yellow pilot lamp extinguishes when the selected temperature is reached.

RK 31 H: Tip switch

RK 52 H/CH and upwards:  
temperature knob

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
Do not place the unit into a sink.

### *Units with drain*

- Ball cock  
Turn the lever to the outflow direction.


- Outlet hose, closing cap

When unscrewing the closing cap, make sure that the hose end is in a higher position than the liquid level in the oscillating tank.

### *Units without drain*

Pour out the liquid over the rear left tank corner. Avoid splashing from below and into the socket.

## 3. Recommendations

### 3.1. Cleaning agents

- Only use approved and especially developed cleaning agents, e.g. TICKOPUR, TICKOMED, STAMMOPUR in order to obtain a high cleaning efficiency. Tap water without any additives does not clean sufficiently.
- When using the cleaning agents you must always strictly observe the safety instructions in the product information.
- Keep cleaning agents away from children and from persons not instructed for correct use by means of the product information.
- Do neither ingest nor inhale the cleaning agents. Avoid any contact of your eyes and/or skin with the cleaning agents.
- Replace used liquids and disinfectants, do not “freshen them up” with an added dose.
- Layers of dirt on the tank bottom or in the insets (abrasives and polishing agents, for example, as well as other deposits) reduce the cleaning efficiency. Empty the oscillating tank and clean it regularly. Replace the cleaning liquid.
- The cleaning and contact liquids in the ultrasonic oscillating tank may not exceed a maximum operating temperature of 100 °C.
- Distilled or deionized water without additives is unsuitable for ultrasonic cleaning directly in the oscillating tank. They lack the necessary cavitation properties to form a uniform ultrasonic field. Irregular cavitation can damage the oscillating tank or the objects being cleaned. (You may use distilled or deionized water in the inset beakers and non-perforated insert tubs.)
- Be careful when working with aggressive cleaning agents in insets or non-perforated insert tanks; avoid splashing into the contact liquid or onto stainless steel surfaces; if necessary, replace the contact liquid immediately, clean the surfaces and wipe them dry.

**The manufacturer does not take any responsibility for damage of equipment or damages at the parts to be cleaned caused by inappropriate cleaning chemicals.**

## 3.2. Cleaning and care

- Be sure to pull the plug before cleaning the unit.
- Do not shower the units or immerse them in water; instead, wipe them with a damp cloth, in- and outside.
- Rinse out the stainless steel tank frequently and wipe it dry.
- Do not use steel wool, scrapers or graters.
- If rings should appear in the stainless steel tank after a longer period of operation, remove them with a customary stainless steel cleaner without scouring additives.
- Metal parts remaining on the stainless steel surface, as well as particles of rust, can penetrate the passive protective layer of the stainless steel, thereby “activating” it and causing it to rust. This rust provokes pittings on the stainless steel surface. Therefore, remove all remaining metal parts such as screws, metal chips and the like; remove minor rust spots immediately with a soft cloth and customary stainless steel cleaners without scouring additives.

**SONOREX units do not require maintenance.**

## 3.3. Treatment of medical instruments

- If the ultrasonic cleaning unit is used in connection with disinfectants (medical device, class IIa) in order to support or to accelerate disinfection, then the ultrasonic unit is considered an accessory for a medical device and must be classified itself as medical device, class II!
- When disinfecting and cleaning contaminated medical instruments in an ultrasonic bath, the hygienic safety is important after the treatment. In case of improper and irregular disinfection and cleaning a micro-biological contamination is possible caused by settlement of microorganism, especially at the tank rim or in the drain section. Cross infection might be a result. Therefore, the tank and housing surfaces must be regularly disinfected and cleaned. According to the hygienic plan the user must carry out regular disinfection and cleaning with a surface disinfecting agent certified by VAH/DGHM\*, or another efficiently classified agent.
- The instruments may not be placed on the tank bottom. The basket avoids damages at the instruments and at the tank bottom. Do not staple instruments, overloading reduces the cleaning efficiency. Pliers and scissors must be opened completely or detached, if necessary. Instruments must be completely covered with cleaning liquid, air has to escape from cavities and hoses.

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Remark: Supplementary application information has to be respected.

## 4. Disfunctions

- Oscillations are too weak and non-uniform, the oscillation noise is too loud; the cleaning effect is insufficient:
  - Has the liquid been properly degassed?
    - ⇒ sonicate for 15 minutes.
  - Is the unit overloaded?
    - ⇒ remove some parts.
  - Non-uniform noise (wobbling) is not a defect
    - ⇒ change the fluid level.
- Heating defect:  
No objection to operate the unit without the heating.
- Minor erosion marks at the tank bottom?
  - ⇒ common wear and tear effect.  
The unit is faultless.
- The mains fuse can be released.
  - ⇒ The unit is equipped with an overvoltage arrestor which protects it from voltage peaks; the mains fuse can be released during overcharge.
- Repairs should only be carried out by authorised service.
- Replace defective parts only by original **SONOREX** parts.

Please inform us about disfunctions in writing.

## 5. Dosing table for SONOREX – compact ultrasonic cleaners

You can order this dosing table (DIN A 4-size) free of charge or download it as PDF-file from the following internet address::

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Values have been rounded up mathematically.

number printed in normal type: water

number printed in **bold** type: concentrate

Fill the oscillating tank to 1/3 with water and add the appropriate quantity of concentrate. Then add water up to the filling mark and sonicate for approx. 15 minutes for degassing.

or

According to the dosing table, fill the oscillating tank with water, add the appropriate quantity of concentrate, mix the solution for a short while and sonicate for approx. 15 minutes for degassing.

### Calculation of other volumes

Example:

- 10 litres solution ready for use
- 2,5 % dosage of the concentrate

$$\frac{10 \text{ l} \times 2,5 \%}{100 \%} = 0,25 \text{ l concentrate}$$

$$10 \text{ l} - 0,25 \text{ l} = 9,75 \text{ l water}$$

## 6. Technical data

**SONOREX** compact units are RFI-proofed and C € marked.

Compliance with the limits according to EN 55011, EN 61000-6-1/08.2002 and EN 60601-1-2.

Nominal voltage: 230V~ 50/60 Hz, (115 V on request)  
cable length 2 m  
High frequency: 35 kHz  
Protection class: IP 32  
Oscillating tank: stainless steel

### **Environmental conditions in accordance with EN 61 010-1 (IEC 1010-1)**

Degree of soiling: 2 according to IEC 60664-1:11.02  
Excess voltage category: II  
Permissible ambient temperature: 5 - 40 °C  
Permissible relative humidity up to 31 °C: 80 %  
Permissible relative humidity up to 40 °C: 50 %  
Dew not permitted.

### **Information for the use as a medical product**

Denotation: ultrasonic cleaning device  
UMDNS-nomenclature (ECRI): 14-263  
Purpose determination: ultrasonic-disinfection and cleaning of medical and dental instruments  
Classification (Medical Devices Directive 93/42/EEC, appendix IX): class I; active, not invasive, not implantable medical device  
Type, model, serial number, year of manufacture: see type designation plate on rear of device

### **Information according to IEC 601-1:**

Protection class: I  
Protection type (B, BF, CF): not applicable, no application part available

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