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Vielen Dank

General Information

This high quality product has been designed and developed to require a minimum of maintenance. Please follow the guidelines below and read this manual in its entirety to ensure safe operation of the unit.

Warning

- Please always disconnect mains plug before cleaning the device.
- For indoor use only.

Safety Information

- To avoid eye damage, don't look directly into the light source when turned on.

- Versions with UV light release intensive ultraviolet radiation during operation. A major part of the light is invisible. The beams can damage your health. Protect your eyes and skin during operation.

Assessment of the risks posed by optical radiation

The assessment of eye safety by optical radiation is carried out in accordance with the standard DIN EN 62471 (VDE 0837-471:2009-03): Photobiological safety of lamps and lamp systems. In this system, the UV LED specified in this data sheet is assigned to risk group 3: lamps/luminaires pose a hazard even in the case of fleeting or short-term irradiation. Protect your eyes during operation with UV protection goggles and/or a UV protection filter.

Package Contents

Cold Light Source LED.70	1
Power Cord	1
User Manual	1

Function and Operation

Technical Data

Luminous flux (light guide 1 m × Ø 9 mm)	approx. 1.500 lm
Color temperature	6.900 K
Color rendering index (CRI)	≥ 75
LED life span	50.000 h (L70, according to manufacturer)
Light guide access	15 mm (Standard)
Risk group	2
Power consumption	approx. 100 W
Operating voltage	100 - 240 V AC, 50/60 Hz
Protection class	I
IP code	IP20
Dimensions (l × w × h)	283 mm × 230 mm × 116 mm
Material	anodized aluminium
Color	black and silver
Weight	approx. 5,2 kg

Versions

LED.70 [EU]	<i>EU connector</i>	100-011724
LED.70 [US]	<i>US connector</i>	100-011842
LED.70 [UK]	<i>UK connector</i>	100-011843
LED.70 [CH]	<i>CH connector</i>	100-011844

Getting Started

The LED lamp is a plug & play device very easy to use.


Please check before use whether...

...the power cord is connected.


...the device is switched on (the yellow LED light glows).

...the brightness is set to the desired value.

Description of Buttons and Switches

- a) **Light Guide Access**
Connection of the light guide with a standard diameter of 15mm
- b) **On/Off Switch and Brightness Control**
For the adjustment of the light (standby until 100%)
- c) **Operation Indicator**
Yellow LED lights when the device is turned on
- d) **Main Switch**
- e) **Connector for Power Cord**
- f) **Assembly of Handle**
Optionally available

Technical Assistance

Troubleshooting

The information in the table below can help to resolve the possible operating problems quickly.

Problem	Causes	Solution
No function at the device	No network connection	Please check whether the power adaptor is connected.
	Electronic defective	Please contact your distributor.
	LED defective	Please contact your distributor.
LED gets dark	Thermal overload of the system	Please ensure an ambient temperature from 5° to 35°C.
	LED life span exceeded	Please contact your distributor.

Storage, Maintenance, Repair

The device is maintenance-free. Some points, however, you should take into consideration.

1. Please unplug the device before cleaning.
2. For the cleaning of the housing please use a smooth, semi-wet cloth, if necessary with a smooth household cleaning fluid.
3. No moisture should entering into the housing interior.
4. To avoid damages during storage or transport please retain the packaging.

Repairs are only allowed by trained and experienced personnel. In case of improper handling the warranty expires.

Technical Service

Web Resources

Visit our StarLight website at starlight-oe.com for:

- Contact details of technical customer service.
- StarLight product overview and price list as well as infos about accessories and spare parts.
- Additional product information and special offers.

Contact

For information or technical assistance contact your local distributor or directly StarLight Opto-Electronics GmbH & Co. KG.

Warranty

StarLight Opto-Electronics GmbH & Co. KG warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of purchase. If a defect is present, StarLight will, at its option, repair, replace, or refund the purchase price of this product at no charge to you. Precondition is the device is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, or misapplication, or from ordinary wear and tear.

For your protection, items being returned must be insured against possible damage or loss. This warranty shall be limited to the replacement of defective products.

It is expressly agreed that this warranty will be in lieu of all warranties of fitness and in lieu of the warranty of merchantability.

Equipment Disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste. Instead, it's your responsibility to correctly dispose of your equipment at lifecycle end by handling it over to an authorized facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.

For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

Thank you

Notizen / Notes