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pUC19-Marker *ready-to-use* (pUC19/Msp I)

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Instructions for use

pUC19 marker *ready-to-use* (pUC19/Msp I)

The marker was manufactured using pUC19-DNA and restriction endonuclease *Msp* I. Made from natural plasmids: fragment sizes can be validated!

DNA-ladder for analysis of small DNA fragments. The concentration of the individual fragments is equimolar so that rough quantification of the analysed DNA is possible.

The ready-to-use pUC19-marker has already been dissolved in 1x concentrated gel loading buffer ROTI®Load DNA (with glycerol, Art. No. 0100.1) and can be used directly. It has a concentration of 0.1 µg DNA/µl marker-solution.

The following fragment sizes (in bp) are included:
501, 489, 404, 331, 242, 190, 147, 111, 110, 67, 34, 34, 26.
pUC19 has a total length of 2686 bp.

Content:
pUC19 marker *ready-to-use*, 0.1 µg/µl

Application:
Our pUC19 marker *ready-to-use* is diluted in 1x ROTI®Load DNA (with glycerol), Art.No.: 0100.1, and can be used directly. The concentration of the marker solution is 1 µg DNA /10 µl.

<u>1x ROTI®Load DNA (with glycerol):</u>	
Tris/HCL, pH 7.5	10 mM
Sodium acetate	5 mM
EDTA	4 mM
Glycerol	8 % (w/v)
Bromophenol blue	0.02 % (w/v)
Xylene cyanol blue	0.015 % (w/v)
Bacteriocidal stabiliser	0.02 % /w/v)

Sample application/concentration:
Standard loading for mini to midi gels per lane is:

- with bands visible in UV-light after Et.-Br.-staining: 0.5 – 1.0 µg
- with detection after Et.-Br.-staining with signal-integrated camera systems: 0.1– 0.5 µg

Storage:

Optimal storage temperature is -20 °C. Please avoid repeated (>10 times) thawing and freezing. We recommend freezing in aliquots. Marker aliquots currently in use can be stored for short term (some weeks) at 4 °C (or at room temperature).


Fig.: 1 % agarose (NEEO, Art. No. 2267) + 2 % Synergel™ (Art. No. 0184)

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe
P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0
Fax: +49 (0) 721/ 5606-149
info@carlroth.com • www.carlroth.com ip 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

pUC19 marker *ready-to-use*

X901.1 500 µl
X901.2 4 x 500 µl