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
pBR322-Plasmid-DNA

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
Instructions for use

pBR322-Plasmid-DNA

High-purity plasmid DNA from pBR322 (Acc. No.: V01119) with a total length of 4361 base pairs.

OD _{260/280} :	>1,70
genomic. DNA:	<2 %
oc-Form:	<3 %
ccc-Form:	>95 %
base pairs:	4361


1,2 % agarose gel

Design: (gene name)

Tetracyclin resistance (*tet*): *tet* / 82-1276
Ampicillin resistance (*bla*): *amp* / 4153-3293
Origin of replication (*ori*): *ori* / 3122-2534

The coordinates are listed clockwise in the plasmid map.

The circular sequence was numbered in such a way that 0 lies in the middle of the singular EcoRI-restriction site. The copy number of pBR322 is limited by the gene product of *rop* to 20 copies per cell. However, by adding chloroamphenicol (final concentration 170 µg/ml) to log-phase cultures and then incubating for 8 hours, the copy number can be increased considerably. The total sequence can be obtained from the EMBLdata bank (ebi.ac.uk/embl/) under accession number (L01119).

The **pBR322-plasmid-DNA** was produced in host bacteria low in nuclease, and isolated through ion-exchange chromatography. It is available at an extremely high percentage in "supercoiled" form (ccc-form >95 %) and being a high-purity DNA, (OD_{260/280} >1.70) it is ideal for different applications:

- controlling transformation efficiency
- comparison with plasmid-ccc-forms of unknown size
- cloning DNA-fragments

- restriction (production of DNA-linear markers)
- calibrating DNA-precipitating or DNA-preparing procedures or equipment
- blocking of DNA-chips with plasmid-DNA of a defined sequence.

DNA is delivered lyophilised.

pBR322-Plasmid-DNA 50 µg X912.1

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe
P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0
Fax: +49 (0) 721/ 5606-149
info@carlroth.com • www.carlroth.com ip 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.