1 kbp DNA-Leiter

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1 kbp DNA-ladder

1 kbp DNA-ladder (Y014) with terminal sticky *EcoRI*-restriction sites. Deproteinised and lyophilised.

Universal marker for electrophoretic separation of fragments between 500 bp and 10 kb, optimally suited for size determination of inserts, plasmids, and restricted DNA. While a large running distance is necessary to separate larger bands in most markers, in this 1 kbp DNA-ladder the mass of the larger fragments was reduced. This allows obtaining clear and distinct bands with short running distances (app. 70 mm in 1.2 % agarose gels) even for the large marker fragments. The terminal sticky *EcoRI*-restriction sites allow simple labelling of the marker fragments via end-labelling techniques.

The following fragment sizes (in kbp) are included:
10, 8, 6, 5, 4, 3, 2.5, 2, 1.5, 1, 0.5

Content:

Y014.1 - 1 kbp DNA-ladder (50 µg, lyophilised)
0100.1 - 1 x sample buffer (ROTI®Load with glycerol) (500 µl) (Art. No. 0100.1)
Y014.2 - 4 x Y014.1

Application:

Depending on the application, the marker can be solubilised directly in the supplied 1x sample buffer ROTI®Load with glycerol (Art. No. 0100.1) or in TE-buffer (ROTI®Stock 100 x TE, Art. No. 1052.1). The lyophilized DNA is dissolved for 15 minutes at room temperature in an appropriate volume (see below) of sample buffer under occasional stirring.

Tip: For direct detection of the marker DNA in the gel under UV- or blue light we recommend to solubilize the marker DNA in 415 µl sterile water or 1xTE buffer. Then add 85 µl ROTI®Load DNASTAIN 1 (6x conc., Art. No. 5783.1). ROTI®Load DNASTAIN contains a fluorescent dye which may be excited by UV- and blue light, emitting at 525 nm.

Sample application/concentration:

Standard loading for mini to midi gels per lane is:

- with bands visible in UV-light after ethidium bromide staining:
0.5 – 0.7 µg
- with detection after ethidium bromide staining with signal-integrated camera systems:
0.1– 0.5 µg

Storage:

Optimal storage temperature is -20 °C.
Please avoid repeated thawing and freezing.
We recommend storage of aliquots.


Fig.: 1.2 % agarose (NEEO, Art. No. 2267). The 1 kbp DNA-ladder is highly compatible with the 100 bp DNA-ladders (equimolar: T834.1 or equalised: T833.1), resulting in a broad and reliable range of DNA-marker bands for all standard applications in the laboratory. An example of a mixture of both markers (0.4 µg each in a 1.2 % agarose gel) is depicted in figure B.

1 kbp DNA-ladder

Y014.1 50 µg + sample buffer
Y014.2 4 x 50 µg + sample buffer

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe
P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0
Fax: +49 (0) 721/ 5606-149
info@carlroth.com • www.carlroth.com sse 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.