## Anwendungen:

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## Ligase T4

5 U/μl, for molecular biology

The ligase T4 is an ATP-dependent recombinant enzyme isolated from *Escherichia coli*. The ligase is used for molecular cloning, site-directed mutagenesis, nick repair in duplex DNA, RNA or DNA/RNA hybrids and ligation-mediated PCR. T4 DNA ligase catalyses the formation of a phosphodiester bond between adjacent 5' phosphate and 3' hydroxyl end groups in duplex DNA or RNA.

**Unit Definition:** A (Weiss) unit of Ligase T4 catalyses the conversion of 1 nmol <sup>32</sup>P from pyrophosphate to norite adsorbable material in 30 minutes at 37 °C.

**Storage buffer:** 10 mM Tris-HCl (pH 7.4), 0.1 mM EDTA, 50 mM KCl, 1 mM DTT, 50% (v/v) glycerol

**10x T4 ligation buffer:** 660 mM Tris-HCl (pH 7.5), 100 mM DTT, 100 mM MgCl<sub>2</sub>  
The 10x T4 ligation buffer does not contain ATP. This must be added separately. The final concentration of ATP in the reaction should be 0.5 - 1.0 mM.

### Applications:

- Cloning of PCR fragments
- Cloning of restriction fragments
- Self-circulation of linear DNA
- Site-directed mutagenesis
- Nick repair in duplex DNA, RNA or DNA/RNA hybrids

### Protocol:

- Pipette the following reagents into a sterile nuclease-free reaction tube (on ice)

Component	Volume
Vector DNA	x μl (20 – 50 ng)
Insert DNA	y μl (3 -10 molar to vector)
10x T4 ligation buffer	2 μl
ATP solution (25 mM)	0,4 – 0,8 μl
T4 DNA ligase (5U/μl)	1 μl
Nuclease-free water	up to 20 μl

- Mix gently
- Incubate sticky ends for 0.5 - 1 hour at 20-25 °C and blunt ends for 1 - 2 hours at 20-25 °C
- Cool down on ice and transform 1 - 5 μl into 50 μl *E. coli* competent cells

### Further information:

- Thaw and resuspend 10x T4 ligation buffer and ATP solution at room temperature
- For blunt end ligations higher amounts of vector and insert DNA should be used

- The efficiency of the blunt end ligation is improved by the addition of PEG. The recommended PEG concentration is 5 % (w/v). It should be noted that a higher PEG concentration reduces the transformation efficiency.
- For sticky end ligation it is recommended to heat the vector and the insert DNA prior to ligation.
- It is recommended to use 3-10 times the amount of insert in relation to the vector. Calculation using the following formula:

$$\text{ng (insert)} = \frac{\text{ng (vector)} * \text{kb (insert)} * \text{molar ratio (insert:vector)}}{\text{kb (vector)}}$$

- Inactivate ligase at 65 °C (10 min) or 70 °C (5 min)
- The Kit contains:  
Ligase T4 (5 U/μl)  
Ligation buffer (10x)  
ATP solution (25 mM)

Ligase T4	500 U	plastic	3729.1
	2500 U	plastic	3729.2

### Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe  
P.O. Box 100121 • 76231 Karlsruhe  
Phone: +49 (0) 721/ 5606-0  
Fax: +49 (0) 721/ 5606-149  
info@carlroth.com • www.carlroth.com ip 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.