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Operating instructions for type DS450, Article No.: 11042500 / 11042502
Read these operating instructions carefully before installing and starting up the system!

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Preface

Dear Sir or Madam,

The stakpure team thanks you for the confidence you have placed in us.

With your selection of this mixed bed ion exchanger system you have chosen an innovative and high-quality product that will provide a long service life.

Please study the information that is given in these operating instructions on installation and operation before you start to install and start this water purification system.

Please also only appoint informed and trained personnel to operate the system and only use stakpure original accessories, replacement parts and consumables.

The above requests are particularly important as we, the manufacturer, cannot accept liability for any damages to the instrument, connected peripherals, building or persons that occur because of improper installation, operation or use of the system, or of the use of parts of any type of other manufacturers.

We wish you every success with and because of your new water purification system!

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
1 Table of content

2	Parts list	3
6	General notes and safety precautions	5
7	Mounting and assembling instructions	5
8	Note on the waste disposal of equipment	5
12	Change of cartridge	7
13	Technical Data	7
14	EEC – Declaration of conformity	7

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6 General notes and safety precautions

- Read these notes and precautions for your own well-being and proper functioning of the cartridge!
- The deionized water that is produced is not fit for drinking.
- Use only the original raw water hose, as the flow restrictor integrated in it prevents damage to the cartridge.
- **Never insert a shut-off device in the cartridge outlet nor kink the outlet hose, as the plastic cartridge could burst under the resulting pressure build up. Ensure free water flow at all times.**
- The use of this system in any other way than that described in these operating instructions invalidates the guarantee.
- An earthing contact socket outlet is required for the electrical connection.
- Under some circumstances, improper handling could result in the leakage of ion exchange resin from an ion exchange cartridge. **CAUTION: Danger of slipping** on leaked out resin!
Proceed as follows: As the resin material can cause slight irritation on contact with skin and/or eyes, immediately turn off the raw water tap to stop further leakage. On skin contact, immediately wash contacted skin areas with plenty of water. On contact with eyes, remove contact lenses if his is possible without difficulty, then immediately and carefully wash eyes with plenty of water for about 10 minutes. When irritation still occurs, seek medical advice immediately.
- Wear protective gloves to collect up the leaked out resin, fill it into a clearly labelled plastic container and return it to the authorized specialist supplier for waste disposal.

7 Mounting and assembling instructions

- 7.1 Draw the cartridge out from the bracket. Mount it near to the raw water tap using the enclosed screws and dowels. Ensure that the DS 450 mixed bed ion exchanger is freely accessible.

CAUTION: Before starting step 1, please make sure that there are no electrical wires or water-carrying pipes in the vicinity of the planned drilling holes.

- 7.2 After mounting the bracket, use slight pressure to insert the cartridge with attached conductivity meter in it.
- 7.3 Connect the raw water hose, free of kinks and coils, to the water tap and to the R 3/4" thread at the bottom of the cartridge.

CAUTION: The straight connector of the raw water hose contains the water restrictor and must be connected to the water tap. The angled connector is to be fitted to the R 3/4" thread of the cartridge.

- 7.4 Fit the pure water hose to the hose nipple of the conductivity meter.

CAUTION: Ensure that the pure water hose is neither blocked in any way nor kinked, as the plastic cartridge could then burst under the resulting build-up of pressure.

- 7.5 Connect the power plug to the appropriate socket with earth contact: 115V/230V, 50/60 Hz.
- 7.6 Carefully open the water tap and watch for any leaks. Should there be leakage at the union nut of the conductivity meter or at the hose connections, turn the water tap off and tighten the union nut or hose connections by hand until the water leak is stopped.
- 7.7 The system is now ready to produce demineralized water. Keep a watch on the conductivity meter, as effluent water is first to be taken for use when the pointer of the conductivity meter is in the green section or at a maximum of 20 µS/cm.

8 Note on the disposal of waste equipment

According to your state government requirements and the 2002/96 EC and 2006/66/EC directives, equipment that is to be scrapped can be brought to authorized collection points for recycling. Alternatively, it can be returned to us for proper recycling/waste disposal. In the

case of return for regeneration, repair, incorrect delivery or double delivery, please use the original cardboard box/packaging whenever possible. Send fall and knock protected.

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12 Change of cartridge

12.1 Close the water tap and remove the electric plug from the socket, then disconnect the raw water hose from the cartridge.

CAUTION: When opening, water may leak out of the cartridge. Make sure, that this does not impair any other systems in the vicinity.

12.2 Remove the cartridge with conductivity meter from the bracket.

12.3 Unscrew the conductivity meter from the exhausted cartridge.

12.4 Fit the conductivity meter on the fresh cartridge

and fit the combination in the bracket. Connect the raw water hose to the fresh cartridge. Only hand tighten the hose nipple.

12.5 Plug in the power plug. The system is now ready for use again. Please note: When no water flows through the cartridge, reionization effects could cause the conductivity meter to show incorrect values. Please do not use the water produced until the flow of it through the measuring cell brings the conductivity meter display to a correct value.

13 Technical Data DS 450

stakpure mixed-bed ion exchanger

Max. flow rate
Capacity at 10°GSG (total dissolved solids, 10° dH GSG=12.5° e)
Pure water quality
Max. water temperature

Electric connection

Height x diameter, incl. conductivity meter
Depth
Weight
Article no.

Replacement cartridge (2 pcs.)

Height x diameter

	DS 450	DS 450
Max. flow rate	50 l/h	50 l/h
Capacity at 10°GSG (total dissolved solids, 10° dH GSG=12.5° e)	425 l	425 l
Pure water quality	0.1 – 20 µS/cm	0.1 – 20 µS/cm
Max. water temperature	30° C	30° C
Electric connection	230 V / 50-60 Hz, 10 VA	115 V / 50-60 Hz, 10 VA
Height x diameter, incl. conductivity meter	100 x 600 mm	100 x 600 mm
Depth	115 mm	115 mm
Weight	3 kg	3 kg
Article no.	11042500	11042502
Replacement cartridge (2 pcs.)	11042700	11042700
Height x diameter	100 x 435 mm	100 x 435 mm


14 EEC – Declaration of conformity

Designation: Mix-bed ion exchanger
Type: Stakpure DS 450
Article No.: 11042500 / 11042502

Applicable EEC Directives:
EEC- Directive low voltage
EEC- Directive electromagnetic compatibility
EMV 2004/108/EG

Applied standard:
55011 und 61000

Niederahr, June 2014

Manufacturer/CE-Representative: Leo Trumm

stakpure GmbH

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D 56414 Niederahr

Telefon: 02602 10673-0

Telefax: 02602 10673-200

info@stakpure.de

www.stakpure.de