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Diaminobenzidin (DAB)

CN75

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**Substrate of horseradish peroxidase.
For detection of endogenous peroxidase and for immunohistochemical stainings.**

I. Introduction

Popular substrate of peroxidases. DAB is used for biochemical and histological detection of endogenous peroxidases like katalase and cytochrome oxidase. Additionally, it is the most common substrate for horseradish peroxidase (HRP) in histochemical and cytological assays. Suitable for staining of sections and Western-Blots. Also suitable for detection of activity of endogenous peroxidases in Mitochondria and in Peroxisomes.

II. Mechanism

DAB is oxidised by HRP, forming a brown precipitate, which is insoluble in aqueous and organic solvents. This precipitate can be detected in visible light and does not bleach during long-term storage.

III. Application

Reconstitution: 5 % in distilled water (add few drops of HCl if DAB does not dissolve completely.)

III.1 Detection of Peroxidase

Stock solution: 1 % (10 mg/ml) in 50 mM Tris buffer, or phosphate buffer (e.g. PBS), pH approx. 7.3

Working concentration: 0.05 - 0.1 % (0.5 - 1 mg/ml)

Incubation Buffer: Tris- or phosphate buffered systems like e.g. 0.05 M Tris-Cl or 0.1 M Na-Phosphate buffer, pH 7.2 bis 7.6, PBS, TBS
Optional: Add 50 – 100 mM imidazole

1. Add 3 ml DAB stock solution (1 %) to 27 ml incubation buffer (final conc. 0.1 %) and mix well.
2. Add 10 µl of 30 % fresh hydrogen peroxide and mix well.
3. Incubate sections or Blots in the staining solution up to the colour intensity required (approx. 5 - 15 mins.)
4. Wash in water or pure incubation buffer (3 x 1 min)
5. Counterstain if required
6. Embed

III.2 Detection of Catalases

Working concentration: 1 - 2 % (10 - 20 mg/ml)

Incubation Buffer: Highly alkaline buffer, like e. g. glycine-NaOH or 50 mM Veronal acetate or Tricine-buffer. pH 9 – 10.5

Note:

- DAB has to be dissolved in water. Solubilisation in buffer leads to cloudy solutions, presumably due to formation of the free base of DAB.
- The pH of the incubation buffer is essential! pH <7.0 diminishes staining intensity, pH >7.0 enhances background staining (peroxidase-staining).
- Do not titrate buffer solutions for DAB with NaOH, if DAB shall be used with blood. In the presence of NaOH, the iron of the hemoglobine starts oxidising DAB at once, resulting in a precipitating brown dye.
- Always use freshly prepared working solution. For oxidation high concentrations of fresh hydrogen peroxide are needed.
- Addition of Imidazole enhances sensitivity of the staining.
- Ultra-thin-sections and Western-Blots may be incubated directly in the staining mix. Thicker sections should be incubated for approx. 5 mins. in DAB solution prior to adding the hydrogen peroxide, in order to allow the DAB to diffuse to the target molecules.
- DAB is not suitable for detection of Peroxidase or Catalase of *Caenorhabditis elegans* (Togo *et al.*, 2000, *Eur. J. Biochem.* 267:1307-12).

IV. Signal enhancement

In Western-Blot systems as well as in sections, sensitivity and intensity of staining may be enhanced by a colour reaction of metal ions (Nickel, Copper, Cobalt, Silver) and Imidazole (e.g. Green *et al.*, 1989, *J. Clin. Pathol.* 42:875-80). We recommend enhancement by Nickel-, or Cobalt chloride combined with Imidazole:
0.05-0.1 % NiCl₂ in DAB/0.05 M Tris-Cl buffer, pH 7.6 or
0.02-0.06 % CoCl₂ in DAB/0.1 M Na-Phosphate buf. pH 7.3

- Addition of metal salt, mix well
- Addition of Imidazole ad 0.07 % (10 mM), mix well
- Addition of H₂O₂ ad 0.01 %, mix well
- Incubation up to the colour intensity required (approx. 5 - 15 mins.)

Note:

- Sections stained with DAB/metal ions may be dehydrated and embedded hydrophobously.
- In order to optimize staining, the amount of metal salts may be varied from 0.01 to 0.25 %. Nickel and Cobalt may also be added in parallel, in many cases further enhancing the contrast of stained tissue.
- Ultra-thin-sections and Western-Blots may be incubated directly in the staining mix. Thicker sections should be incubated for approx. 5 mins. in DAB / enhancer solution prior to adding the hydrogen peroxide, in order to allow the DAB to diffuse to the target molecules.

V. Storage

Lyophilisate: +4 °C

Storage of stock solution in aliquots at -20 °C

VI. Further Reagents

ROTI®Stock 10X PBS ready-to-use solution

Art. No. 1058

ROTI®Stock 10X TBS ready-to-use solution

Art. No. 1060

Hydrogen peroxide 30 % Art. No. 8070

Imidazole Art. No. X998

Nickel(II) chloride hexahydrate Art. No. 4489

Cobalt(II) chloride hexahydrate Art. No. T889

 **Warning** H315-H319-H335-H341-H351
P202-P261-P280-P305+P351+P338-P308+P313

Full text of hazard- and precautionary statements
see material safety data sheet section 2.2

Diaminobenzidin tetrahydrochloride

CN75.1	1 g
CN75.2	5 g
CN75.3	10 g