## ROTI® DipSlides PCA/Legionella

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# Product Data Sheet

## ROTI® DipSlides PCA/Legionella

DipSlides and contact agar in one versatile tool. Sterile nutrient broth for total germ count and for detection of legionellas

2493

- ✓ Ready-to-use and sterile
- ✓ Very easy application, reliable assay and easy interpretation
- ✓ Suitable for analysis of each surface and waters from air conditioners, water pipelines etc.
- ✓ For supervision of units for heating of drinking water (Drinking Water Regulations) and for surveillance of ventilation and air conditioning.

### COMPOSITION

PCA: Plate Count Agar (acc. to APHA, ISO 4833:2003) (X930) for total germ count (**amber**)

Legionella Agar: Legionella Agar (acc. to ISO 11731-2) (CL48) for selective detection of legionellas (**black**)

### STORAGE, EXPIRY AND SHIPMENT

ROTI® DipSlides are shipped in special boxes on gel ice. Please make sure that ROTI® DipSlides are stored immediately at 4-15 °C. Store boxes lid-up, in order to avoid moistening of the agar with condensed water.

ROTI® DipSlides may not be frozen. Slides that have been frozen have to be disposed of, since the medium is likely to have taken severe damage and growth of microorganisms is significantly hindered.

Store at 4-15 °C and protected from direct light. In order to avoid drying of the agar, temperature should vary for maximally ±2 °C. Expiry of ROTI® DipSlides is 5 months after production. Expiry date is noted on each package and on each vial with „EXP“.

### COMPOSITION OF AGARS (in g/l, approximated)

For further information, please see the respective data sheets of the agars.

Nutrient	Casein peptone 5.0, yeast extract 2.5, glucose 1.0, Agar 15.0, pH value 7.0 ± 0.2. <b>colour: amber</b>	Legionella Agar (CL48)	Yeast extract 10.0, charcoal activated 2.0 GVPC selective additive CL50 (complies), growth additive (CL49) (complies), agar 13.0, pH value 6.9 ± 0.2. <b>colour: black</b>
----------	--	------------------------	---

### IMPORTANT TO NOTE FOR SAFETY

**Inhalation of legionellas due to aerosols represents highest risk of infection.**

Appropriate precautions have to be implemented in order to prevent infection. Persons who carry out the assessment should have the ability, experience, instruction, information, training and resources to perform the assay competently and safely. In particular, they should know about potential sources and the risks of infection they present when working with *Legionella spp.* The current instructions-for-use gives only general information for test performances, which have to be adopted by each user for his or her special requirements, particularly if specific regulations have to be followed.

**At least note and implement the following precautions:**

- Be cautious when taking samples, particularly of moderately warm water. Prevent formation of aerosols!
- Although there is no evidence that legionellas escape the tubes during incubation, we recommend to not inhale the air coming from the incubator when opening the door. Thoroughly ventilate the room containing the incubator.
- Do not open the tubes for analysis of the agar (see EVALUATION OF RESULTS).
- Always wear gloves and mask (X458.1) when manipulating DipSlides.

### APPLICATION

- Remove carrier from the sterile vial. The screw cap serves automatically as holder. Make sure the agars don't show any sign of contamination or dehydration.
- Press agar for a few seconds onto the surface to be tested. Alternatively, immerse the agar for approx. 10 seconds in the liquid, in order to wet both sides completely. Carefully allow excess sample to run off. In case of samples of high viscosity, we recommend diluting the sample by 1:2 to 1:10. Dilution must be taken into account when calculating the sample contamination. Do not touch the agar (risk of contamination!)
- Replace carrier into the vial and close the cap moderately.
- Incubate vial with carrier at a temperature suitable for the bacteria (35 ± 2 °C for 10 days over all). During incubation (upright, standing on the lid), the vials should be screwed on only loosely, in order to enable entrance of oxygen. In accordance to official regulations or internal protocols other incubation parameters may be chosen.

5. PCA: After 48 hours of incubation count and register colonies.

Legionella Agar: After 3 days of incubation count and register colonies. Further analyses at, for instance, days 5, 7, and 10.

**Please note:** Many legionellas grow very slowly, forming only small colonies even after 10 days of incubation. Well detectable colonies appearing on days 2 or 3 often represent unspecific accompanying flora. Record and repeat the analysis.

Characteristic appearance of legionellas is a white or blue/grey colour with even rims (sometimes rose coloured), which look a bit like being made from porcelain or frosted glass. In case white/grey/bluish colonies grow after 7-10 days, professional microbiological analyses *have to be performed*.

## EVALUATION OF RESULTS

After the incubation time span required colonies are counted and noted as cfu (colony forming units). For this purpose, the carrier should not be removed from the vial, in order to prevent release of the microorganisms. After incubation, tightly close the tube and let it cool for at least 2 hours at 4 °C. Condensed water may then be tapped into the lid, enabling detailed analysis of the agars through the plastic wall.

Any validation required has to be undertaken by the user. We recommend to test each sample liquid or surface at least twice (with two ROTI®DipSlides), in order to validate results. While calculating sample contamination, eventual sample dilution factors have to be taken into account. Evaluation of results is done semi-quantitatively according to the table given below.

**For solutions the following approximation may be applied:**

**Cfu of one side x 40 is approx. cfu / ml solution. Detection limit is 10<sup>2</sup> CFU / ml solution.**

Table 1: Evaluation

Colony no.	cfu / ml *	Result	Interpretation
0	<10 <sup>2</sup>	Negative	Contamination below lower detection limits
1-5	ca. 10 <sup>2</sup>	Weakly positive	Weak contamination of solution/surface
10-50	ca. 10 <sup>3</sup>	Positive	Low contamination of solution/surface
100-500	ca. 10 <sup>4</sup>	Strongly positive	High contamination of solution/surface
>500	>10 <sup>4</sup>	Highly positive	Very high contamination of solution/surface

\* This calculation may be used as approximation to numbers only. For definite quantitation of contaminated solutions or surfaces, sophisticated microbiological procedures have to be performed using Roth nutrient dry media.

## DISPOSAL

Disposal of ROTI®DipSlides bearing cultures has to be done according to all common regulations for disposal of the respective microorganisms, and to regulations of good laboratory practice. E.g. autoclaving for 20 mins. at 121 °C or over night-incubation in appropriate disinfectants like, for instance, hypochlorite.

## MICROBIOLOGICAL TEST

**PCA** (32 ± 2 °C / 18-48 h and. \*: 30 °C / 72 ± 3 h, respectively)

Test Strain	Growth	Colour
* <i>Escherichia coli</i> ATCC 25922	Good	light-beige
<i>Escherichia coli</i> ATCC 8739	Good	light-beige
<i>Staphylococcus aureus</i> ATCC 25923	Good	light-beige
* <i>Staphylococcus aureus</i> ATCC 6538	Good	light-beige
<i>Staphylococcus epidermidis</i> ATCC 12228	Good	light-beige
* <i>Bacillus subtilis</i> ATCC 6633	Good	light-beige

**Legionella-Agar** (35 ± 2 °C / 24-72 h and up to 10 days)

Test strain	Growth	Colour
<i>Legionella pneumophila</i> ATCC 33153	Good	white
<i>Escherichia coli</i> ATCC 25922	Inhibited	-
<i>Staphylococcus epidermidis</i> ATCC 12228	Inhibited	-

Table 2: Overview

Ord. No.	Product	Agar Side 1	Agar Side 2	Application	Inc. temp.	Inc. time
2493.1	ROTI®DipSlide PCA/Legionella	Plate Count (X930)	Legionella Agar (CL48)	Whole germs + legionellas	35 ± 2 °C	10 days*

\* Recommended time points of analysis see above.

**ROTI®DipSlide PCA/Legionellas**

**20 slides**

**2493.1**

### Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe • P.O. Box 100121 • 76231 Karlsruhe  
 Phone: +49 (0) 721/ 5606-0 • Fax: +49 (0) 721/ 5606-149 • info@carlroth.com • www.carlroth.com

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

sse 07/2021