CRONOBACTER-CHROMOGENER-AGAR

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Product Data Sheet

CRONOBACTER CHROMOGENIC AGAR

Selective medium for detection of *Cronobacter* (formerly *Enterobacter*) *sakazakii* in milk powder and powdered infant formulae


ISO 22964:2017

2997

Formulation in g/l:

Caseinpeptone (Pancreas hydrolysate)	7.0
Sodium chloride	5.0
Yeast extract	3.0
Sodium desoxycholate	0.25
Chromogene Mixture	0.15
Sodium thiosulfate	1
Ammonium iron(III) citrate	1
Agar	15.0
Final pH	7.3 ± 0.2

Store at 2-8 °C

Cronobacter
(*Enterobacter*) *sakazakii*

PREPARATION

Suspend 32.4 g of the medium in 1 l distilled water. Solubilise under heating with frequent agitation and boil for approx. 1 min until complete dissolution. Autoclave for 15 mins. at 121 °C. Let cool to 45-50 °C, homogenise gently, and pour into Petri dishes. The prepared medium is amber in colour and should be stored in the dark at approx. 10 °C.

USES

Cronobacter Chromogenic-Agar is a selective medium for the detection of *Cronobacter* (*Enterobacter sakazakii*) in milk powder and powdered infant formula. The ISO normative 22964 recommends this medium for the isolation of *Cronobacter* spp.

Cronobacter is a Gram-negative, ubiquitous rod-bacterium which is more or less found everywhere. In healthy adults, its pathogenic relevance is marginal. It is, however, currently considered a putative pathogen in un-weaned babies, causing severe meningitis and necrotic enterocolitis with mortality rates of 40 to 80 %. Hence, it is vital to review the manufacturing process of dairy products particularly meant for infants, guaranteeing the absence of the bacterium in the final product.

After incubation in standard or enrichment broth, a loopful (approx. 10 µl) is streaked onto the surface of the Cronobacter Chromogenic Agar plate. The plate is incubated at 44 °C ± 1 °C for 24 h ± 2 h, followed by examination for the presence of presumptive *Cronobacter* by analysing pigment formation. The high incubation temperature of 44 °C efficiently inhibits Gram-positive bacteria.

The chromogenic mixture contains the chromogenic substrate α-D-glucopyranosid, that is cleaved by the bacterial α-D-glucosidase, resulting in a blue to greenish colouring of the *Cronobacter* colonies (see below). Typical colonies are small to medium sized green to blue-green colonies. Other enterobacteria grow in light to reddish colonies. Non-typical colonies are colourless to reddish/violet coloured and often slightly transparent.

Presumptive *Cronobacter* colonies should be confirmed on CASO Agar (Art. No. X937), on which they present a yellow colour, including biochemical confirmation.

MICROBIOLOGICAL TEST

The following results were obtained from type cultures after incubation at a temperature of 44 °C ± 1 °C and observed after 24 ± 2 hours. Temperatures have to be reviewed quite stringently: at ≥45 °C growth of *Cronobacter* is inhibited, while at ≤43 °C growth of the accompanying bacteria are not efficiently reduced.

Microorganisms	Growth	Colony Colour
<i>Cronobacter</i> / <i>Enterobacter sakazakii</i> ATCC 29544	Good	green-blue – turquoise/cyan
<i>Escherichia coli</i> ATCC 25922	Good	colourless/transparent – reddish/violet
<i>Escherichia coli</i> ATCC 8739	Good	colourless/transparent – reddish/violet
<i>Staphylococcus aureus</i> ATCC 25923	Inhibited	-
<i>Enterococcus faecalis</i> ATCC 29212	Inhibited	-

ISO normative 22964:2017 Horizontal method for detection of *Cronobacter* spp.

Guillaume G.O., et al. (2005) Cultural Method for Detection of *Enterobacter sakazakii* in environmental samples. *Journal of Food Protection*, 68:64-9.

CRONOBACTER CHROMOGENIC AGAR

100 g

2997.1

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe • P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0 • Fax: +49 (0) 721/ 5606-149 • info@carlroth.com • www.carlroth.com

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

see 06/2021