
UTI-CHROMOGENER AGAR

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Samra et al. (1998) Evaluation of use of a new chromogenic agar in detection of urinary tract pathogens. *J. Clin. Microbiol.* 36:990-4.

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Product Data Sheet


UTI CHROMOGENIC AGAR

Medium for detection and differentiation of microorganisms causing urinary tract infections

3089

Approximate formula in g/l:

Nutrients (mixture of extracts, glucose)	13.0
Peptone	16.0
Tryptophane	2.0
Chromogene Mixture	0.5
Agar	16.0
Final pH	7.2 ± 0.2

PREPARATION

Suspend 47.5 g of the medium in 1 l distilled water. Solubilise under heating with frequent agitation and boil for approx. 1 min. Autoclave for 15 mins. at 121 °C. Let cool to 45-50 °C, homogenise gently, and pour into Petri dishes. Store in the dark at approx. 10 °C. Colour: amber, slightly opalescent.

USES

UTI Chromogenic-Agar is a medium for the detection and differentiation of bacteria causing infection of the urinary tract. Due to the specially designed chromogenic mixture, the relevant species may be immediately differentiated and identified on one plate. Once incubation is completed, the colonies are analysed for pigment formation. The chromogenic mixture contains two major components: a) The chromogenic substrate X-β-glucosid, that is cleaved by the β-D-glucosidase found in *Enterobacteriaceae*, resulting in a blue colouring of *Enterobacter*, *Enterococcus*, and *Klebsiella* colonies (see below). b) A substrate cleaved by the β-galactosidase present in Coliforms, causing colonies of *E. coli*, *Enterobacter*, and *Klebsiella* to grow in a salmon to pink colour. Mixture of both colours causes *Enterobacter* and *Klebsiella* to form dark violet colonies. Due to the addition of tryptophane to the media formulation, followed by activation of bacterial tryptophane desaminase, also colonies of *Proteus spp.*, *Morganella spp.*, and *Providencia spp.* may be identified through their light-brown colouring. *Staphylococcus aureus* grows well in natural (cream) coloured colonies.

Please note: Microorganisms with atypical enzyme patterns may give anomalous reactions. For instance, 45 % of *Enterobacter cloacae* population do not contain a proper β-glucosidase, therefore resulting in pink colonies not distinguishable from *E. coli*. For confirmation of *E. coli* colonies we recommend performing an Indole tests.

MICROBIOLOGICAL TEST

The following results were obtained from type cultures after incubation at a temperature of 35 °C±2 °C and observed after 18-24 hours.

Microorganisms	Growth	Colony Colour
<i>Escherichia coli</i> ATCC 25922	Good	salmon-coloured - pink
<i>Enterobacter aerogenes</i> ATCC 13048	Good	dark blue - dark violet
<i>Klebsiella pneumoniae</i> ATCC 13883	Good	dark blue - dark violet
<i>Proteus mirabilis</i> ATCC 13315	Good	light brown – yellowish-brown
<i>Staphylococcus aureus</i> ATCC 25923	Good	cream-coloured (natural colony colour)
<i>Enterococcus faecalis</i> ATCC 19433	Good	light blue

Samra *et al.* (1998) Evaluation of use of a new chromogenic agar in detection of urinary tract pathogens. *J. Clin. Microbiol.* 36:990-4.

Store at 2-8 °C

UTI CHROMOGENIC AGAR

100 g

3089.1

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe • P.O. Box 100121 • 76231 Karlsruhe
 Phone: +49 (0) 721/ 5606-0 • Fax: +49 (0) 721/ 5606-149 • info@carlroth.com • www.carlroth.com

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

sse 07/2021