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CompactDry SL

4817

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CompactDry SL

4817

NISSUI ready-to-use, sterile
Ready-made plates with anhydrous chromogenic
nutrient medium for detection of salmonellae

A. Introduction

Compact Dry plates made by NISSUI are suitable for detection and enumeration of microorganisms from foods and feeds, cosmetics and pharmaceutical samples, as well as for hygiene monitoring and contamination control in food-producing areas. Due to the chromogenic additives, bacterial colonies depict specific colouring, making them easy to interpret. The nutrient medium used for Compact Dry SL plates is deduced from the formulation of Triple Sugar Iron Agar. Salmonellae may easily be identified by the green colouring of the colonies, accompanied by colour shift of the plate to yellow due to alcalization of the medium. Salmonellae producing hydrogen sulphide grow in black colonies on yellow medium. The presence of *Salmonella* in the sample is detected by the combination of three different test principles:

Carl Roth GmbH + Co. KG
Schoemperlenstraße 3-5
76185 Karlsruhe
Postfach 100121
76231 Karlsruhe
Telefon: +49 (0) 721/ 5606-0
Telefax: +49 (0) 721/ 5606-149
E-Mail: info@carlroth.de
Internet: www.carlroth.de

gh 01/2020

- Alkalisiation of the medium by *Salmonella*'s lysine decarboxylase ability (medium color will change blue-purple to yellow)
- Greening colonies caused by decomposition of chromogenic substrate with specific enzyme of *Salmonella*. Black colonies are generated by hydrogen sulfide producing *Salmonella*.
- Motility of *Salmonella*.

B. Storage and Usage

CompactDry plates are stable at room temperature and at 4 °C. Store protected from moisture and light. Do not use the product after its expiry date. Quality of the product is not warranted after being expired. Do not use the product that contains any foreign materials, discoloured or dehydrated, or its container is damaged. After opening the aluminum bag, any plates unused should be put back into the aluminum bag to be sealed with tape to avoid light and moisture, and use up as soon as possible.

C. Sample preparation

From waters or liquid foods

Add 9 times volume of Buffered Peptone Water to 1 vol. liquid specimen. Filtrate the liquid sample through membrane filter, and put the filter into Buffered Peptone Water.

From solid foods

Take 25 g of solid specimen into the sterilized homogenizer bag. Add 225 ml of sterile Buffered Peptone Water into the bag, and homogenize for about one (1) minute

From surfaces

Wir recommend using the CompactDry Swab (Art. No. 4809). Use the swab to wipe the surface, put into the device with wiping solution. Add 9 ml of Buffered Peptone Water to 1 ml of this solution.

D. Application

1. Incubate the samples in the closed homogenizer bag for 20-24 hours at 35-37 °C in the incubator for preenrichment.
2. Rub the bag for homogenization.
3. Open the plate and gently drop 0.1 ml (3 drops) of enriched specimen onto the dry sheet (approx. 1 cm far from the edge of plate). This enriched culture will stay at dropped point. Diffusion of this dropped specimen shall not reach to the edge of plate.
4. Gently drop 1 ml of sterilized water at the opposite point where the specimen dropped. The water diffuses automatically and evenly into the sheet and transforms the dried sheet into a gel within seconds.
5. Put the cap back onto the plate and note date and sample on the memorandum section.
6. Turn over the capped plate and incubate at appropriate temperature.
7. After incubation count the number of coloured colonies underneath the plate. White paper placed under the plate helps with counting.

Incubation:

Please use the incubation time/temperature according to the national food analysis recommended for total viable count.

General recommendation: 42 °C for 24 hours

E. Interpretation

Salmonella Positive

- 1.) Black to green isolated or fused colonies.
 - 2.) Colour of sheet around the colonies is changed to yellow.
 - 3.) Colonies have wandered all over the plate (motility).
- In case of all three parameters are positive, consider sample as *Salmonella* contaminated. If a large quantity of *Salmonella* is present, no isolated colonies are formed. Rather there may be several spots with fused black or green colonies, and whole plate sheet becomes yellow.

Putative *Salmonella* positive or negative colonies have to be confirmed by microbiological identification and confirmation tests.

Salmonella Negative

There is no colour with the medium. If it occurs, the sheet colour would change to red or reddish purple due to fermentation of lactose or sucrose by coliforms. No black or green colonies are observed.

Please note: The sheet colour might change to yellow caused by large colony numbers of *Pseudomonas* or *Proteus*. But yellow portion is small and limited because of their less motility. Other Gram-negative bacteria are inhibited in growth or don't show any colour shift.

F. General comments

- Some colonies may be only weakly coloured. Since false-negative results have to be avoided, those colonies have to be counted as positives.
- Make sure CompactDry plates are well closed during incubation, in order to avoid drying.
- High concentrations on plates will cause the entire growth area to become red/pink. In this case dilute the sample.
- After use please follow the current disposal regulations.
- The growth area is 20 cm². The back of the plate has a grid carved of 1 x 1 cm to make the colony counting easier. In case of any difficulties to count colonies due to large number of colonies grown, total viable count can be obtained by multiplying 20 by an average number of colonies per grid counted from several grids.
- Compact Dry plates may grow up to 300 colonies per plate. In those cases we recommend appropriate dilution of the samples. We recommend using the CompactDry Dilution Set (Art. No. 4823).
- Compact Dry plates are produced at an ISO 9001 certified site.

CompactDry SL 4817.1	40 plates
4817.2	240 plates