## CompactDry EC

4826

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AOAC approval No. 110402

MicroVal approval No. 0805-004LR according to ISO 4832 (2006) and No.0805-005LR according to ISO 16649-2-2:2001

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4826

ready-to-use, sterile

Ready-made plates with anhydrous chromogenic nutrient medium for detection of coliforms in general, particularly *E. coli*

### A. Introduction

Compact Dry plates made by NISSUI are suitable for detection and enumeration of microorganisms from foods and feeds, cosmetics and pharmaceutical samples, as well as for hygiene monitoring and contamination control in food-producing areas. Due to the chromogenic additives, bacterial colonies depict specific colouring, making them easy to interpret.

The nutrient medium used for Compact Dry EC plates is related to the Coliform chromogenic Agar. Colonies of *E. coli* show a bright blue colour, while other coliforms grow in red to pink coloured colonies. Addition of red and blue colonies results in total count of coliforms.

### B. Storage and Usage

CompactDry plates are stable at room temperature and at 4 °C. Store protected from moisture and light.

Do not use the product after its expiry date. Quality of the product is not warranted after being expired. Do not use the product that contains any foreign materials, discoloured or dehydrated, or its container is damaged. After opening the aluminum bag, any plates unused should be put back into the aluminum bag to be sealed with tape to avoid light and moisture, and use up as soon as possible.

### C. Sample preparation

#### From waters or liquid foods

One ml of specimen (dilute if necessary) is used for inoculation.

#### From solid foods

Add buffer solution to the sample and homogenize. One ml of this homogenate (dilute if necessary) is used for inoculation.

#### From surfaces

We recommend using the CompactDry Swab (Art. No. 4809). Use the swab to wipe the surface, put into the device with wiping solution. One ml of this specimen (dilute if necessary) is used for inoculation.

### D. Application

1. Open the cap and drop 1 ml of specimen on the middle of the Compact Dry plate. Specimen diffuses automatically and evenly into the sheet and transforms the dried sheet into a gel within seconds.
2. Put the cap back onto the plate and note date and sample on the memorandum section.
3. Turn over the capped plate and incubate at appropriate temperature.
4. After incubation count the number of coloured colonies underneath the plate. White paper placed under the plate helps with counting.

#### Incubation:

Please use the incubation time/temperature according to the national food analysis recommended for total viable count.

General recommendation: 37 °C for 24 hours

### E. Interpretation

*E. coli* forms blue colonies, other coliforms grow in red. Added numbers of red and blue colonies are the total coliform group count.

### F. General comments

- *E. coli* 0157 forms pink/red purple colonies.
- High concentrations on plates will cause the entire growth area to become white/pink. In this case dilute the sample.

- Some colonies may be only weakly coloured. Since false-negative results have to be avoided, those colonies have to be counted as positives.
- Make sure CompactDry plates are well closed during incubation, in order to avoid drying.
- High concentrations on plates will cause the entire growth area to become red/pink. In this case dilute the sample.
- After use please follow the current disposal regulations.
- The growth area is 20 cm<sup>2</sup>. The back of the plate has a grid carved of 1 x 1 cm to make the colony counting easier. In case of any difficulties to count colonies due to large number of colonies grown, total viable count can be obtained by multiplying 20 by an average number of colonies per grid counted from several grids.
- Compact Dry plates may grow up to 300 colonies per plate. In those cases we recommend appropriate dilution of the samples. We recommend using the CompactDry Dilution Set (Art. No. 4823).
- Compact Dry plates are produced at an ISO 9001 certified site.

**AOAC approval No. 110402**

**MicroVal approval No. 0805-004LR according to ISO 4832 (2006) and No.0805-005LR according to ISO 16649-2-2:2001**

**NordVal approval No. 036 and 037**

## CompactDry EC

**4826.1 40 plates**

**4826.3 120 plates**

**4826.2 240 plates**

### Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe

P.O. Box 100121 • 76231 Karlsruhe

Phone: +49 (0) 721/ 5606-0

Fax: +49 (0) 721/ 5606-149

info@carlroth.com • www.carlroth.com

ip 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.