MAXIMAL-WIEDERBELEBUNGSLÖSUNG

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For homogenization and dilution of microbiological samples

Approximate formula in g/l:

Casein peptone 1.0
Sodium Chloride..... 8.5
pH 7.0 ± 0.2

PREPARATION

Dissolve 9.5 g of the medium in one liter of deionized or distilled water. Mix well. Heat to dissolve the medium completely. Dispense and sterilize at 121 °C for 15 minutes.

In order to determine carbohydrate fermentation patterns, add 1.8 ml of 1 % phenol red to reconstitute the dry medium. After dispensing into test tubes with Durham gas collecting vials for gas detection, sterilize at 121°C for 15 minutes. Aseptically add sterile carbohydrate solution (Dextrose) to yield 1 % final concentration. Distribute carbohydrate in the tube by rotating gently.

USES

Maximum Recovery Solution is an isotonic diluent used for maximum recovery of microorganisms, and for the growth of bacterial cultures, mainly marine bacteria. ISO 6887 recommends this medium as a diluent for the preparation of initial suspension for microbiological samples. The low concentration of peptone does not cause a multiplication of the organisms within 1-2 hours of dilution of the sample. It is also used for carbohydrate fermentation tests in many food and environment studies, amongst others.

Inoculate tubes with a sample and incubate at 3 ±2°C for 18-24 hours. The fermentation of carbohydrate produces acid, causing a drop in the pH and a change of colour to yellow; gas production is indicated by gas bubbles in the Durham tubes.

MICROBIOLOGICAL TEST

The following results were obtained in the performance of the given medium from type cultures after incubation at a temperature of 35±2 °C and observed after 18-24 hours.

Microorganisms	Growth
<i>Escherichia coli</i> ATCC 25922	Good
<i>Salmonella typhimurium</i> ATCC 14028	Good
<i>Staphylococcus aureus</i> ATCC 25923	Good

Coccolin L, Manzano M, Cantur C., Comi G. (2001) *App. Environ Microbiolog.* 67:5113-21

Destoumieux-Garzon D. *et al.* (2001) *J. Biol Chem.* 276. 50:47070-7

ISO 6887 Microbiology – general guidance for the preparation of dilutions for microbiology examinations.

MAXIMUM RECOVERY SOLUTION

500 g

5766.1

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe • P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0 • Fax: +49 (0) 721/ 5606-149 • info@carlroth.com • www.carlroth.com

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

sse 07/2021