ROTI® DipSlides D/Eneutral

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Product Data Sheet

ROTI® DipSlides D/Eneutral

DipSlides and contact agar in one versatile tool. Sterile nutrient medium for detection of microorganisms from a sanitised environment or from liquids in presence of preservatives. Contains 5 neutralisers.

8023

- ✓ Ready-to-use and sterile
- ✓ Very easy application, reliable assay and easy interpretation
- ✓ Suitable for analysis of each surface and a multitude of liquids like water, milk, water-soluble paints, cosmetics, urine, etc.
- ✓ For microbiological laboratories, clinics and for hygiene control in food production

COMPOSITION

Both sides (D/Eneutral): D/E Neutralizing Agar (Art. No. 8850) for detection and enumeration of microorganisms from a sanitised environment (**violet-red**)

STORAGE, EXPIRY AND SHIPMENT

ROTI®DipSlides are shipped in special boxes on gel ice. Please make sure that ROTI®DipSlides are stored immediately at 4 -15 °C. Store boxes lid-up, in order to avoid moistening of the agar with condensed water.

ROTI®DipSlides may not be frozen. Slides that have been frozen have to be disposed of, since the medium is like to have taken severe damage and growth of microorganisms is significantly hindered.

Store at 4-15 °C and protected from direct light. In order to avoid drying of the agar, temperature should vary for maximally ± 2 °C. Expiry of ROTI®DipSlides is 5 months after production. Expiry date is noted on each package and on each vial with „EXP“.

COMPOSITION OF AGARS (in g/l, approximated)

For further information, please see the respective data sheets of the agars.

D/Eneutral (8850)	Dextrose 10, Lecithin 7, Natriumthiosulfat 6, Caseinpepton 5, Polysorbat 80 5, Hefeextrakt 2,5; Natriumbisulfid 2,5; Natriumthioglycolat 1, Bromkresolpurpur 0,02, Agar 15, pH-Wert 7,6 $\pm 0,2$ Colour: violet-red
-------------------	--

APPLICATION

1. Remove carrier from the sterile vial. The screw cap serves automatically as holder.
2. Press agar gently onto the surface that has to be tested. Alternatively, immerse the agar for several seconds in the liquid, in order to wet both sides completely. Carefully allow excess sample to run off. In case of samples of high viscosity, we recommend diluting the sample by 1:2 to 1:10. Dilution must be taken into account when calculating the sample contamination. Please note: In case the DipSlides are used as contact agar, two different surfaces may be tested. We strongly recommend, however, NOT to use one DipSlide for testing two different liquids, since the risk of cross-contamination is very high! Do NOT touch the agar (risk of contamination!)
3. Replace carrier into the vial and close the cap moderately.
4. Incubate vial with carrier at a temperature suitable for the bacteria (e.g. 35 ± 2 °C for 18 - 24 h). During incubation, the vials should be screwed on only loosely, in order to avoid condensation of water. *Please note:* In accordance to official regulations or internal protocols other incubation parameters may be chosen.
5. After the appropriate incubation time span count and register colonies.

EVALUATION OF RESULTS

After the incubation time span required colonies are counted and noted as cfu (colony forming units). For this purpose, the carrier should not be removed from the vial, in order to prevent release of the microorganisms. After incubation, tightly close the tube and let it cool for at least 2 hours at 4 °C. Condensed water may then be tapped into the lid, enabling detailed analysis of the agars through the plastic wall.

Any validation required has to be undertaken by the user. We recommend to test each sample liquid or surface at least twice (with two ROTI®DipSlides), in order to validate results. While calculating sample contamination, eventual sample dilution factors have to be taken into account. Evaluation of results is done semi-quantitatively according to the table given below.

For solutions the following approximation may be applied:

Cfu of one side x 40 is approx. cfu / ml solution. Detection limit is 10² CFU / ml solution.

Table 1: Evaluation

Colony no.	cfu / ml *	Result	Interpretation
0	<10 ²	Negative	Contamination below lower detection limits
1-5	ca. 10 ²	Weakly positive	Weak contamination of solution/surface
10-50	ca. 10 ³	Positive	Low contamination of solution/surface
100-500	ca. 10 ⁴	Strongly positive	High contamination of solution/surface
>500	>10 ⁴	Highly positive	Very high contamination of solution/surface

* This calculation may be used as approximation to numbers only. For definite quantitation of contaminated solutions or surfaces, sophisticated microbiological procedures have to be performed using Roth nutrient dry media.

DISPOSAL

Disposal of ROTI®DipSlides bearing cultures has to be done according to all common regulations for disposal of the respective microorganisms, and to regulations of good laboratory practice. E.g. autoclaving for 20 mins. at 121 °C or incubation in 70 % ethanol.

MICROBIOLOGICAL TEST

D/Eneutral (35 ± 2 °C / 18 - 24 h)

Test Strain	Growth
<i>Escherichia coli</i> ATCC 25922	gut
<i>Bacillus subtilis</i> ATCC 6633	gut
<i>Pseudomonas aeruginosa</i> ATCC 27853	gut
<i>Staphylococcus aureus</i> ATCC 25923	gut
<i>Salmonella typhimurium</i> ATCC 14028	gut

PLEASE NOTE!

Application and evaluation of results of ROTI®DipSlides may only be carried out by trained personnel. The current instructions-for-use gives only general information for test performances, which have to be adopted by each user for his or her special requirements, particularly if specific regulations have to be followed. Any validation required has to be undertaken by the user.

Table 2: Overview

Ord. No.	Product	Agar Side 1	Agar Side 2	Application	Inc. temp.	Inc. time
8023.1	ROTI®DipSlide D/Eneutral	D/E Neutralizing Agar (8850)	D/E Neutralizing Agar (8850)	whole germs after disinfection	35 ±2 °C	18 - 24 h

Our whole product range of ROTI®DipSlides for many different applications may be found in our current catalogue or in our webshop.

ROTI®DipSlide D/Eneutral

20 slides

8023.1

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe • P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0 • Fax: +49 (0) 721/ 5606-149 • info@carlroth.com • www.carlroth.com

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

sse 07/2021