COLIFORME CHROMOGENER AGAR (ISO)

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Product Data Sheet

COLIFORMS CHROMOGENIC AGAR (ISO)


For detection and enumeration of *E. coli* and other coliforms in water

Acc. to ISO 9308 / ISO 11133 / for Microbiology

8849

Approximate formulation in g/l:

Sodium chloride (NaCl)	5.0
di-sodium hydrogen phosphate	2.7
Sodium dihydrogen phosphate.....	2.2
Yeast extract.....	2.0
Casein, enzymatically digested	1.0
Sodium pyruvate.....	1.0
Tryptophane	1.0
Sorbitol.....	1.0
Chromogene mixture.....	0.4
Sodium heptadecyl sulphate (Tergitol-7).....	0.15
Agar	10
pH value.....	6.8 ± 0.2

A: *E. coli* ATCC 25922 (blue-dark violet)

B: *Enterobacter aerogenes* ATCC 13048 (red-rose)

PREPARATION

Suspend 26.45 g of the medium in one litre of distilled water.

Mix well and heat under frequent agitation. Boil for 1 minute until complete dissolution.

Do not overheat! DO NOT AUTOCLAVE!

Let the medium cool to 45 - 50 °C and pour in petri dishes. The prepared medium is amber and should be stored protected from light at 8 - 15 °C.

USES

The Coliforms chromogenic Agar is recommended by ISO 9308-1 (CCA) as a selective medium for detection of *E. coli* and other coliforms from water and food samples. The qualitative ingredients enable a fast growth of the colonies, inhibiting the Gram positive bacteria by the contained Tergitol 7. The medium also contains tryptophans, enabling the performance of an indol test in order to determine *E. coli*. Coliforms are able to metabolise the chromogenic component salmon-gal by the formation of the enzyme β-D-galactosidase, and therefore grow in red colonies. Since *E. coli* forms additionally the enzyme β-D-glucuronidase, it metabolises the X-glucuronide and the colonies typically grow in dark blue-violet.

The TWVO recommends to filter the sample first using a membrane, and place this membrane on top of the CCA plate. Incubate at 36 ± 2°C for 21 ± 3h.

For evaluating the results, the red colonies are counted first and noted (Coliforms, not *E. coli*). Colourless colonies are classified by an oxidase test as coliform (oxidase negative) or not coliform (oxidase positive).

Then the dark blue-violet colonies are counted and noted as the count of *E. coli*. The whole germ count of coliform germs is the result of all oxidase negative colonies plus all red colonies and plus all dark blue-violet colonies.

Please note: Some strains of *Shigella* are able to form the glucuronidase, too, and hence grow in light blue colonies.

The strains of *E. coli* that do not form glucuronidase, grow in colourless colonies, e.g. *E. coli* O157:H7.

MICROBIOLOGICAL TEST

The following results were obtained in the performance of the given medium from type cultures after incubation at a temperature of 36 ± 2 °C and observed after 21 ± 3 hours.

Microorganisms	Growth	Colony colour
<i>Escherichia coli</i> ATCC 25922	Good	Blue-dark violet
<i>Escherichia coli</i> ATCC 8739	Good	Blue-dark violet
<i>Enterobacter aerogenes</i> ATCC 13048	Good	Red-pink
<i>Pseudomonas aeruginosa</i> ATCC 10145	Good	Colourless
<i>Enterococcus faecalis</i> ATCC 19433	Inhibited	---

Acc. to:

- 1.) DIN EN ISO 9308-1:2014-12 Wasserbeschaffenheit
- 2.) DIN EN ISO 7218:2014-09
- 3.) Byamukama et al (2000). *J. Appl Environ Microbiol.* 66(2): 864-868
- 4.) Geissler et al (2000). *J. Appl- Microbiol.* 88 : 280 - 285
- 5.) Ossmer et al (1999) Poster Presentation Congreso de la Sociedad, Espanola de Microbiologia, Granada, Spain.
- 6.) USEPA: 40 CFR Part 141: Federal Register/Vol. 67, No. 209 (2002): Rules and Regulations
- 7.) Lange, Strathmann und Oßmer (2013). *J. Appl. Microbiol.* 57(6): 547-553
- 8.) ICFMH-WPCM (2012) Reference Strain Catalogue

COLIFORMS CHROMOGENIC AGAR (ISO)

500 g

8849.1

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe • P.O. Box 100121 • 76231 Karlsruhe

Phone: +49 (0) 721/ 5606-0 • Fax: +49 (0) 721/ 5606-149 • info@carlroth.com • www.carlroth.com

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

sse 06/2021