## LB-Agars mit Antibiotika

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jh 09/2022

## LB Agars with Antibiotics

Agars for selective growth of recombinant *Escherichia coli* in genetic engineering. Premixed agar formulation including antibiotics.

8852, 8853, 8855, 8858, 8860

### Approximate formulation in g/l:

Art. No.	Nutrient Medium	Tryp- tone*	Yeast extract	NaCl	Antibiotic	Agar- Agar	pH value
8852	LB-Agar (Lennox) / Amp100	10 g/l	5 g/l	5 g/l	Ampicillin 100 mg/l	15 g/l	7,0±0,2
8853	LB-Agar (Lennox) / Kana50	10 g/l	5 g/l	5 g/l	Kanamycin 50 mg/l	15 g/l	7,0±0,2
8855	LB- Agar (Luria/Miller) / Amp100	10 g/l	5 g/l	10 g/l	Ampicillin 100 mg/l	15 g/l	7,0±0,2
8858	LB- Agar (Luria/Miller) / Amp50	10 g/l	5 g/l	10 g/l	Ampicillin 50 mg/l	15 g/l	7,0±0,2
8860	LB-Agar (Luria/Miller) / Kana50	10 g/l	5 g/l	10 g/l	Kanamycin 50 mg/l	15 g/l	7,0±0,2

\* Pancreatically digested peptone from casein.

### PREPARATION

Suspend 35 g (Lennox) or 40 g (Luria/Miller), respectively, in one litre of distilled or deionized water. Mix well and heat with frequent agitation. Boil for one minute until fully dissolved.

LB agars with **ampicillin**: **Do not overheat. Do not autoclave.**

LB agars with **kanamycin** or **chloroamphenicol**: **Do not overheat. Do not autoclave.**

Cool to 45-50 °C in a water bath. Mix well and dispense into sterile petri dishes.

The hardened dishes should be stored bottom up at 8-15 °C and may be used for approx. 2 weeks, as long as sterility is guaranteed. The colour of the agar is amber, slightly opalescent.

### USES

Use of LB Agars with antibiotics is recommended for the selective propagation of *E. coli* recombinant with plasmids that carry a selection marker in terms of resistance against ampicillin, kanamycin or chloroamphenicol. The transformed *E. coli* are plated **directly** onto the prepared selective agar. Further addition of antibiotics during preparation of the agar dish is not necessary. IPTG and X-β-Gal may be added for induction of protein expression (blue/white selection).

LB agars according to the Lennox modification (less salt) are recommended for growth of *E. coli* strains that require less salt concentration.

LB (Lysogeny **B**roth) media is a group of nutritionally rich growth solutions for the propagation and maintenance of *E. coli* bacteria. Originally, the first LB medium has been developed by Guiseppe Bertani in 1951 in order to optimise growth of *Shigella* and plaque formation of P1, P2 and P3 phages<sup>1</sup>. Since then, the media were, and still are, widely used for growth of *E. coli* in general and recombinant *E. coli* in particular, as well as for propagation of plasmids and expression proteins *in vitro*. During the following years, some modifications were implemented including the reduction of sodium chloride<sup>2,3</sup>, amongst other things in order to provide optimised conditions for salt sensitive bacterial strains.

Tryptone provides nitrogen, vitamins, minerals and amino acids essential for growth. Yeast extract is a rich source of vitamins, particularly the B-group. Sodium chloride supplies essential electrolytes for transport and osmotic balance. In LB agars, high-pure agar-agar is added as solidifying agent.

## MICROBIOLOGICAL TEST

The following results were obtained in the performance of the given medium from type cultures after incubation at a temperature of 35±2 °C and observed after 18-24 hours.

### General QC assays

Microorganisms	Growth
<i>Escherichia coli</i> ATCC 25922	Inhibited
<i>Escherichia coli</i> ATCC 8739	Inhibited

### Specific QC assays

Art. No.	Agar medium	Microorganisms	Growth
8852	LB-Agar (Lennox) / Amp100	<i>Escherichia coli</i> DH5 alpha + pUC19	Good
8853	LB-Agar (Lennox) / Kana50	<i>Escherichia coli</i> DH5 alpha + pHSG 298	Good
8855	LB- Agar (Luria/Miller) / Amp100	<i>Escherichia coli</i> DH5 alpha + pUC19	Good
8858	LB- Agar (Luria/Miller) / Amp50	<i>Escherichia coli</i> DH5 alpha + pUC19	Good
8860	LB-Agar (Luria/Miller) / Kana50	<i>Escherichia coli</i> DH5 alpha + pHSG 298	Good

Acc. to:

- 1.) Bertani, G. Studies on lysogenesis. I. The mode of phage liberation by lysogenic *Escherichia coli*. *J. Bacteriol.* 62:293-300 (1951)
- 2.) Lennox, E. S. Transduction of linked genetic characters of the host by bacteriophage P1. *Virology* 1:190-206 (1955)
- 3.) Luria, S. E.; Adams, J. N.; Ting, R. C. Transduction of lactose-utilizing ability among strain of *E. coli* and *S. dysenteriae* and the properties of the transducing phage particles. *Virology* 12:348-390 (1960)
- 4.) Miller, J. H. Experiments in molecular genetics. *Cold Spring Harbor Laboratory*, Cold Spring Harbor, N.Y. (1972)
- 5.) Atlas, R.M., Parks L.C. Handbook of Microbiological Media. *CRC Press*, Inc. London (1993)
- 6.) Sambrook and Russell. Molecular cloning: a laboratory manual, 3rd ed. *Cold Spring Harbor Laboratory Press*, Cold Spring Harbor, N.Y. (2001)

**LB Agar (Lennox) / Amp100**  
 500 g 8852.1  
 1 kg 8852.2  
 2.5 kg 8852.3

**LB Agar (Lennox) / Kana50**  
 500 g 8853.1  
 1 kg 8853.2

**LB Agar (Luria/Miller) / Amp100**  
 500 g 8855.1  
 1 kg 8855.2  
 2.5 kg 8855.3

**LB Agar (Luria/Miller) / Amp50**  
 500 g 8858.1  
 1 kg 8858.2  
 2.5 kg 8858.3

**LB Agar (Luria/Miller) / Kana50**  
 500 g 8860.1  
 1 kg 8860.2

### Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe • P.O. Box 100121 • 76231 Karlsruhe  
 Phone: +49 (0) 721/ 5606-0 • Fax: +49 (0) 721/ 5606-149 • info@carlroth.com • www.carlroth.com

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

jh 09/2022