## LACTOSE-TTC-AGAR MIT NATRIUM-HEPTADECYLSULFAT (BASIS)

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# Product Data Sheet


## TTC-CHAPMAN-AGAR (WITH SODIUM HEPTADECYLSULPHATE) (BASE)

For isolation and enumeration of coliform bacteria from water samples by the membrane filtration method  
CP65 DEV / for Microbiology

### Formulation in g/l:

Lactose .....	20.0
Peptone (ex meat) .....	10.0
Yeast extract .....	6.0
Beef extract .....	5.0
Bromothymol blue .....	0.05
Sodium heptadecyl sulphate (Tergitol 7) .....	0.1
Agar .....	15.0
Final pH .....	7.2±0.1

Figure: A) *E. coli* ATCC 25922 B) *Klebsiella* spp. C) non-lactose fermenting species (*Salmonella*) D) uninoculated TTC Agar

### PREPARATION (EN ISO 9308-1:2000 – B.1)

Suspend 56.2 g of the medium in 1 l distilled water. Solubilise under heating with frequent agitation and boil for few seconds. Pour aliquots of 200 ml each and sterilise for 15 minutes at 121 °C in an autoclave. Let cool to approx. 50 °C and add 0.5 ml each of a sterile, 1 % TTC solution (CP72.1, reconstitute one vial (50 mg) in 5 ml sterile, distilled water). Mix carefully and pour into Petri dishes (at approx. 5 mm thickness). Let cool. The Petri dishes have to be stored in the dark at 3-5 °C. Do not use longer than 7 days. Longer storage result in decreasing colour reaction in the border area.

### USES

TTC-Chapman-Agar with Sodium Heptadecyl sulphate is prescribed by European ISO Standard 9308-1:2000 for analyzing water samples to determine the presence of lactose-fermenting bacteria on membrane filter plates. This detection is the first step towards identifying coliform bacteria. The water samples are filtered through sterile membrane filters (cellulose ester, 0.45 µm, Art. No. CY49.1) and the filters incubated on TTC-Chapman-Agar with Sodium Heptadecyl sulphate. Sodium heptadecyl sulphate (Tergitol 7) inhibits Gram positive bacteria and the spreading of *Proteus*. Once incubation is completed, the colonies are analyzed for pigment formation.

Caused by dehydrogenases of the respiratory chain of living cells, the added TTC is reduced to formazane, causing a basic dark-red colouring of all colonies. Lactose-fermenting bacteria produce acids during fermentation which cause a drastic pH change in the medium, indicated by the pH-indicator bromothymol blue as a yellow staining of the entire medium or by a yellow corona around the colonies. Hence, coliform bacterial colonies (e.g. *Escherichia*, *Enterobacter*, *Klebsiella*, *Citrobacter*) are stained yellow or orange on yellow medium. Non-lactose fermenting (β-galactosidase-negative) bacteria (e.g. *Pseudomonas*, *Salmonella*, *Shigella*, *Proteus*) are stained purple to violet on blue medium. All yellow or orange-red colonies with a yellow corona have to be subcultivated on tryptone soy agar (Art. No. CP70.1) and tryptophan culture broth (Art. No. CP71.1) for species confirmation. Coliform bacteria are Gram negative, bacillary bacteria which do not represent a taxonomic group, but are characterized for being able to produce energy from lactose under acid and aerosis. They are prevalent particularly in the intestines of warm-blooded animals. They can, however, survive for a time outside their natural environment and are therefore extremely suitable for analyzing water samples for faecal contamination. However, some types of coliform bacteria do in fact live in the soil and in surface waters with the result that the respective contamination of the water samples may be ascribed to other errors during treatment of water. European ISO Standard 9308-1:2000 defines reference procedures for the detection and enumeration of coliform bacteria in all types of water meant for human use. Two procedures are given – a standard test which is especially suitable for water low in germs and a quick-test for short term but less sensitive analysis.

### MICROBIOLOGICAL TEST

The following results were obtained in the performance of the given medium from type cultures after incubation at a temperature of 35 ± 2°C and 44°C and observed after 24-48 hours.

Microorganisms	44 ° C	35 ° C	Colony Colour	Medium Colour
<i>Escherichia coli</i> ATCC 25922	Good	Good	Yellow with orange center	Yellow
<i>Citrobacter</i> spp.	Nil	Good	Yellow with orange Center	Yellow
<i>Klebsiella</i> spp.	Nil	Good	Red to yellow	Yellow
<i>Enterobacter</i> ATCC 13048	Nil	Good	Red, dark yellow or green-yellow with orange center	Yellow
Non lactose-fermenting species	Nil	Good	Light violet to purple	Blue

**TTC CHAPMAN AGAR (WITH SODIUM HEPTADECYL SULPHATE) (BASE)**

**500 g**

**CP65.1**

### Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe • P.O. Box 100121 • 76231 Karlsruhe  
Phone: +49 (0) 721/ 5606-0 • Fax: +49 (0) 721/ 5606-149 • info@carlroth.com • www.carlroth.com

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

sse 06/2021