Produkt-Datenblatt

TRYPTON-GALLE-AGAR (TBA)

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Product Data Sheet

TRYPTONE BILE SALTS AGAR (TBA)

For detection and enumeration of *E. coli* in foods and water samples.

Agar for quick-testing of indol production

Acc. EN ISO 9308-1:2000

CP69

Formulation in g/l:

Tryptone.....	20.0
Bile salts.....	1.5
Agar	15.0
Final pH.....	7.2±0.1

PREPARATION (EN ISO 9308-1:2000 – B.4)

Suspend 36.5 g of the medium in 1 l distilled water. Solubilise under heating with frequent agitation. Let cook for 1 minute. Pour aliquots of 200 ml each and sterilise for 15 minutes at 121 °C in an autoclave. Let cool to approx. 50 °C and pour into Petri dishes (at approx. 5 mm thickness). Let cool. Double layer agar plates for speedy testing of *E. coli* are produced by pouring a thin layer (approx. 1 mm) of warm TSA-agar (Art. No. CP70.1) over ready-to-use TBA-plates. Prepare freshly for each test.

USES

Tryptone Bile Salts Agar is prescribed by European ISO Standard 9308-1:2000 for analyzing water samples to determine the presence of *Escherichia coli* with a quick-test. After filtering the water samples through sterile membrane filters (cellulose ester, 0.45 µm, Art. No. CY49.1), preincubate the membrane on a Tryptone Soy Agar (Art. No. CP70.1) plate and transfer it to a Tryptone Bile Salts Agar plate after 4-5 hours. Transfer the membrane to a filter paper saturated with indole reagent after incubation. The colonies which appear red under UV-light should be counted as *Escherichia*. *Escherichia coli* (originally *Bacterium coli* commune) are Gram negative, rod-shaped, facultative anaerobic bacteria belonging to the enterobacteriaceae family. The coliform bacteria group was named after them. They can be found in the physiological intestinal flora of warm-blooded animals. Because they are facultative pathogenic bacteria, they cause infection in other organs such as the urinary passage or lungs. Some strains are even intestinally pathogenous. *E. coli* bacteria are typical faecal contamination in drinking water and are regarded as the main organism for these tests. *Escherichia* transform lactose under acid production via β-galactosidase (yellow colonies with a yellow surrounding area on TTC-Chapman-Agar (Art. No. CP65)). Furthermore, they are negative in the oxidase test and positive in the quick-test in indole production with dimethyl-amino-benzyldehyde.

European ISO Standard 9308-1:2000 defines reference procedures for the detection and enumeration of coliform bacteria in all types of water meant for human use. Two procedures are given – a standard test which is especially suitable for water low in germs and a quick-test for short term but less sensitive analysis.

MICROBIOLOGICAL TEST

The following results were obtained from type cultures after incubation at a temperature of 44 ± 0.5 °C and observed after 19-20 hours.

Microorganisms	Growth	Colony Colour	Indole
<i>Escherichia coli</i> ATCC 25922	Good	Red	+
<i>Klebsiella pneumoniae</i> ATCC 13833	Nil	-	-

TRYPTONE BILE SALTS AGAR (TBA)

500 g

CP69.1

Carl Roth GmbH + Co. KG

Schoenperlenstraße 3-5 • 76185 Karlsruhe • P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0 • Fax: +49 (0) 721/ 5606-149 • info@carloth.com • www.carloth.com

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

sse 07/2021