Produkt-Datenblatt

TRYPTON-SOJA-AGAR (TSA)

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Product Data Sheet

TRYPTONE SOY AGAR (TSA)

For detection and enumeration of *E. coli* and other coliform bacteria in foods and water samples. Non-selective agar for subculturing prior to oxidase tests

Acc. EN ISO 9308-1:2000

CP70

Formulation in g/l:

Casein, tryptic digest	15.0
Soy peptone.....	5.0
Sodium chloride	5.0
Agar	15.0
Final pH.....	7.2±0.1

PREPARATION (EN ISO 9308-1:2000 – B.3)

Suspend 40 g of the medium in 1 l distilled water. Solubilise under heating with frequent agitation.

Pour aliquots of 200 ml each and sterilise for 15 minutes at 121 °C in an autoclave.

Let cool to approx. 50 °C and pour into Petri dishes (at approx. 5 mm thickness). Let cool.

Double layer agar plates for speedy testing of *E. coli* are produced by pouring a thin layer (approx. 1 mm) of TSA-agar over ready-to-use TBA-plates (Art. No. CP69.1). Prepare freshly for each test.

USES

Tryptone Soy Agar is prescribed by European ISO Standard 9308-1:2000 for analyzing water samples to determine the presence of coliform bacteria. After filtering the membranes and culturing the samples on TTC-Chapman-Agar with Sodium Heptadecylsulphate (Art. No CP65.1), subculture all yellow to orange-coloured colonies with a yellow corona on Tryptone Soy Agar. Analyze the colonies with an oxidase test on filter paper after incubation. Oxidase-negative colonies should be counted as coliform bacteria. For a quick *E. coli* test, place the filtered membrane directly on the TSA-agar plates or on a double layer TSA/TBA combination plate. Coliform bacteria are Gram negative, rod-shaped bacteria which do not represent a taxonomic group, but are characterized for being able to produce energy from lactose under acid and aeration. They are prevalent particularly in the intestines of warm-blooded animals. They can, however, survive for a time outside their natural environment and are therefore extremely suitable for analyzing water samples for faecal contamination. However, some types of coliform bacteria do in fact live in the soil and in surface waters with the result that the respective contamination of the water samples may be ascribed to other errors during treatment of water.

European ISO Standard 9308-1:2000 defines reference procedures for the detection and enumeration of coliform bacteria in all types of water meant for human use. Two procedures are given – a standard test which is especially suitable for water low in germs and a quick-test for short term but less sensitive analysis.

MICROBIOLOGICAL TEST

The following results were obtained in the performance of the given medium from type cultures after incubation at a temperature of 36 ± 2 °C and observed after 18-24 hours.

Microorganisms	Growth
<i>Escherichia coli</i> ATCC 25922	Good
<i>Klebsiella pneumoniae</i> ATCC 13833	Good

TRYPTONE SOY AGAR (TSA)

500 g

CP70.1

1 kg

CP70.2

Carl Roth GmbH + Co. KG

Schoenperlenstraße 3-5 • 76185 Karlsruhe • P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0 • Fax: +49 (0) 721/ 5606-149 • info@carlroth.com • www.carlroth.com

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

sse 07/2021