## ORANGENSERUMAGAR

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# Product Data Sheet

## ORANGE SERUM AGAR

For isolation, detection and enumeration of acid-tolerant pathogens in fruit juices and fruit juice concentrates, particularly citrus products

HP09

### Approximate formula in g/l:

Casein peptone .....	10.0
Orange extract.....	5.0
Glucose.....	4.0
Potassium phosphate.....	3.0
Yeast extract.....	3.0
Agar .....	15.0
Final pH .....	5.5 ± 0.2

### PREPARATION

Suspend 40 g of the medium in one liter of distilled or deionised water. Heat with frequent agitation and boil for one minute. Dispense into appropriate containers and sterilise at 118 °C for 15 minutes in an autoclave.

**Do not overheat!**

### USES

Orange Serum Agar is recommended for proliferation of the acid-tolerant, lactic acid producing micro flora spoiling citric products, such as *Bacillus*, *Lactobacillus*, molds, etc, based on the medium suggested in 1951 <sup>(1)</sup>. By adding orange extract to the formulation, the medium has been optimized for production control in the fruit juice industry. The low pH value of the medium, in addition to the acid formation resulting from glucose fermentation, promotes growth of *Lactobacillus*, *Leuconostoc*, yeast and molds.

The low pH of citric juices and products limits the growth of microorganisms to acid tolerant pathogens.

In general, use the plate-pouring method. 0.1 - 1 ml of the appropriate dilution of the test sample is applied into a petri dish. Subsequently, approx. 20 ml of the cooled agar (45 - 55 °C) is added. Agar and samples are mixed by gently swirling and the plate is allowed to cool and solidify in a level position. The petri dishes are aerobically incubated at 35 ± 2 °C (*Lactobacillus*) or 30 ± 2 °C (other microorganisms) and daily observed for 40 - 48 hours. In case presence of fungi is suspected, prolong incubation to 5 days.

Since Orange Serum Agar is not a selective medium, we recommend species identification on specialised agar media like MRS-Agar (Art. No. X924) for *Lactobacillus*, and Yeast-Glucose Agar with chloramphenicol (Art. No. CL00) or Sabouraud 2 % Glucose Agar (Art. No. AE22) for yeasts and molds.

### MICROBIOLOGICAL TEST

The following results were obtained in the performance of the given medium from type cultures after incubation at a temperature of 35 ± 2 °C and 30 ± 2 °C and observed after 40 - 48 hours.

Microorganisms	Growth at 30 °C	Growth at 35 °C
<i>Aspergillus brasiliensis</i> ATCC 16404*	Good	<i>n.d.</i>
<i>Lactobacillus fermentum</i> ATCC 9338	<i>n.d.</i>	Good
<i>Saccharomyces cerevisiae</i> ATCC 9763	Good	<i>n.d.</i>

\* Renamed by the ATCC from *A. niger*, January 2011

Acc. to: 1.) Hays, G.L. (1951) *Proc. Fla. State Hort. Soc.* 54:135-7.

**Warning** H315-H317-H411-EUH208 P280-P302+P352-P332+P313

**ORANGE SERUM AGAR**

**500 g**

**HP09.1**

### Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe • P.O. Box 100121 • 76231 Karlsruhe  
Phone: +49 (0) 721/ 5606-0 • Fax: +49 (0) 721/ 5606-149 • info@carlroth.com • www.carlroth.com

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

sse 07/2021