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## ROTI®Store Cryoröhrchen

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# Instructions for use

## ROTI®Store cryo-vials

For preparing deep-freeze cultures of microorganisms

### A. Product description

ROTI®Store is a product for preparing deep-freeze cultures. It enables safe storage and recultivation of microorganisms. These bind to glass beads of exceptional surface quality and are protected by a special cryo-medium when freezing. They may be stored at -20 °C for maximal one week. For long term-storage we recommend a minimum of -80 °C. In order to recultivate simply remove a glass bead and bring it into contact with desired culture medium.

Also suitable for conservation of marine microorganisms, even during field work (Siebert et al., *Appl. Environ. Microbiol.* 2004 (70) 5912-5).

#### **Carl Roth GmbH + Co. KG**

Schoemperlenstraße 3-5  
76185 Karlsruhe  
Postfach 100121  
76231 Karlsruhe  
Telefon: +49 (0) 721/ 5606-0  
Telefax: +49 (0) 721/ 5606-149  
E-Mail: info@carlroth.de  
Internet: www.carlroth.de


jh 07/2020

### B. Pack contents and storage

ROTI®Store contains 50 cryo-vials with 25±5 glass beads per vial in a special cryo-medium. The non-inoculated ROTI®Store cryo-vials can be stored between +4 °C to +25 °C (optimal storage: +4 to +8 °C). They should not be constantly exposed to sunlight.

### C. Inoculation

- You will require a fresh overnight culture or colonies resuspended in culture medium with an OD<sub>660</sub> of at least 2 \*
- Check the vial prior to use for possible contamination (cloudiness). Please label the cryo-vial with a waterproof marker (e.g. ROTILABO®LabMarker, Roth Art. No. 3083.2) and add approx. 0.5 ml culture solution.
- Close the vial tightly and invert it carefully several times.
- Remove the supernatant above the pellet with a sterile pipette (Roth Art.No. EA69.1) and dispose of it correctly.
- The inoculated and tightly closed ROTI®Store cryo-vial can now be frozen and stored at -20 °C for one week or at -80 °C (or lower) for years.

\*Direct inoculation of ROTI®Store cryo-vials using picked colonies is possible, it does result,

however, in much less efficient recultivation rates. In these cases, we recommend to resuspend the colony well in the cryo-medium. Slightly agitate the tube 15-20 mins. at 4 °C and remove the supernatant only then.

### D. Recultivation

Remove the ROTI®Store cryo-vial from the deep-freezer and transport it in a cold cryo-block (e.g. ROTILABO®ice box, Roth Art.No. K088.1). The vial should not be allowed to thaw.

Remove a bead under sterile conditions with a sterile pincette or inoculation needle (e.g. Roth Art.No. EA88.1) and then close the cryo-vial.

The bead may now be added directly to the liquid medium or rolled out onto a solid culture medium.

Ø: 11 mm  
length: 47 mm

### ROTI®Store cryo-vials

<b>P730.1</b>	50 vials with 25 glass beads per vial
<b>P730.2</b>	250 vials with 25 glass beads per vial
<b>P730.3</b>	5 vials with 25 glass beads per vial