## KARTOFFEL-EXTRAKT-GLUCOSE AGAR

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sse 06/2021

## POTATO-EXTRACT-GLUCOSE AGAR


Recommended by the Harmonized Method (Ph. Eur. 6.0), Ph. Eur. 6.3, and 7.0  
For the cultivation and identification of yeasts and molds

Ph. Eur.  
X931

### Formulation in g/l:

Potato Infusion..... 4.0  
Dextrose ..... 20.0  
Agar ..... 15.0  
pH ..... 5.6 ± 0.2

*Aspergillus brasiliensis*  
ATCC 16404

### PREPARATION

Suspend 39 g of the medium in one liter of distilled or deionized water. Mix well and heat with frequent agitation. Boil for one minute. Sterilize at 118 °C for 15 minutes. Cool to 40-45 °C and pour into Petri dishes. The medium is amber and slightly opalescent due to the potato infusion.

### USES

Potato-Extract-Glucose Agar is used for the isolation, enumeration and identification of yeast and molds. This medium is recommended by the Association of Public Health (APHA).

Potato-Extract-Glucose Agar can be used in the analysis of dairy products, bottled drinks, frozen food, and other types of food. It can also be used in the identification of fungi and yeasts in parallel with their cellular morphology or in methods of microcultivation in slides.

When the medium is to be used for enumeration of molds and yeasts, add to the medium, sterilized and cooled to 45-50 °C, approximately 14 ml of a sterilized 10 % solution of tartaric acid to obtain a pH of 3.5. Do not heat the medium after adding the acid, because the agar may hydrolyze and not solidify.

### MICROBIOLOGICAL TEST

The following results were obtained in the performance of the given medium from type cultures after incubation at 35±2 °C and observed after 18-24 hours.

Microorganisms	Growth
<i>Aspergillus brasiliensis</i> ATCC 16404*	Good
<i>Candida albicans</i> ATCC 10231	Good
<i>Saccharomyces cerevisiae</i> ATCC 9763	Good
<i>Trychophyton mentagrophytes</i> ATCC 9533	Good

\*Renamed by the ATCC from *A. niger*, January 2011

<b>POTATO-EXTRACT-GLUCOSE AGAR</b>	500 g	X931.1
	1 kg	X931.2
	2,5 kg	X931.3

### Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe • P.O. Box 100121 • 76231 Karlsruhe  
Phone: +49 (0) 721/ 5606-0 • Fax: +49 (0) 721/ 5606-149 • info@carlroth.com • www.carlroth.com

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

sse 06/2021