
SALMONELLA-SHIGELLA AGAR

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Product Data Sheet


SALMONELLA-SHIGELLA AGAR

For the isolation of pathogenic enterobacteria
X934

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Beef Extract	5.0
Lactose	10.0
Sodium Citrate	8.5
Ferric Citrate	1.0
Neutral Red.....	0.025
Peptone.....	5.0
Bile Salts	8.5
Sodium Thiosulphate	8.5
Brilliant Green	0.00033
Agar.....	13.5
pH-value	7.0 ± 0.2

Salmonella typhimurium
ATCC 14028

Organism	Colony
<i>E. coli</i>	small, red-pink
<i>Enterobacter and Klebsiella</i>	larger than <i>E. coli</i> , mucoid, pale-cream to pink
<i>Proteus and some Salmonellae</i>	colourless, transparent with black center (in case of H ₂ S production)
<i>Shigella and several Salmonellae</i>	clear, colourless, transparent

PREPARATION

Suspend 60 g of the medium in one liter of deionised or distilled water. Mix well. Heat with frequent agitation, and boil for one minute. Do not sterilise in an autoclave. Pour into plates. Avoid freezing. The plates of the medium can be kept for at least a week in refrigeration.

USES

Salmonella-Shigella Agar is a selective and differential medium and is widely used in sanitary bacteriology to isolate Salmonella and Shigella from feces, urine, and fresh and canned foods. Inhibition of Gram positive microorganisms is obtained by the bile salts mixture.

Due to its strong inhibitory power, Salmonella-Shigella Agar can be streaked with a heavy inoculum but other less inhibitory media such as e.g. MacConkey Agar (Art. No. X922) or XLD Agar (Art. No. X941) should be streaked in parallel.

Non-lactose fermenting bacteria produce clear colonies, transparent or colorless, while coliforms are sufficiently inhibited, and form small colonies that vary from rose to red in colour. The H₂S producing bacteria produce colonies with black centers and a clear halo such as Proteus and other species of Salmonella.

MICROBIOLOGICAL TEST

The following results were obtained in the performance of the given medium from type cultures after anaerobic incubation at a temperature of 35 ± 2 °C and observed after 18 – 24 hours.

Microorganisms	Growth	Colony colour
<i>Salmonella enteritidis</i> ATCC 13076	Good	Colourless with black center
<i>Salmonella typhi</i> ATCC 6539	Good	Colourless with black center
<i>Salmonella typhimurium</i> ATCC 14028	Good	Colourless with black center
<i>Shigella flexneri</i> ATCC 12022	Good	Colourless
<i>Enterobacter aerogenes</i> ATCC 13048	Partially inhibited	Cream-pink
<i>Escherichia coli</i> ATCC 25922	Inhibited	-
<i>Enterococcus faecalis</i> ATCC 19433	Inhibited	-

SALMONELLA-SHIGELLA AGAR

500 g

X934.1

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe • P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0 • Fax: +49 (0) 721/ 5606-149 • info@carlroth.com • www.carlroth.com

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

sse 07/2021