KRISTALLVIOLETT-GALLE-LACTOSE AGAR

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
VIOLET RED BILE AGAR WITH LACTOSE

For the detection and enumeration of coliforms in milk, foods and other materials
 Acc. to ISO 4832
 X940

Formulation in g/l:

| | |
|-------------------------|-----------|
| Yeast Extract | 3.0 |
| Bile Salts..... | 1.5 |
| Sodium Chloride..... | 5.0 |
| Neutral Red | 0.03 |
| Peptone ex Gelatin..... | 7.0 |
| Lactose | 10.0 |
| Crystal Violet | 0.002 |
| Agar | 15.0 |
| pH | 7.4 ± 0.2 |

Escherichia coli ATCC 25922 *Salmonella gallinarum* NCTC 9240

PREPARATION

Suspend 41.5 g of the medium in one liter of deionized or distilled water. Mix well and heat with frequent agitation. Boil for one minute until complete dissolution. If desired it can be sterilized at 118 °C for 10 minutes. But in the studies of Hartman, he demonstrated that media prepared only by boiling gave the same results as media sterilized by autoclaving (Hartman, J., Milk and Food Tech. 23:43, 1980).
 Cool to 45-50 °C and pour into plates.

USES

Violet Red Bile Agar with Lactose is a differential and mildly selective medium based on MacConkey Agar for the detection of coliforms in water as well as milk and other food materials, as well as animal feeding stuffs. It may also be used according to the APHA Standard Methods for the examination of Milk Products. ISO 4832 recommends Violet Red Bile Agar with Lactose for colony count of Coliforms from foods and animal feeding stuffs.

Gram-positive organisms are markedly inhibited by the bile salts and the crystal violet. The colonies of lactose-fermenting bacteria are red in color whose size depends on the number of colonies on the plate. Occasionally the cocci of the intestinal tract can develop as small, punctiform red colonies.

The material sample is seeded in small aliquots immediately onto the Agar. If desired, after the plates have solidified and been stored, but before the sample is seeded, another thin layer can be poured on top. Some laboratories are accustomed to this method and dismiss any growth on the lower layer as contamination.

MICROBIOLOGICAL TEST

The following results were obtained in the performance of the given medium from type cultures after incubation at a temperature of 35 ± 2 °C and observed after 18 – 24 hours.

| Microorganisms | Growth | Colony colour |
|--|-----------|--------------------|
| <i>Escherichia coli</i> ATCC 25922 | Good | Purple |
| <i>Enterobacter aerogenes</i> ATCC 13048 | Good | Purple |
| <i>Salmonella gallinarum</i> NCTC 9240 | Good | Colourless |
| <i>Staphylococcus aureus</i> ATCC 6538 | Inhibited | - |
| <i>Enterococcus faecalis</i> ATCC 29212 | Inhibited | - |
| <i>Pseudomonas aeruginosa</i> ATCC 27853 | Good | Colourless - beige |

VIOLET-RED-BILE AGAR WITH LACTOSE

500 g

X940.1

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe • P.O. Box 100121 • 76231 Karlsruhe
 Phone: +49 (0) 721/ 5606-0 • Fax: +49 (0) 721/ 5606-149 • info@carlroth.com • www.carlroth.com

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

sse 06/2021