

Transfection Reagent: Roti[®]-Fect (P001)

Cells / Cell line: Panc-1 (Human epithelioid carcinoma cell line)

Protocols and Results



A)

Cell density	12wells, 2x10 ⁵ cells / ml, 24 h prior to transfection
Medium	x
DNA / RNA	Plasmid with GFP
Transfection complex	1,1 µg Plasmid-DNA + 50 µl serum free medium + 3 µl Roti [®] -Fect plus in serum free medium Formation of transfection complex for 20 mins. at ambient temp.
Transfection	in 1 ml medium + FCS over night
Incubation	in full medium (37 °C, 5 % CO ₂) for 48 hours overall
Detection	FACS and fluorescence microscopy
Result	Rot [®] -Fect is well tolerated by the cell line with very low to no toxicity. (Fugene HD, used in parallel with the same parameters, did show significant toxicity.) Transfection rate was 4.0 %. (Fugene HD 2.0 %)

13.07.2016
Carl Roth GmbH + Co. KG
i.V. Dr. Stefanie Seipp