



INSTRUCTIONS FOR USE

SPONGES

TRADE NAME	Sponges for histology cassettes
IDENTIFICATION	SPONGES FOR HISTOLOGY CASSETTES
INTENDED USE	Device for professional laboratory analysis use with the function of collecting and protecting biological samples to be inserted into histology cassettes for the performance of laboratory analysis
NON AUTOMATED DEVICE	-----
QUANTITATIVE DEVICE	The sponges are used to contain and protect small sized samples. The results of analyses conducted must always be subject to the scrutiny of highly qualified personnel for validation and reporting
SAMPLE	Biological sample or sample of another nature
IN VITRO MEDICAL DIAGNOSTIC DEVICE	Compliant with (EU) regulation 2017/746 for in vitro medical-diagnostic devices
USERS	Highly qualified personnel. For professional laboratory analysis use
STORAGE	Protect from direct sunlight and UV rays. Store in a dry location at room temperature
MATERIAL STABILITY	The characteristics of the material and the resistance values allow the product to remain stable over time. It is established that these types of devices are valid for 5 years from the production date listed on the label.
USAGE INSTRUCTIONS AND RESTRICTIONS	comply with laboratory sampling methods and/or environmental protocols in which the item is used report any defects that prevent usage of the sponges to the product distributor report serious accidents to the proper authorities dispose of sponges according to the applicable law. wear gloves and personal protective equipment when handling potentially contaminating samples. Keep out of children's reach Do not use the item after the expiration date Verify the resistance to specific substances in the Technical/Chemical Resistance section of the current catalogue.
SINGLE USE	② do not re-use. Danger of contamination



Kartell SpA
Via delle Industrie 1
20082 Noviglio (MI)
Tel. 02/900121
Website: www.kartellabware.com



to access instructions in other languages: www.kartellabware.com/en/download/