

**INSTRUCTIONS FOR USE**

**SAMPLE CONTAINERS:**

TRADE NAME	Single use containers for saliva, faeces, urine, cylindrical vessels and metallic trays for histology samples, histology cassettes, rings for the inclusion and Petri dishes.
IDENTIFICATION	<b>SINGE USE CONTAINERS FOR BIOLOGICAL SAMPLES</b>
INTENDED USE	Device for professional laboratory analysis use <b>for the collection of biological samples</b> such as urine, saliva, faeces, histology samples for the purpose of conducting analyses for the measurement of the chemical, physical and microscopic characteristics of the samples obtained. Containers used as instruments for the growth of cell cultures or in daily use in the laboratory as sample containers.
NON AUTOMATED DEVICE	-----
QUANTITATIVE DEVICE	The results of analyses conducted on the samples contained in the device must always be subject to the scrutiny of highly qualified personnel for validation of the report.
SAMPLE	Biological sample or sample of another nature
IN VITRO MEDICAL DIAGNOSTIC DEVICE	Compliant with (EU) regulation 2017/746 for in vitro medical-diagnostic devices
USERS	Highly qualified personnel. For professional laboratory analysis use
STORAGE	Protect from direct sunlight and UV rays. Store in a dry location at room temperature
MATERIAL STABILITY	The plastic materials used allow the product to remain stable over time. It is established that these types of devices are valid for 5 years from the production date listed on the label.
USAGE INSTRUCTIONS AND RESTRICTIONS	<p>comply with laboratory and/or environmental protocols in which the item is used for the execution of the analyses.</p> <p>If the container equipped with a cap, once it has been opened to collect the sample, place the cover in a hygienic place with the internal part facing upwards.</p> <p>Ensure that the internal part of the cover is not touched or contaminated in any way. Immediately cap the container after collecting the sample to be examined</p> <p>Report any defects that prevent use of the containers to the product distributor.</p> <p>Report serious accidents to the proper authorities.</p> <p>Do not approach open flames or heat sources</p> <p>Dispose of the containers according to the applicable law</p> <p>Wear gloves and personal protective equipment when handling potentially contaminating solutions.</p> <p>Keep out of children's reach</p>
SINGLE USE OR REUSABLE DEVICES	② do not re-use. Danger of contamination



Kartell SpA  
 Via delle Industrie 1  
 20082 Noviglio (MI)  
 Tel. 02/900121  
 Website: [www.kartelllabware.com](http://www.kartelllabware.com)



to access instructions in other languages:  
[www.kartelllabware.com/en/download/](http://www.kartelllabware.com/en/download/)