

# EU Declaration of Conformity



Product: Laboratory Equipment

File Reference: TFR926061\_001

Manufacturer: Antylia Scientific Ltd  
Beacon Road  
Stone, Staffordshire  
ST15 0SA  
United Kingdom

EU Representative

Antylia Scientific GmbH  
Futtererstraße 16  
97877 Wertheim  
Deutschland

This declaration of conformity is issued under the sole responsibility of the manufacturer

## Object of Declaration

**Stuart 150 mm Hotplate / stirrers  
04807-xx 'Undergrad' Range**

The object of the declaration described above is in conformity with the relevant Union Harmonisation Legislation

**Electromagnetic Compatibility Directive**  
**Low Voltage Directive**  
**RoHS Directive**

**2014/30/EU**  
**2014/35/EU**  
**2011/65/EU (including CDD 2015/863)**

References to the relevant harmonised standards used or references to the other technical specifications in relation to which conformity is declared:

EN 61326-1:2013

Electrical equipment for the measurement, control and Laboratory use – Part 1: General requirements

EN 61010-1:2010  
/A1:2019/AC:2019-04

Safety requirements for electrical equipment for measurement, control and laboratory use – Part 1: General requirements

EN 61010-2-010:2020  
(Applicable to unit with heating functionality only)

Safety requirements for electrical equipment for measurement, control and laboratory use – Part 2-010: Particular requirements for laboratory equipment for the heating of materials

EN 61010-2-051:2015  
(Applicable to unit with stirring functionality only)

Safety requirements for electrical equipment for measurement, control and laboratory use – Part 2-051: Particular requirements for laboratory equipment for mixing and stirring

Signed for and on behalf of the above manufacturer

Year of Marking  
Place of Issue  
Authorised Representative

2021

Stone, UK

Martin Card

Title

EMEA Quality and Engineering Manager

Signature

Date

21st April 2021

Catalogue Code	Part Code	Description
04807-50	UC150	Analogue hotplate, ceramic, 150mm x 150 mm 230 V a.c.
04807-52	US150	Analogue hotplate, aluminium, 150mm x 150 mm 230 V a.c.
04807-54	UC151	Analogue stirrer, ceramic, 150mm x 150 mm 230 V a.c.
04807-56	US151	Analogue stirrer, stainless steel, 150mm x 150 mm 230 V a.c.
04807-58	UC152	Analogue stirring hotplate, ceramic, 150mm x 150 mm 230 V a.c.
04807-60	US152	Analogue stirring hotplate, aluminium, 150mm x 150 mm 230 V a.c.
04807-62	UC152D	Digital stirring hotplate, ceramic, 150mm x 150 mm 230 V a.c.
04807-64	US152D	Digital stirring hotplate, aluminium, 150mm x 150 mm 230 V a.c.